## HOME BLOOD PRESSURE KIT INSTRUCTIONS – MODEL UA-101

### IMPORTANT INFORMATION

#### PLEASE READ THIS IMPORTANT INFORMATION BEFORE USING YOUR MONITOR.


- Please remember that only a medical practitioner is qualified to interpret your blood pressure measurements. Use of this device should not replace regular medical examinations.
- Have your physician review the procedure for using this device. He or she will want to verify blood pressure measurements before prescribing or adjusting medications.
- Consult your physician if you have any doubt about your readings. Should a mechanical problem occur, contact A&D Medical.
- Clean the monitor and cuff with a dry, soft cloth or a cloth dampened with water and a mild detergent. Never use alcohol, benzene, thinner or other harsh chemicals to clean monitor or cuff.

### PRECAUTIONS

This Blood Pressure monitor is designed to be used at home or in a doctor's/nurse's office, by those who are eighteen (18) years and older, to monitor blood pressure (systolic and diastolic) and pulse rate.

**Caution:** This product contains natural rubber latex which may cause allergic reactions.

### HOME BLOOD PRESSURE KIT COMPONENTS

#### LIST OF COMPONENTS:

|                   |                    |                  |
|-------------------|--------------------|------------------|
| Ear tips          | Left Bladder Tube  | Stethoscope Head |
| Binaurals         | Right Bladder Tube | Carrying Case    |
| Inflation Bulb    | Y-Tube             | Calibration key  |
| Air Release Valve | Cuff D-Ring        |                  |
| Manometer         | Cuff               |                  |

### ASSEMBLY INSTRUCTIONS

- Place cuff on your upper arm.
- Make sure cuff end has been inserted through the D-Ring.
- Insert single end of the Y-Tube into stethoscope head.
- Insert each arm of the Y-Tube into each bottom end of the binaurals.
- Insert air release valve into left bladder tube.
- Insert manometer into right bladder tube.

### TAKING YOUR BLOOD PRESSURE

#### BLOOD PRESSURE TAKING TIPS

- Relax for several moments before measurement.
- Do not smoke or ingest caffeine at least 30 minutes prior to measurement.
- Remove constricting clothing and place cuff on bare arm.
- Unless your physician recommends otherwise, use left arm to measure pressure.
- Do not talk during measurement.
- Do not cross legs and keep your feet flat on the floor.

#### NOW YOU ARE READY. FOLLOW THESE SIMPLE STEPS:

- Sit comfortably with your left arm resting on a flat surface so that the center of your upper arm is at the same height as your heart.
- Lay left arm on the table with palm up. Place the cuff on upper arm with the tube facing downward and toward the inside of your arm. Fasten cuff securely with hook and loop fastener strap.
- The cuff should be snug but not too tight. You should be able to insert two fingers between the cuff and your arm.
- With your palm facing up, place the stethoscope head one-inch above the bend of the elbow on the inside of the arm. This is where the main artery is located.
- Insert the ear tips of the stethoscope into your ears. Adjust by turning ear tips to fit comfortably in your ears.
- Hold the manometer in your left hand and the inflation bulb in your right hand.
- Close the air release valve on the bulb by turning the valve clockwise.
- Inflate the cuff by repeatedly squeezing the bulb with your right hand. Listen to the pulse beat while inflating the cuff.
- Watch the gauge. When you can no longer hear the pulse beat, raise the pressure an additional 30 mmHg.
- After reaching your inflation level, slowly deflate the cuff by opening the air release valve. To open the air release valve, turn it counter-clockwise so that the pressure drops 2-4 mmHg per second. The needle on the manometer should drop an average of one to two marks every second.
- After opening the air release valve, listen carefully for a pulse beat. The moment you hear the first thumping sounds of your pulse beat, note the reading on the gauge. This is your systolic blood pressure measurement.
- When you can no longer hear any sounds, read the gauge. This is your diastolic blood pressure measurement.
- Remove cuff and make note of your blood pressure.

**NOTE:** To be certain your technique was good and the readings were as accurate as possible, wait 10-15 minutes and repeat the procedure.

#### TWO-PARTY METHOD:

- If you plan on taking someone else's blood pressure, have the person whose blood pressure is to be taken to rest and relax at least 15 minutes beforehand.
- Have the person extend his or her left arm toward you, palm up. Be sure the arm rests at the same level as the heart.
- Position the cuff as indicated in the self-taking instructions.
- Attach the manometer to the loop on the cuff, so that it is clearly visible to you.
- Perform Steps 1-13 in Taking Your Blood Pressure. If possible, have your physician watch your technique to be sure it is correct.

### ABOUT BLOOD PRESSURE

#### WHAT IS BLOOD PRESSURE?


Blood pressure is the force exerted by blood against the walls of the arteries. Systolic pressure occurs when the heart contracts; diastolic pressure occurs when the heart expands. Blood pressure is measured in millimeters of mercury (mmHg).

#### WHAT AFFECTS BLOOD PRESSURE?

Blood pressure is affected by many factors: age, weight, time of day, activity level, climate, altitude and season. Certain activities can significantly alter one's blood pressure. Walking will likely raise an individual's blood pressure. Sleeping will likely decrease blood pressure. Not relaxing for several moments before measurements can influence readings as well.

In addition to these factors, beverages containing caffeine or alcohol, certain medications, emotional stress and even tight-fitting clothes can make a difference in the readings.

#### Fluctuation within a day (case: 35 year old male)

#### WHAT CAUSES VARIATIONS IN BLOOD PRESSURE?

An individual's blood pressure varies greatly from day to day and season to season. For hypersensitive individuals, these variations are even more pronounced. Normally, blood pressure rises during work or play and falls to its lowest levels during sleep.

#### ASSESSING HIGH BLOOD PRESSURE FOR ADULTS

The following standards for assessing high blood pressure (without regard to age or gender) have been established as a guideline. Please note that other risk factors (e.g. diabetes, obesity, smoking, etc.) need to be taken into consideration and may affect these figures. Consult with your physician for an accurate assessment.

#### JNC7 Classification Table – for adults within the U.S.

| BP Classification    | Systolic (mmHg)              | and | Diastolic (mmHg)             |
|----------------------|------------------------------|-----|------------------------------|
| Normal               | less than 120                |     | less than 80                 |
| Prehypertension      | 120-139                      | or  | 80-89                        |
| Stage 1 Hypertension | 140-159                      | or  | 90-99                        |
| Stage 2 Hypertension | greater than or equal to 160 | or  | greater than or equal to 100 |

**SOURCE:** The Seventh Report of the Joint National Committee on Prevention, Evaluation and Treatment of High Blood Pressure for adults. National Heart, Lung and Blood Institute - May 2003.

#### WHO Classification Table – for adults outside of the U.S. (e.g. Canada, Mexico)

| BP Classification    | Systolic (mmHg)              | and | Diastolic (mmHg)             |
|----------------------|------------------------------|-----|------------------------------|
| Optimal              | less than 120                |     | less than 80                 |
| Normal               | less than 130                | or  | less than 85                 |
| High-Normal          | 130-139                      | or  | 85-89                        |
| Grade 1 Hypertension | 140-159                      | or  | 90-99                        |
| Grade 2 Hypertension | 160-179                      | or  | 100-109                      |
| Grade 3 Hypertension | greater than or equal to 180 | or  | greater than or equal to 110 |

**SOURCE:** Standards to assess high blood pressure, without regard to age or gender, have been established by the World Health Organization (WHO).

#### WHAT IS HYPERTENSION?

Hypertension (high blood pressure) is the diagnosis given when readings consistently rise above normal. It is well known that hypertension can lead to stroke, heart attack or other illness if left untreated. Referred to as a "silent killer" because it does not always produce symptoms that alert you to the problem, hypertension is treatable when diagnosed early.

#### CAN HYPERTENSION BE CONTROLLED?

In many individuals, hypertension can be controlled by altering lifestyle and minimizing stress, and by appropriate medication prescribed and monitored by your doctor. The American Heart Association recommends the following lifestyle suggestions to prevent or control hypertension:

- Don't smoke.
- Reduce salt and fat intake.
- Exercise routinely.
- Have regular physical checkups.
- Maintain proper weight.
- Monitor your blood pressure at periodic intervals.

#### WHY MEASURE BLOOD PRESSURE AT HOME?

It is now well known that, for many individuals, blood pressure readings taken in a doctor's office or hospital setting might be elevated as a result of apprehension and anxiety. This response is commonly called "white coat hypertension."

In any case, self-measurement at home supplements your doctor's readings and provides a more accurate, complete blood pressure history. In addition, clinical studies have shown that the detection and treatment of hypertension is improved when patients both consult their physicians and monitor their own blood pressure at home.

### MAINTENANCE INSTRUCTIONS

#### BASIC CARE OF YOUR HOME BLOOD PRESSURE KIT

With proper care, your Home Blood Pressure Kit will provide many years of reliable service. To insure a well functioning unit, follow these basic rules:

- Handle the pressure manometer with care. Do not drop or jar. The needle should indicate zero when cuff is fully deflated.
- Always deflate cuff before storing.
- Avoid hitting objects against surface of stethoscope.
- Store all items, including instruction booklet, in the convenient storage case provided.

A&D Medical has checked your blood pressure unit to assure reliability and accuracy prior to shipment. However, as with any sensitive instrument subject to repeated use, we recommend your blood pressure unit be checked periodically for proper calibration to ensure accurate blood pressure readings.

### CHECKING FOR ACCURACY

#### HOW TO CHECK YOUR HOME BLOOD PRESSURE KIT FOR ACCURACY

- Detach manometer from bladder tube. Hold manometer in a vertical position.
- If the needle on the dial rests within the zero indicator on the face plate, your unit is guaranteed to be accurate.
- If the needle rests outside the zero indicator, it should be recalibrated.
- To recalibrate, insert the calibration key into the hose connector end of the manometer. Turn the calibration key in either direction without going all the way around the dial to get the needle back to the zero indicator.
- If the unit cannot be recalibrated using this method, please contact A&D Medical at 1-888-726-9966.

VISIT OUR WEBSITE AT [WWW.ANDMEDICAL.COM](http://WWW.ANDMEDICAL.COM)  
FOR ADDITIONAL SUPPORT.

### CONTACT INFORMATION

Our products are designed and manufactured using the latest scientific and technological methods, and offer accurate, easy to use, home monitoring and treatment options. Our full product line includes:

- Digital Blood Pressure Monitors
- Manual Blood Pressure Monitors
- Blood Pressure Cuffs
- Wireless Products
- Personal Health Scales
- Activity Monitors
- Digital Thermometers
- Stethoscopes

This A&D blood pressure product is covered by a two year warranty.

Distributed by:  
**A&D Engineering, Inc.**  
1756 Automation Parkway  
San Jose, CA 95131  
1-888-726-9966  
[www.andmedical.com](http://www.andmedical.com)

LifeSource is a trademark of A&D Medical. ©2012 A&D Medical. All rights reserved. Specifications subject to change without notice. I-MAN-UA-101-062012-REV-a Printed in China

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