TALKING PREMIUM+

BLOOD PRESSURE MONITOR
UPPER ARM
UA-1030T

Instruction Manual
ORIGINAL

Greetings

Congratulations on purchasing a state-of-the-art A&D blood pressure monitor. Designed for ease of use and accuracy, this monitor will facilitate your daily blood pressure regimen. We recommend that you read through this manual carefully before using the monitor for the first time.

Intended Use

- The monitor is designed for use on adults only. Do not use on newborns or infants.
- Environment for use. The monitor is intended for use in the home healthcare environment.
- This monitor is designed to measure blood pressure and pulse rate of people for diagnosis.

Precautions

- Precision components are used in the construction of this monitor. Extremes in temperature, humidity, direct sunlight, shock or dust should be avoided.
- Clean the monitor and cuff with a dry, soft cloth or a cloth dampened with water and a neutral detergent. Never use alcohol, benzene, thinner or other harsh chemicals to clean the monitor or cuff.
- Avoid tightly folding the cuff or storing the hose tightly twisted for long periods, as such treatment may shorten the life of the components.
- The monitor and cuff are not water resistant. Prevent rain, sweat and water from soiling the monitor and cuff.
- Measurements may be distorted if the monitor is used close to televisions, microwave ovens, cellular telephones, X-ray or other devices with strong electrical fields.
- When using the monitor, confirm that the monitor is clean.
- Used equipment, parts and batteries are not treated as ordinary household waste, and must be disposed of according to the applicable local regulations.

Precautions

- When the AC adapter is used, make sure that the AC adapter can be readily removed from the electrical outlet when necessary.
- Do not modify the monitor. It may cause accidents or damage to the monitor.
- To measure blood pressure, the arm must be squeezed by the cuff hard enough to temporarily stop blood flow through the artery. This may cause pain, numbness or a temporary red mark to the arm. This condition will appear especially when measurement is repeated successively. Any pain, numbness, or red marks will disappear with time.
- Do not apply the cuff on an arm in which another medical device is attached. The equipment may not function properly.
- People who have a severe circulatory deficit in the arm must consult a doctor before using the monitor, to avoid medical problems.
- Do not self-diagnose the measurement results and start treatment by yourself. Always consult your doctor for evaluation of the results and treatment.
- Do not apply the cuff on an arm with an unhealed wound.
- Do not apply the cuff on an arm receiving an intravenous drip or blood transfusion. It may cause injury or accidents.
- Do not use the monitor where flammable gases such as anesthetic gases are present. It may cause an explosion.
- Do not use the monitor in highly concentrated oxygen environments, such as a high-pressure oxygen chamber or an oxygen tent. It may cause a fire or explosion.
- Take care to avoid accidental strangulation of babies or infants with the hose and cable.
- Do not twist the air hose during measurement. This may cause injury due to continuous cuff pressure.
- Wireless communication devices, such as home networking devices, mobile phones, cordless phones and their base stations, walkie-talkies can affect this blood pressure monitor. Therefore, a minimum distance of 30 cm (12") should be kept from such devices.

Precautions

- Measuring blood pressure too frequently may cause harm due to blood flow interference. Check that the operation of the monitor does not result in prolonged impairment of blood circulation, when using the monitor repeatedly.
- Clinical testing has not been conducted on newborn infants and pregnant woman.
- If you have had a mastectomy, please consult a doctor before using the monitor.
- Do not let children use the monitor by themselves and do not use the monitor in a place within the reach of infants. It may cause accidents or damage.
- There are small parts that may cause a choking hazard if swallowed by mistake by infants.
- Do not touch the batteries, the DC jack and the patient at the same time. That may result in electrical shock.
- Unplug the AC adapter when not in use during the measurement.
- Use of accessories not detailed in this manual may compromise safety.
- Should the battery short-circuit, it may become hot and potentially cause burns.
- Allow the monitor to adapt to the surrounding environment before use (about one hour).
- Do not inflate without wrapping the cuff around the upper arm.

Contents


Parts Identification	E-6
Symbols	E-7
Using the Monitor	E-10
Taking your Blood Pressure	E-16
Measurements	E-17
About Memory	E-21
What is an Irregular Heartbeat.....	E-22
About Cuff Inflation Meter	E-22
AHA Classification Indicator	E-23
About Blood Pressure	E-24
Troubleshooting	E-26
Technical Data	E-27
Warranty	E-32

Symbols

Symbols that appear on the display


Symbols	Function/Meaning	Recommended Action
	Appears while measurement is in progress. It blinks when the pulse is detected.	Measurement is in progress. Remain as still as possible.
	Irregular Heartbeat symbol (I.H.B.)Appears when an irregular heartbeat is detected during measurement	_____
	TriCheck mode Takes three consecutive measurements automatically and displays the average values of the three measurements.	_____
	Appears when a body or arm movement is detected.	The reading may yield an incorrect value. Take another measurement. Remain still during measurement.
	Appears during measurement when the cuff is applied loosely.	The reading may yield an incorrect value. Apply the cuff correctly, and take another measurement.
	Previous measurements stored in memory	_____
	Average data	_____
	The battery power indicator during measurement.	_____
	The battery is low when it blinks.	Replace all batteries with new ones when the mark blinks.

Symbols

Symbols	Function/Meaning	Recommended Action
Err	Unstable blood pressure due to movement during measurement.	Take another measurement. Remain very still during measurement.
	The systolic and diastolic values are within 10mmHg of each other.	Apply the cuff correctly, and take another measurement.
	The pressure value did not increase during the inflation.	
Err CUF	The cuff is not applied correctly.	
E	PUL DISPLAY ERROR The pulse is not detected correctly.	
Err E	Blood pressure monitor internal error	Remove the batteries and press the START button, and then install the batteries again. If the error still appears, contact the dealer.
Err S		
SYS	Systolic blood pressure in mmHg	_____
DIA	Diastolic blood pressure in mmHg	_____
PUL/min	Pulse per minute	_____
AM	Data taken in the AM	_____
PM	Data taken in the PM	_____
	Pressure settings Indicates the pressure value as set by the user.	_____

Symbols

Symbols that are printed on the device case

Symbols	Function/Meaning
	Standby and Turn the device on.
	Battery installation guide
	Direct current
	Type BF: Device, cuff and tubing are designed to provide special protection against electrical shocks.
SN	Serial number
	Refer to instruction manual/booklet
	Polarity of DC jack
IP	International Protection Symbol
	Keep dry

Symbols that are printed on the AC adapter


Symbols	Function/Meaning
	Alternating current
	Direct current
	Polarity of DC jack
	Class II device
	For indoor use only
	UL/cUL certification device label
	Energy efficiency level VI certification device label
IP	International protection symbol

Using the Monitor

Installing/Changing the Batteries

You must install 4 type A (1.5 volt) batteries (alkaline batteries recommended), or use the AC Adapter and attach the cuff to the monitor before using it. To install batteries (or replace them if the "Low Battery" symbol appears on display), proceed as follows:

1. Remove battery compartment cover by gently pushing down on arrow and sliding cover forward.
2. Put in bottom row of batteries first. Place the batteries in compartment with positive (+) and negative (-) terminals matching those indicated in the compartment. Be sure batteries make contact with compartment terminals.
3. Replace cover by sliding it into the compartment and gently pressing into place.

Note: Rechargeable batteries are not recommended for use with this monitor. Please dispose of batteries following local regulations.

CAUTION


- Insert the batteries as shown in the battery compartment. If installed incorrectly, the device will not work.
- When  (LOW BATTERY mark) blinks on the display and the device announces that the battery needs to be replaced, replace all batteries with new ones. Do not mix old and new batteries. It may shorten the battery life, or cause the device to malfunction. Replace the batteries two seconds or more after the device turns off. If  (LOW BATTERY mark) appears even after the batteries are replaced, make a blood pressure measurement. The device may then recognize the new batteries.
-  (LOW BATTERY mark) does not appear when the batteries are drained.

Using the Monitor

- The battery life varies with the ambient temperature and may be shorter at low temperatures. Generally, four new LR6 batteries will last approximately for one year when used twice for measurement each day.
- Use the specified batteries only. The batteries provided with the device are for testing monitor performance and may have a limited life.
- Remove the batteries if the device is not to be used for a long time. The batteries may leak and cause a malfunction.

Connecting the AC adapter to the monitor

Gently insert the AC adapter plug into a 120V AC outlet. Connect the AC adapter plug into the jack on the back of the monitor


Disconnecting the AC adapter from the monitor:

1. Turn the unit off by pressing the START button.
2. Disconnect the plug from the jack of the monitor quickly.
3. Gently remove the AC adapter from the outlet

Note: If monitor has batteries and an AC adapter connected, the measurements stored in memory may be lost if AC adapter is unplugged from the wall first. To ensure that measurements remain in memory, make sure the unit is off and that the jack plugged into the monitor is removed quickly before unplugging adapter from the wall.

Connecting the Air Hose

Insert the air connector plug into the air socket firmly.


Using the Monitor

Adjusting the Clock

The UA-1030T blood pressure monitor has a built-in clock that gives you the time and date and can tell you when a blood pressure reading stored in memory was taken. We recommend that you adjust the clock prior to use.

To adjust the clock, follow these simple steps:

1. Press the set button (Ⓢ) located to the right of the display.
2. Use the arrow ▲▼ buttons to set the year.
3. After reaching the desired year, press the set button (Ⓢ) again.
4. Repeat the process for the month, day and hour.
5. The next setting will be for pressure.

Set Button

Adjusting the Pressure

1. After setting the year, month, day and hour, the next mode is the pressure setting.
2. Default is **AUTO**, which stands for automatic.
3. For even more rapid readings, you can adjust the default to about 30 mmHg or more above your expected systolic pressure.
4. Use the arrow ▲▼ buttons to set the pressure.
5. After reaching the desired pressure, press the set button (Ⓢ) again.
6. The next setting will be for TriCheck™ Mode.

Using the Monitor

Note: Pressing  will toggle between volume setting → language setting → instruction setting → original state (standby) .

Adjusting the Speaker Volume

1. Press the  button to enter the volume setting mode.
2. There will be a confirmation tone and the current volume will be displayed by the pressure bar (volume off, 1, 2, 3, 4, 5, 6). If it is off , the LCD will display “ oFF ”.
3. Change the volume value with the ▲ ▼ buttons.

Selecting the Language

1. Press  to enter the language setting mode. L1, L2 or L3 will be displayed.
L1: English
L2: French
L3: Spanish
2. Change the language with the ▲ ▼ buttons.

Instructions On / Off

1. Press  to enter the instruction setting mode.
2. Turn the audio ON and OFF with the ▲ ▼ buttons.

Using the Monitor

Tricheck™ Mode

The UA-1030T blood pressure monitor has an optional feature to automatically take three consecutive readings one minute apart and then provide you with the average. This complies with the American Heart Association's recommended measuring method.

1. After setting the year, month, day, hour and pressure, the next mode is the  Tricheck™ Mode.
2. Default is off.
3. To turn on, use the arrow ▲ ▼ buttons.
4. To turn off again, use the arrow ▲ ▼ buttons.
5. After reaching the desired setting, press the set  button again.
6. This will end your settings session.
7. To adjust any values, push the set  button and begin again.

Using the Monitor

Selecting the Correct Cuff Size


Using the correct cuff size is important for an accurate reading. If the cuff is not the proper size, the reading may yield an incorrect blood pressure value.

- The arm size is printed on each cuff.
- The index ▲ and proper fit range, on the cuff, tell you if you are applying the correct cuff. (Note: the SmoothFit Cuff, UA-392A, does not have the index symbol)
- If the index ▲ points outside of the range, contact your local dealer to purchase a replacement cuff. The arm cuff is a consumable. If it becomes worn, purchase a new one.

Catalog Number	Recommended Cuff Size	Arm Size
UA-289A	Small Cuff	16-24cm
UA-290A	Medium Cuff	23-37cm
UA-291A	Large Cuff	31-45cm
UA-392A	SmoothFit™ Cuff	23-37cm
UA-420A	Wide Range Cuff	22-42 cm

Applying the arm cuff

1. Wrap the cuff around the upper arm, about 1-2 cm above the inside of the elbow, as shown in the figure to the right.
2. Place the cuff directly against the skin. Clothing may cause a faint pulse, and result in a measurement error. Constriction of the upper arm, caused by a tightly rolled up shirtsleeve, may prevent accurate readings.

Taking your Blood Pressure

This blood pressure monitor is designed to detect the pulse and to inflate the cuff to a systolic pressure level automatically.

Tips for Taking Blood Pressure

Before Your Measurement:

- For 30 minutes prior to taking your blood pressure
 - » Do not exercise
 - » Do not drink coffee, caffeinated soda or alcohol
 - » Do not smoke
- Sit quietly for 5-10 minutes
- Ensure you are using the right size cuff

During Your Measurement:

- Do not talk
- Sit with your back straight and supported
- Uncross your legs and place feet flat on floor
- Rest arm on a table so that the cuff is at heart level
- Measure two times a day, in the morning and evening

Notes for Accurate Measurement

- This monitor bases its measurements on the heartbeat. If you have a very weak or irregular heartbeat, the monitor may have difficulty determining your blood pressure.
- Should the monitor detect a condition that is abnormal, it will stop the measurement and display an error symbol. Refer to the section "Symbols" for the description of the symbols.
- If you have emotional stress, the measurement will reflect this stress as a higher (or lower) than normal blood pressure reading and the pulse reading will usually be faster than normal.
- An individual's blood pressure varies constantly, depending on what you are doing and what you have eaten. What you drink can have a very strong and rapid effect on your blood pressure.

Measurements

1. Press the START button. All of the display segments are displayed. Zero (0) is displayed blinking briefly. The display changes, as indicated in the figure at the right, as the measurement begins.

The cuff starts to inflate. It is normal for the cuff to feel very tight. A pressure bar indicator is displayed, on the left edge of the display, during the inflation.


Note: If you wish to stop inflation at any time, press the START button again

2. When inflation is complete, deflation starts automatically and ♥ blinks, indicating that the measurement is in progress. Once the pulse is detected, the mark blinks with each pulse beat.

Note: If an appropriate pressure is not obtained, the device starts to inflate again automatically. To avoid re-inflation, see "Measurement with Set Pressure".

3. The display measurement is complete, the systolic and diastolic pressure readings and pulse rate are displayed. The device announces the measurement results. The cuff exhausts the remaining air and deflates completely


4. Press the START button to turn the device off. After one minute of non-operation, the device will turn off automatically.
5. The reading is then stored into memory. See "About Memory" for more information.

Measurements

Measurement using TriCheck

1. Press the START button. All of the display segments are displayed. Zero (0) is displayed blinking briefly and the first measurement starts.
2. When the measurement is complete, the systolic and diastolic pressure readings and pulse rate are displayed, then a one-minute countdown starts for the second measurement.
3. After one minute, the second measurement starts.
4. When the measurement is complete, the average readings of the first and second measurements are displayed, then a one-minute countdown starts for the third measurement.
5. After one minute, the third measurement starts.
6. When the measurement is complete, the average readings of the three measurements are displayed and stored in memory.

Notes: During the measurement, the TriCheck  symbol is displayed.

To cancel the measurement, press the START button. In this case, no data is stored in memory.

When the cuff fit error symbol appears during the first measurement, cancel the measurement, reapply the cuff correctly and start a new measurement.


After measurement, the average data of the three measurements is stored in memory. No data will be stored, when canceling the measurements, or before completion of the third measurement.

Measurements

Measurement with Set Pressure

The UA-1030T is designed to detect the pulse and to inflate the cuff to a systolic pressure level automatically. If re-inflation occurs repeatedly, use the following methods. During the blood pressure measurement, re-inflation may occur. A fixed pressure value can be set to avoid re-inflation.

1. At step 5 of the clock adjustment procedure (see page E-12), press the set button to go to the pressure setting mode. The current setting blinks.
2. Press the ▲ or ▼ button to select a pressure value about 30mmHg or more above your expected systolic pressure from the following.

AUTO	: Automatic pressurization (default value)
180	: Pressure value of 180mmHg (fixed)
210	: Pressure value of 210mmHg (fixed)
240	: Pressure value of 240mmHg (fixed)

3. Press the set button to go to the TriCheck setting mode.


Press the START button to turn the device off. After three minutes of non-operation, the device will turn off automatically. The next measurement will be performed with the new pressure value.

Measurements

Measurement with desired Systolic Pressure

Use this method when re-inflation occurs repeatedly even if the pressure value is set to 240 in the procedure above or when the results are not displayed even if the pressure decreases to 20mmHg or less.

1. Place the cuff on the arm (preferably the left arm).
2. Press and hold the START button until a number about 30 to 40mmHg higher than your expected systolic pressure appears.
3. When the desired number is reached, release the START button to start measurement. Continue to measure your blood pressure as described on the E-17.

About Memory

Note: This device stores the last 90 measurements in memory. The device announces the memory data values as they are displayed.

1. Press the ▼ or ▲ button. The average of all measurements and the number of data are displayed. (If no data, "0" is displayed. Press the ▼, ▲ or START button to turn the device off.)

2. Each time the ▼ button (or the ▲ button to display the data in the reverse order) is pressed, the memory data is displayed as follows.

Average of all AM (morning) measurements taken between 4:00 and 9:59. (In the example, 10 measurements. If no data, "--" is displayed.)


Average of all PM (evening) measurements taken between 18:00 and 1:59. (In the example, 9 measurements. If no data, "--" is displayed.)

Most recent data (See No. 35). Three seconds after the data number display, the measurement data is displayed.

Last data (No.1). Three seconds after the data number display, the measurement data is displayed.

3. After the last data is displayed, press the ▼ button to return the average of all measurements.

4. Press the START button to turn the device off. After one minute of non-operation, the device will turn off automatically.



What is an Irregular Heartbeat?

The UA-1030T blood pressure monitor provides a blood pressure and pulse rate measurement even when an irregular heartbeat occurs. The Irregular Heartbeat symbol (♥) will appear in the display window in the event an Irregular Heartbeat has occurred during measurement. An irregular heartbeat is defined as a heartbeat that varies from the average of all heartbeat intervals during the blood pressure measurement. It is important that you relax, remain still and refrain from talking during measurements.

Note: We recommend contacting your physician if you see this symbol frequently.

About Cuff Inflation Meter

The Cuff Inflation Meter is located on the left side of the display screen to tell you when the blood pressure monitor is inflating and deflating the cuff. The Cuff Inflation Meter moves up during inflation and moves down during deflation.

Inflation in Progress	Inflation Complete	Deflation/ Measurement in Progress
 The display shows a vertical bar on the left side that is approximately 1/3 full. The main display shows a large number '118' with 'mmHg' to its right. Below this, it shows '825' and '10:24' with 'AM' to the left of the time. A battery icon is in the bottom left corner.	 The display shows a vertical bar on the left side that is completely full. The main display shows a large number '150' with 'mmHg' to its right. Below this, it shows '825' and '10:24' with 'AM' to the left of the time. A battery icon is in the bottom left corner.	 The display shows a vertical bar on the left side that is approximately 1/3 full and has a downward-pointing arrow next to it. The main display shows a large number '90' with 'mmHg' to its right. Below this, it shows '825' and '10:24' with 'AM' to the left of the time. A heart symbol is in the bottom left corner. A battery icon is in the bottom left corner.

AHA Classification

Each segment of the bar indicator corresponds to the AHA blood pressure classification.

Example

Hypertension Stage 2	Hypertension Stage 1	Elevated
		

AHA Classification Indicator

-  ← Hypertension Stage 2
-  ← Hypertension Stage 1
-  ← Elevated
-  ← Normal

: The indicator displays a segment, based on the current data, corresponding to the AHA classification.

Blood Pressure Classification	Systolic (mmHg)		Diastolic (mmHg)
Normal	< 120	and	< 80
Elevated	120 – 129	and	< 80
Hypertension Stage 1 (High Blood Pressure)	130 – 139	or	80 – 89
Hypertension Stage 2 (High Blood Pressure)	≥ 140	or	≥ 90
HYPERTENSIVE CRISIS (consult your doctor immediately)	≥ 180	and / or	≥ 120

About Blood Pressure

What is Blood Pressure?

Blood pressure is the force exerted by blood against the walls of the arteries. Systolic pressure occurs when the heart contracts; diastolic pressure occurs when the heart expands. Blood pressure is measured in millimeters of mercury (mmHg).

What affects Blood Pressure?

Blood pressure is affected by many factors: age, weight, time of day, activity level, climate, altitude and season. Certain activities can significantly alter one's blood pressure. Walking will likely raise an individual's blood pressure. Sleeping will likely decrease blood pressure. Not relaxing for several moments before measurements can influence readings as well. In addition to these factors, beverages containing caffeine or alcohol, certain medications, emotional stress and even tight-fitting clothes can make a difference in the readings.

What causes Variations in Blood Pressure?

An individual's blood pressure varies greatly from day to day and season to season. For hypersensitive individuals, these variations are even more pronounced. Normally, blood pressure rises during work or play and falls to its lowest levels during sleep.

About Blood Pressure

What is Hypertension?

Hypertension (high blood pressure) is the diagnosis given when readings consistently rise above normal. It is well known that hypertension can lead to stroke, heart attack or other illness if left untreated. Referred to as a “silent killer” because it does not always produce symptoms that alert you to the problem, hypertension is treatable when diagnosed early.

Can Hypertension be Controlled?

In many individuals, hypertension can be controlled by altering lifestyle and minimizing stress, and by appropriate medication prescribed and monitored by your doctor. The American Heart Association recommends the following lifestyle suggestions to prevent or control hypertension:

- Do not smoke
- Exercise regularly
- Reduce salt and fat intake
- Have regular physical checkups
- Maintain proper weight

Why measure Blood Pressure at Home?

It is now well known that, for many individuals, blood pressure readings taken in a doctor’s office or hospital setting might be elevated as a result of apprehension and anxiety. This response is commonly called “white coat hypertension.” In any case, self-measurement at home supplements your doctor’s readings and provides a more accurate, complete blood pressure history. In addition, clinical studies have shown that the detection and treatment of hypertension is improved when patients both consult their physicians and monitor their own blood pressure at home.

Troubleshooting

Problem	Probable Cause	Corrective Action
Nothing appears on the display, even when the power is turned on.	Battery is drained.	Replace all batteries with new ones.
	Battery terminals are not in the correct position.	Reinstall the batteries with negative and positive terminals matching those indicated on the battery compartment.
The cuff does not inflate.	Battery voltage is too low.  (LOW BATTERY mark) blinks. If the batteries are drained completely, the mark does not appear.	Replace all batteries with new ones.
The device does not measure. Readings are too high or too low.	The cuff is not applied properly.	Apply the cuff correctly.
	Readings are too high or too low.	Make sure you remain very still and quiet during measurement.
	The cuff position is not correct.	Sit comfortably and still. Place your arm on a table with your palm facing upward and the cuff at the same level as your heart.
	_____	If you have a very weak or irregular heart beat, the device may have difficulty in determining your blood pressure.
Other	The value is different from that measured at a clinic or doctor's office.	Refer to "Why Measure Blood Pressure at Home?".
	_____	Remove the batteries. Place them back properly and take another measurement.

Note: If the actions described above do not solve the problem, contact the dealer. Do not attempt to open or repair this product, as any attempt to do so will make your warranty invalid.

Technical Data

Model	UA-1030T
Type	Oscillometric
Display	Digital character height Pressure (20.0mm)/Pulse (11.5mm) displayed simultaneously
Memory	90 readings
Measurement range	Pressure: 0 - 299mmHg Systolic pressure: 60 - 279mmHg Diastolic pressure: 40 - 200mmHg Pulse: 40 - 180 beats per minute
Accuracy	Pressure: ± 3 mmHg Pulse: $\pm 5\%$
Rating range	DC 6V 3W
Power source	4 x 1.5 volt batteries (AA, LR6 or R6P) or 120 volt AC adapter (TB-233)
Classification	Internally powered ME equipment (Supplied by batteries) /Class II (Supplied by adapter) Continuous operation mode.
Clinical test	According to ISO81060-2: 2013 In the clinical validation study, K5 was used on 85 subjects for determination of diastolic blood pressure.
EMD	IEC 60601-1-2: 2014
Operating Environment	50°F to 104°F (10°C to 40°C) 15%RH to 85%RH / 800hPa to 1060hPa
Transport/Storage Environment	-4°F to 140°F (-20°C to 60°C) 10%RH to 95%RH / 700hPa to 1060hPa
Dimensions	Width: 140mm; Height: 60mm; Length: 105mm

Technical Data

Weight	10.6 oz. (300g) without batteries
Ingress protection	Device: IP20
Number of measurements	1000 measurements, when AA alkaline batteries are used, with pressure value of 180mmHg at room temperature of 23°C.
Applied part.	Cuff Type BF 
Useful life	Device: 5 years (when used six times a day) Cuff: 2 years (when used six times a day)
AC adapter	5 years (when used six times a day)

Accessory AC adapter

The adapter is to connect the device to a power source at home. Please contact your local A&D dealer for purchasing. The AC adapter is required to be inspected or replaced periodically.

TB-233	Input: 120V~ 60Hz 0.15A Output: 6V  500mA      IP21
--------	---

Note: Specifications are subject to change without prior notice. IP classification is the degrees of protection provided by enclosures in accordance with IEC 60529. This device is protected against solid foreign objects of 12 mm diameter and greater such as a finger. This device is not protected against water.

Technical Data

EMD Technical Data

Battery-operated or AC Adapter-operated Blood Pressure Monitor

Medical Electrical Equipment needs special precautions regarding EMD and needs to be installed and put into service according to the EMD information provided in the following.

Portable and mobile RF communication equipment (e.g. cell phones) can affect Medical Electrical Equipment.

The use of accessories and cables other than those specified may result in increased emissions or decreased immunity of the unit.

Table 1 - EMISSION Limits -

Phenomenon	Compliance
Conducted and radiated RF EMISSION CISPR 11	Group 1, Class B
Harmonic distortion IEC 61000-3-2	Class A
Voltage fluctuations and flicker IEC 61000-3-3	Compliance

Table 2 - IMMUNITY TEST LEVELS : Enclosure Port -

Phenomenon	IMMUNITY TEST LEVELS
Electrostatic discharge IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Radiated RF EM fields IEC 61000-4-3	10 V/m 80 MHz - 2.7 GHz 80 % AM at 1 kHz
Proximity fields from RF wireless communications equipment IEC 61000-4-3	See table 4
Rated power frequency magnetic fields IEC 61000-4-8	30 A/m 50 Hz or 60 Hz

Table 3 - IMMUNITY TEST LEVELS : Input a.c. power Port -

Phenomenon	IMMUNITY TEST LEVELS
Electrical fast transients / bursts IEC 61000-4-4	±2 kV 100 kHz repetition frequency

Technical Data

Phenomenon	IMMUNITY TEST LEVELS
Surges Line-to-line IEC 61000-4-5	± 0.5 kV, ± 1 kV
Conducted disturbances induced by RF fields IEC 61000-4-6	3 V 0.15 MHz - 80 MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz
Voltage dips IEC 61000-4-11	0 % U_T ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°
	0 % U_T ; 1 cycle and 70 % U_T ; 25/30 cycle Single phase: at 0°
Voltage interruption IEC 61000-4-11	0% U_T ; 250/300 cycle

NOTE U_T is the AC mains voltage prior to application of the test level.

Table 4 - Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment -

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380 - 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 - 470	GMRS 460 FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
710	704 - 787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800 - 960	GSM 800/900 TETRA 800 iDEN 820 CDMA 850 LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						

Technical Data

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
1720	1700 - 1990	GSM 1800 CDMA 1900 GSM 1900 DECT LTE Band 1, 3, 4, 25 UMTS	Pulse modulation 217 Hz	2	0.3	28
1845						
1970						
2450	2400 - 2570	Bluetooth WLAN 802.11 b/g/n RFID 2450 LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
5240	5100 - 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5500						
5785						

Warranty

LIMITED WARRANTY

A&D Medical

For purchasers within the US only:

Product	Consumer Warranty Term
Monitor	5 year
Cuff	2 year

For outside of US, please contact local distributor or dealer.

Limited Warranty:

A&D Medical (“A&D”) warrants to the first purchaser (“You”) that the A&D product You purchased (the “Product”) will be free from defects in material, workmanship and design for the applicable Warranty Term stated above from the date You purchased the Product under normal use. This Limited Warranty is personal to You and is not transferable. If the Product is defective, then You return the Product to A&D in accordance with the procedure set forth below. A&D’s warranty obligation is limited to the repair or replacement, at A&D’s option, of the defective Product that has been returned by You within the warranty period. Such repair or replacement will be at no charge to You. The repaired or replacement Product is warranted here-under for the longer of the remainder of the original warranty period or 90 days from the date of shipment of the repaired or replacement Product.

To obtain a warranty service, please contact us in **US at 1-888-726-9966** or in **Canada at 1-800-461-0991** for return address, shipping and handling fee, and other instructions for processing warranty. Please ensure you have satisfactory proof of the date of Your purchase and a description of the defect. Returns will not be accepted unless a Return Material Authorization (RMA) Number has been issued from A&D Customer Service Representative.

This Limited Warranty does not cover, and A&D will not be liable for (i) any shipment damage, (ii) any dam-age or defect due to misuse, abuse, failure to use reasonable care, failure to follow written instructions enclosed with the Product, accident, subjecting the Product to any voltage other than the specified voltage, improper environmental conditions, or modification,

Warranty

alteration or repair by anyone other than A&D or persons authorized by A&D, or (iii) expendable or consumable components.

THIS LIMITED WARRANTY IS THE ONLY WARRANTY PROVIDED BY A&D; THERE ARE NO OTHER EXPRESS WARRANTIES. If A&D cannot reasonably repair or replace the Product, A&D will refund the amount You paid for the Product (not including taxes), less a reasonable charge for usage. To receive a refund you must have returned the Product and all associated materials to A&D. The above remedy of repair, replacement or refund is your only and exclusive remedy. IN NO EVENT SHALL A&D BE LIABLE FOR ANY DAMAGES, INCLUDING, WITHOUT LIMITATION, LOST PROFITS, LOST INFORMATION OR REPLACEMENT COSTS, ARISING OUT OF YOUR USE OF OR INABILITY TO USE THE PRODUCT, INCLUDING, WITHOUT LIMITATION, ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, EVEN IF A&D HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. Some states do not allow the exclusion of incidental or consequential damages, so that the above exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that may vary from state to state.

No distributor, dealer or other party is authorized to make any warranty on behalf of A&D or to modify this warranty, or to assume for A&D any liability with respect to its products.

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A&D Company, Limited
A&D Engineering, Inc.
A&D Electronics (Shenzhen) Co., Ltd.