

AUTOMATIC BLOOD PRESSURE MONITOR

TM-2657WP Series (TM-2657WP / TM-2657WP-RS / TM-2657WP-BT / TM-2657WP-BLE)

Instruction Manual ORIGINAL




Table of Contents

1.	Introduction, Warnings and Precautions	4
	1.1 Operational principles	4
	1.2 Warning definitions	4
	1.3 Precautions for use	5
	1.4 Precautions for safe measurement	8
	1.5 Unpacking	9
2.	Features	10
3.	Abbreviations and symbols	. 10
4.	Specifications	. 12
	4.1 Model configuration	12
	4.2 External dimensions	12
5.	Part Names	13
6.	About Blood Pressure	16
	Before use	
	7.1 Monitor installation	17
	7.2 Power connection	17
	7.3 Security slot	17
	7.4 Attaching the instruction panel	18
	7.5 Pre-inspection	18
8.		
•.		
0	Softing the clock	
	Setting the clock	
	Printer	21
	Printer 10.1 Installing the printer paper	.21 21
	Printer	21
10.	Printer 10.1 Installing the printer paper	21 21 22
10.	Printer	21 21 22
10.	Printer 10.1 Installing the printer paper 10.2 Selecting the print format Adding logo or other patterns	21 21 22 23
10. 11.	Printer 10.1 Installing the printer paper 10.2 Selecting the print format Adding logo or other patterns 11.1 Creating a bitmap pattern	21 21 22 23 23 24
10. 11.	Printer 10.1 Installing the printer paper 10.2 Selecting the print format Adding logo or other patterns 11.1 Creating a bitmap pattern 11.2 Sending bitmaps	21 21 22 23 23 24
10. 11.	Printer	21 22 23 23 24 24
10. 11.	Printer 10.1 Installing the printer paper 10.2 Selecting the print format Adding logo or other patterns 11.1 Creating a bitmap pattern 11.2 Sending bitmaps Changing functions 12.1 Procedure to Change Function Settings	21 22 23 23 24 24 24 24 27
10. 11.	Printer 10.1 Installing the printer paper 10.2 Selecting the print format Adding logo or other patterns 11.1 Creating a bitmap pattern 11.2 Sending bitmaps Changing functions 12.1 Procedure to Change Function Settings 12.2 Display time	21 22 23 23 24 24 24 24 27
10. 11.	Printer 10.1 Installing the printer paper 10.2 Selecting the print format Adding logo or other patterns 11.1 Creating a bitmap pattern 11.2 Sending bitmaps Changing functions 12.1 Procedure to Change Function Settings 12.2 Display time 12.3 Applied pressure	21 21 22 23 23 24 24 24 27 27 27 27
10. 11.	Printer 10.1 Installing the printer paper 10.2 Selecting the print format Adding logo or other patterns 11.1 Creating a bitmap pattern 11.2 Sending bitmaps Changing functions 12.1 Procedure to Change Function Settings 12.2 Display time 12.3 Applied pressure 12.4 IHB	21 21 22 23 23 24 24 24 27 27 27 27
10. 11.	Printer 10.1 Installing the printer paper 10.2 Selecting the print format Adding logo or other patterns 11.1 Creating a bitmap pattern 11.2 Sending bitmaps Changing functions 12.1 Procedure to Change Function Settings 12.2 Display time 12.3 Applied pressure 12.4 IHB 12.5 Print quality	21 21 22 23 23 24 24 24 27 27 27 27 27 27 27
10. 11.	Printer 10.1 Installing the printer paper 10.2 Selecting the print format Adding logo or other patterns 11.1 Creating a bitmap pattern 11.2 Sending bitmaps Changing functions 12.1 Procedure to Change Function Settings 12.2 Display time 12.3 Applied pressure 12.4 IHB 12.5 Print quality 12.6 ID and name printing	21 21 22 23 23 24 24 24 27 27 27 27 27 27 27
10. 11.	Printer 10.1 Installing the printer paper 10.2 Selecting the print format Adding logo or other patterns 11.1 Creating a bitmap pattern 11.2 Sending bitmaps Changing functions 12.1 Procedure to Change Function Settings 12.2 Display time 12.3 Applied pressure 12.4 IHB 12.5 Print quality 12.6 ID and name printing 12.7 Mean arterial pressure (MAP) printing	21 21 22 23 23 24 24 24 27 27 27 27 27 27 28
10. 11.	Printer 10.1 Installing the printer paper 10.2 Selecting the print format Adding logo or other patterns 11.1 Creating a bitmap pattern 11.2 Sending bitmaps Changing functions 12.1 Procedure to Change Function Settings 12.2 Display time 12.3 Applied pressure 12.4 IHB 12.5 Print quality 12.6 ID and name printing 12.7 Mean arterial pressure (MAP) printing 12.8 Measurement value printing	21 21 22 23 23 24 24 24 27 27 27 27 27 27 28 28
10. 11.	Printer 10.1 Installing the printer paper 10.2 Selecting the print format Adding logo or other patterns 11.1 Creating a bitmap pattern 11.2 Sending bitmaps Changing functions 12.1 Procedure to Change Function Settings 12.2 Display time 12.3 Applied pressure 12.4 IHB 12.5 Print quality 12.6 ID and name printing 12.7 Mean arterial pressure (MAP) printing 12.8 Measurement value printing 12.9 Graph printing	21 21 22 23 23 24 24 24 27 27 27 27 27 28 28 29
10. 11.	Printer 10.1 Installing the printer paper 10.2 Selecting the print format Adding logo or other patterns 11.1 Creating a bitmap pattern 11.2 Sending bitmaps Changing functions 12.1 Procedure to Change Function Settings 12.2 Display time 12.3 Applied pressure 12.4 IHB 12.5 Print quality 12.6 ID and name printing 12.7 Mean arterial pressure (MAP) printing 12.8 Measurement value printing 12.9 Graph printing 12.10 Comment printing	21 21 22 23 23 24 24 24 27 27 27 27 27 28 28 29 29
10. 11.	Printer 10.1 Installing the printer paper 10.2 Selecting the print format Adding logo or other patterns 11.1 Creating a bitmap pattern 11.2 Sending bitmaps Changing functions 12.1 Procedure to Change Function Settings 12.2 Display time 12.3 Applied pressure 12.4 IHB 12.5 Print quality 12.6 ID and name printing 12.7 Mean arterial pressure (MAP) printing 12.8 Measurement value printing 12.9 Graph printing 12.0 Comment printing 12.11 Bitmap/logo printing	21 21 22 23 23 24 24 24 27 27 27 27 27 27 27 28 8 8 29 29 29

Ар	pendix: EMD information	
17.	Warranty	54
	16.2 Standards	53
	16.1 Performance specifications	52
16.	Performance specifications and standards	
15.	Accessories and options list	
	14.8 Error codes	49
	14.7 Troubleshooting	49
	14.6 Disposing of the component parts	48
	14.5 Checking the number of measurements	47
	14.4 Replacing the arm cuff cover	45
	14.3 Periodic inspection	45
	14.2 Cleaning	43
	14.1 Inspection and safety management	42
14.	Maintenance	
	13.8 <i>Bluetooth</i> [®] Low Energy (External input/output unit : TM-2657-04)	40
	13.7 <i>Bluetooth</i> ® (External input/output unit : TM-2657-05)	38
	13.6 Bluetooth® (External input/output unit : TM-2657-04 & TM-2657-05)	37
	13.3 D-SUB 9-PIN MALE – (EXTERNAL INPUT/OUTPUT UNIT)	35
	13.2 Mini-DIN 8 pin female – (External input/output unit : only TM-2657-01)	34
	13.1 External input/output unit	34
13.	Transmission specifications	
	12.23 Bluetooth® Airplane mode	33
	12.22 Bluetooth® connection timing	33
	12.21 ICT printing	32
	12.20 Time format	32
	12.19 Date format	32
	12.18 Blood pressure result output	32
	12.17 Stop bit (D-Sub)	32
	12.16 Stop bit (Mini-DIN)	31
	12.15 Transmission speed (D-Sub)	31

1. Introduction

Indications for Use

- The TM-2657WP Series Automatic Blood Pressure Monitor is intended to be used by patients to measure systolic and diastolic blood pressure and pulse.
- The intended users are general adults, or 13 and older, with common knowledge about blood pressure measurement, who can perform a measurement on either their right or left arm.
- This device is designed to be used at outpatient clinics or general hospitals. It can also be used at health facilities, fitness gyms and other public facilities for blood pressure management of the visitors.

Notes

- Do not attempt to evaluate the blood pressure measurement results. Always consult with a doctor for evaluation of the
 results and treatment, especially when the results are greatly different from your ordinary values. Self-diagnosis and selftreatment from such results can be dangerous.
- Do not attempt to use this device on newborns or infants. Using this device on small children could cause injury to them. This device is designed for measuring persons 13 or older.
- Facilities with the device installed should employ at least one person who has good knowledge of blood pressure measurement and can give advice to users about how to pose for measurement or general information about blood pressure. The person should also have basic knowledge about maintenance of the monitor and know procedures to request training for maintenance if necessary.

1.1 OPERATIONAL PRINCIPLES

The cuff pressure is raised to approximately 30 mmHg higher than the anticipated systolic pressure and then gradually depressurized. Pulsations occur in the cuff pressure that matches the heart rate. These pulsations have an undulating pattern. They start small and then gradually increase with depressurization. After the maximum amplitude (MAP) is reached, they decrease. An oscillometric blood pressure monitor analyzes the amplitude waveform data of these pulsations to determine the systolic and diastolic blood pressures.

1.2 WARNING AND PRECAUTIONS

To prevent accidents due to inappropriate handling, this product and its manual contain the following warning signs and marks. The meaning of these warning signs and marks are as follows.

Warning Definitions

.	
	An imminently hazardous situation which, if not avoided, will result in death or serious injury.
	A potentially hazardous situation which, if not avoided, could result in death or serious injury.
	A potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practice.

Symbol Examples

Ý	The symbol \triangle indicates "Caution." The nature of the caution required is described inside or near the symbol, using text or a picture. The example on the left indicates caution against electrical shock.
\otimes	The symbol \otimes indicates "Do not." The prohibited action is described inside or near the symbol, using text or a picture. The example on the left indicates "Do not disassemble."
•	The symbol ● indicates mandatory action. The mandatory action is described inside or near the symbol, using text or a picture. The example on the left indicates general mandatory action.

Other

Precautions for each operation are described in the instruction manual. Read the instruction manual before using the device.

1.3 PRECAUTIONS FOR USE

In order to use the TM-2657WP Series Automatic Blood Pressure Monitor safely and correctly, carefully read the following precautions before using the monitor. The following content summarizes general matters regarding the safety of patients and operators, in addition to safe handling of the monitor.

When installing and storing the monitor

	Keep the monitor away from areas where flammable anesthetics or flammable gases are present, high- pressure oxygen chambers, and oxygen tents. Using the monitor in these areas may cause an explosion.

	Consider the following when using and storing the monitor. If the monitor is stored in an environment beyond the specified temperature or humidity, it may not perform to its capabilities.
	Avoid locations where the monitor may be splashed by water.
	• Avoid locations with high temperature, high humidity, direct sunlight, dust, salt and sulfur in the air.
	• Avoid locations where the monitor may be tilted, vibrated, or impacted (including during transportation).
	Avoid locations where chemicals are stored or gas occurs.
•	 Installation site: A location with a temperature between +10°C and +40°C (+50°F and +104°F) and a humidity between 15%RH and 85% RH (no condensation).
	 Storage site: A location with a temperature between -20°C and +60°C (-4°F and +140°F) and a humidity between 10% RH and 95% RH.
	• A location with an electrical outlet that can supply sufficient power (frequency, voltage, current) for the monitor.
	· Avoid locations where removal and insertion of AC power cable is prohibited.
	 The surface temperature of the cuff may become 46°C (115°F) when used in a 40°C (104°F) environment.

NOTE: Please be aware that the rubber feet may discolor the top of the stand.

Before using the monitor

A

Confirm for proper operation before use, if the packaging is damaged, unintentionally opened and exposed to environmental conditions outside of those specified.

	 Make sure that the electrical outlet is properly grounded and supplies the specified voltage and frequency (100-240V~ 50-60 Hz, more than 85VA). Connect the monitor to a grounded, 3-prong outlet. If a grounded, hospital-grade, 3-prong outlet is not available, connect the ground wire to an outlet with a contact terminal and ground it. Using the monitor with an incorrect outlet may cause an electrical shock.

Use the monitor safely and correctly.
 Connect all cables correctly and securely.
 Do not place objects on the monitor or power cable.
 Using other devices in conjunction with this monitor may cause incorrect diagnosis or safety problems. When used, check for safety.
 Always use accessories and consumables approved by A&D.
Carefully read the instruction manuals provided with optional items.
The precautions for these items are not listed in this manual.
• For safe and correct use of this monitor, always perform a pre-inspection (an inspection before use).
 If the monitor is covered with condensation, allow it to dry before switching the power on.
 If the monitor has not been used for an extended period, check that the monitor operates normally and safely before using it.
 The pressure of the cuff may cause a patient's arm to become numb.

1. Introduction

When using the monitor.

^

CONTRAINDICATIONS

\otimes	 Do not use a mobile phone near the monitor. It may cause a malfunction. Do not use the monitor in a moving vehicle as this may result in inaccurate measurements. 	

	Always check the condition of the monitor, its parts and the patient for safety.	
	• If a problem is found with the monitor, its parts or the patient, stop using the monitor, check the status of the patient and take appropriate actions.	
	• Frequent measurements can cause injury to the patient by interfering with blood flow.	
U	• Check the condition of the patient on a regular basis if measurements are performed frequently for a long time. There is a risk of causing damage by interfering with blood flow.	
	 To ensure accurate measuring, we recommend measuring blood pressure after being in a relaxed state for at least five minutes. 	

After using the monitor.

\otimes	Do not forcibly pull out the cables. Hold the connector with your hand when disconnecting the cables.	
0	 Use the specified procedure to return switches to their state before usage, then switch the power off. Clean the accessories and arrange them before storage. Keep the monitor clean and in proper operating condition so that it can be used without problem for the next operation. 	

If you suspect there is a problem with the monitor, perform the following actions.

	Ensure the safety of the patient.	
0	• Stop the operation of the monitor, switch the power off, and then disconnect the power cable from the outlet.	
	 If the air in the cuff is not released by pressing the Switch, press the FAST STOP switch. 	
	Label the monitor with a sign that says "Out of order" or "Do not use" and then contact A&D immediately.	

When performing a maintenance inspection.

	 For your safety, before performing a maintenance inspection, switch the power off and disconnect the power cable from the outlet. 	
•	 Always perform a pre-inspection and maintenance inspection to ensure safe and correct operation. The organization that installs the monitor (hospital, clinic, or other) is responsible for use, maintenance, and management of medical electrical devices. Neglecting pre-inspection and maintenance inspection can result in accidents. 	
\otimes	Never disassemble or modify the monitor (medical electrical device).	



When maintaining the monitor, use a dry, soft cloth. Do not use rags soaked in volatile liquids such as thinner and benzene.

0

Be aware that strong electromagnetic waves can cause malfunctions.

This monitor complies with EMD-standard IEC60601-1-2:2014+A1:2020. However, to prevent electromagnetic interference with other devices, do not use mobile phones close to within 30 cm of near the monitor.
 If this monitor is located near strong electromagnetic waves, noise may enter in waveforms and malfunctions may occur. If unexpected malfunctions occur during use of this monitor, inspect the electromagnetic environment and take appropriate actions.
The following are examples of general causes and countermeasures.
 Use of mobile phones Radio waves may cause unexpected malfunctions.
 Instruct visitors to rooms or buildings with medical electrical devices not to use mobile phones or small wireless devices.
 High frequency noise is being introduced from other devices via the electrical outlet.
 Check for the source of noise, and then perform countermeasures, such as using a noise cancellation device on this line.
 If the noise source is a device that can be stopped, stop using it.
- Use another electrical outlet.
• Effects from static electricity are suspected (discharges from devices or the surrounding area)
- Before using the monitor, ensure that the operator and patient have discharged static electricity.
 Humidify the room.
 If lightning is occurring nearby, the monitor may receive excessive voltage. In such cases, power the monitor using the following method.
 Use an uninterruptible power supply.

Environmental considerations

CAUTION Before disposing of this monitor, remove the lithium battery from the monitor.		
	 When any serious incident occurs in relation to this device, report to the manufacturer and the competent authority in your country. 	
0	 This device is not intended to diagnose heart arrhythmias. If the Irregular Heart Beat indicator illuminates frequently and is unrelated to patient movement during blood pressure measurement, seek medical attention. 	

FCC Caution

This device complies with Part 15 of FCC rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operations.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This transmitter must not be located or operated in conjunction with any other antenna or transmitter.

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment and meets the FCC radio frequency (RF) Exposure Guidelines in Supplement C to OET65. This equipment has very low levels of RF energy that it deemed to comply without maximum permissive exposure evaluation (MPE). But it is desirable that it should be installed and operated keeping the radiator at least 20 cm or more away from person's body (excluding extremities: hands, wrists, feet and ankles).

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- · Reorient or relocate the receiving antenna.
- · Increase the separation between the equipment and receiver.
- · Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

1. Introduction

1 4 PRECAUTIONS FOR SAFE MEASUREMENT


The following lists precautions related to measurement. Always consult with a doctor for evaluation of the results and treatment. Self-diagnosis and self-treatment from results can be dangerous.

0	 Do not measure on an arm receiving an intravenous drip or blood transfusion. This may cause an accident. Do not perform measurement if the arm has external injuries. Not only will the wound worsen, there is a risk of spreading disease. If the arm cuff cover is soiled with blood, dispose of the cover. There is a risk of spreading disease. 	



Items that may be contaminated must be disposed of as medical waste.

A CAUTION



Measurement is intended for arms with circumferences of

18 to 42 cm (7.1 to 16.5").

• The arm of the patient is wet. Wet arms may cause accidents or electrical shock.

NOTES

- Blood pressure measurement may cause subcutaneous bleeding. This subcutaneous bleeding is temporary and disappears with time.
- Consult a doctor before use if you have had a mastectomy or lymph node resection.
- · If thick clothing is worn, correct measurement is not possible. Measure when the patient is wearing a sleeveless or thin shirt
- If the patient rolls up their sleeve and the sleeve pinches their arm, correct measurement is not possible.
- Measurement is not possible with patients with peripheral hypoperfusion, very low blood pressure, or low body temperature (since blood flow to the measurement location is low).
- · Measurement is not possible for the patient who has frequent arrhythmia recurrences.
- · Measurement locations are restricted to the right and left upper arms. Other locations cannot be measured.
- Insert the arm into the arm insertion section up to the top of the shoulder.
- Adjust the height of the cuff so that is the same height of heart using chair. If height of cuff and patient' heart are different, correct measurement is not possible. If the patient does not feel well, stop measurement immediately and take appropriate actions.
- · Hold the optional air cushion chair firmly when you sit on it.
- · Measurement cannot be performed with the following patients.
 - Patients who have just exercised: Blood pressure after exercise is higher than ordinary blood pressure.
 - Measure after the patient has rested for several minutes and has taken deep breaths.
 - Patients with shaking arms: If the patient's body moves, correct measurement is not possible. Wait until the shaking stops, and then perform measurement. (This includes shaking from the cold or muscle movements after moving heavy objects.)
- · Consult the doctor for any of the following situations.
 - The application of the cuff on any limb with intravascular access or therapy, or an arteriovenous (A-V) shunt.
 - Simultaneous use with other medical monitoring equipment on the same limb.
 - The blood circulation of the patient needs to be checked.


1.5 UNPACKING



This monitor is a precision device and must be handled carefully. If it receives a strong impact, it may be damaged.

NOTE: This monitor has been shipped in specially designed packaging to prevent damage during shipping. Check the monitor for damage when unpacking it.

Before using the monitor, check the main unit and each standard accessory for damage. For optional items, see "15. ACCESSORIES AND OPTIONS LIST".



Confirm that all of the parts are included to ensure that the medical device is ready for safe use and as intended.

2. Features

- Measurement can be performed using either the right or left arm.
- The arm cuff is inflated around the arm by pressing the button and deflation speed is automatically controlled. No special adjustment is required. All you have to do is insert your arm into the arm insertion section to the shoulder and press the button. The rest of the procedure is done automatically for a quick and easy measurement of blood pressure.
- The printer is equipped with a cutter to automatically cut the printed paper.

Options

An optional external input/output unit can be connected to a computer for data management or automation as necessary.

3. Abbreviations and Symbols

Abbreviation/Symbol	Meaning
~	Alternating current
mmHg	Blood pressure unit (millimeters of mercury)
/min.	Heartbeats per minute
	Displayed when measurement is not possible
SYS	Systolic blood pressure
MAP	Mean arterial blood pressure
DIA	Diastolic blood pressure
PUL	Pulse
((\C)))	Irregular Heartbeat symbol (IHB)
\bigcirc	Time
0	Power off (disconnected from the power source)
	Power on (connected to the power source)
\Diamond	START button
\bigcirc	STOP button
LOT	Batch code
REF	Catalog number
SN	Serial number
Exx	Error code display (xx=00 to 99)
<u>^</u>	Displays extent of electric shock protection: B-type applied part
(Follow Instructions for use
20XX	Date of manufacture
G→	RS-232C serial interface
<u>×</u>	WEEE label
	Manufacturer
(MEASUREMENT IN PROGRESS)	Displays the measurement status. "MEASUREMENT IN PROGRESS".
TAKE MEASUREMENT AGAIN	Displays the measurement status. "TAKE MEASUREMENT AGAIN"
FAST STOP	"FAST STOP" for rebooting the device.
Please do not pull print paper during printing.	Caution: "Please do not pull printer paper during printing."
(The printer paper is automatically cut.)	Caution: "The printer paper is automatically cut."
POWER	"POWER" switch.


Abbreviation/ Symbol	Meaning
SELECT	Used to change functions.
	Used to change function setting.
COUNT	Used to display the number of measurements to date.
Paper	Describes how to change printer paper.
((*))	To indicate generally elevated, potentially hazardous, levels of non-ionizing radiation, or to indicate equipment or systems e.g. in the medical electrical area that include RF transmitters or that intentionally apply RF electromagnetic energy for diagnosis or treatment.
MD	Medical device
IP	International protection symbol
UDI	Unique device identifier
Ĩ	Refer to the instruction manual.
X	Temperature limitation
<u>s</u>	Humidity limitation
6.9	Atmospheric pressure limitation
<u><u></u><u></u><u></u><u></u></u>	This way up
ų	Fragile, handle with care
	Stacking limit by 3
¢	Universal recycling symbol
\$	Handle with care
#	Keep dry

4. Specifications

4.1 MODEL CONFIGURATION


Printer	0
Measurement status LED	0
Time, Date format	12 hour, MM/DD/YYYY

4.2 EXTERNAL DIMENSIONS




Unit: mm

Front



No.	Name	Description
1	Printer cover button	Opens the printer cover.
2	Printer paper opening	Opening where printer paper comes out.
3	Printer cover	Holds down the printer paper.
4	Power cable	AC power cable.
5	Systolic blood pressure display	Displays the systolic blood pressure measurement value. When a measurement error occurs, the error code is displayed.
6	Diastolic blood pressure display	Displays the diastolic blood pressure measurement value. Displays the pressure during measurement.
7	Pulse display	Displays the pulse measurement value.
8	Clock display	Displays the current time. (12 or 24 hour)
		Displays the measurement status.
9	Measurement status LED	"MEASUREMENT IN PROGRESS"
		"TAKE MEASUREMENT AGAIN"
10	button	If this button is pressed in standby mode, blood pressure measurement is started.
10		If this button is pressed during blood pressure measurement, blood pressure measurement is stopped.
11	Arm cuff cover	Inner cover of the cuff.
12	Front cuff retainer	Holds the arm cuff cover.
13	FAST STOP button	If this button is pressed, the power is switched off and measurement is stopped.


Rear



No.	Name	Description
1	SELECT button	Used to change functions.
2	▲ button	If pressed when the number of measurements to date is displayed, the number of measurements is printed. Used to change functions.
3	COUNT button	Displays the number of measurements to date. (See "12.5. Checking the number of measurements")
4	Bitmap SD socket cover	Use only for maintenance.
5	Arm cuff cover	Inner cover of the cuff.
6	Rear cuff retainer	Holds the arm cuff cover.
7	Armrest	Location to rest the arm during measurement.
8	External input/output unit	The optional external input/output unit.
9	POWER switch	Switches the power on and off. Once the power is switched on, the monitor will be in the standby mode.
10	Cover for pressure inspection area	Used to check pressure accuracy.
11	AC INPUT connector	Location to insert the power cable.
12	Security slot	Can be used with a security cable to secure the monitor. (For theft prevention)


External input/output units (optional)

TM-2657-01 External input/output unit RS 2ch




No.	Name	Description
1	Mini-DIN 8 pin female	RS-232C
2	D-Sub 9 pin male	RS-232C

TM-2657-04 External input/output unit RS+Bluetooth® Low Energy



No.	Name	Description
-	Bluetooth [®] low energy	Bluetooth® Ver. 4.2
1	D-Sub 9pin Male 🕬 😳 🕬	RS-232C

TM-2657-05 External input/output unit RS+Bluetooth®



No.	Name	Description
-	Bluetooth®	<i>Bluetooth</i> [®] Ver.2.1 class1 SPP HDP correspondence
1	D-Sub 9pin Male 🕬 😳 🕬	RS-232C


Note: For details on EXTERNAL INPUT/OUTPUT UNIT (TM-2657-01, TM-2657-04, TM-2657-05), contact your local A&D dealer.

6. About Blood Pressure

Blood pressure variations

Blood pressure is highly sensitive and changes subtly with each beat to match the condition of the heart. It may vary by 30 to 50 mmHg in response to various conditions.

That's why it's important not to focus on a single measurement, but instead measure every day at the same time to learn your average blood pressure and blood pressure trends. This blood pressure information will be important when visiting a doctor. Consult with a doctor to determine the meaning of your results.



What types of high blood pressure are there?

There are 2 types of high blood pressure: essential hypertension and secondary hypertension. Secondary hypertension is caused by disease that raises blood pressure. When kidney inflammation or pregnancy toxicosis causes high blood pressure, treat the problem and the blood pressure will fall naturally.

In the case of essential hypertension, the cause is not clear, but the blood pressure is high. The combination of long periods of stress, high salt intake, obesity and genetic problems can cause essential high blood pressure. Of these causes, genetics play a large factor. If both or one parent has high blood pressure, the occurrence rate of high blood pressure is 60% and 30%, respectively, indicating a genetic component.

What is IHB (Irregular Heartbeat)?

IHB appears when an irregular heartbeat is detected. The mark is printed when a very slight vibration like shivering or shaking is detected.

When the monitor detects an irregular rhythm during the measurements, the IHB symbol will appear on the display with the measurement values.

Ensure your arm is positioned correctly, then take another measurement.

If the symbol continues to appear, we recommend you to consult with your physician.

NOTE We recommend contacting your physician if you see this ((()) IHB symbol frequently.

When is the IHB mark printed?

The IHB mark is printed in the measurement data in the following two cases.

- · When a beat varies during measurement.
- · When the arm or monitor is moved during measurement.

See the precautions at the beginning of this manual and install the monitor in an appropriate location using a safe and correct method.

7.1 MONITOR INSTALLATION

Attaching the armrest

Place the monitor on a stand so that measurement can be performed in an appropriate posture. The patient's heart and the cuff should be at the same height and the patient should be relaxed.

To prevent theft, we recommend using a chain to connect the security slot and stand. (See "6.3. Security slot")




7.2 POWER CONNECTION



To avoid the risk of electric shock, the monitor must be used with an electrical outlet that is properly grounded and supplies the specified voltage and frequency (100-240V~ 50-60 Hz, more than 85VA).
Connect the monitor to a grounded, 3-prong outlet.

Use the 3-prong power cable provided with the monitor to connect between the AC INPUT connector and an electrical outlet.



7.3 SECURITY SLOT

The monitor can be secured to a table or other secure location by passing a security cable through the hole of the protruding tab on the monitor to secure it.

7. Before Use


7.4 ATTACHING THE INSTRUCTION PANEL


See the illustration below to attach the instruction panel to the rear side of the monitor.



Make sure to attach the instruction panel to the main unit before use. The instruction panel contains precautions that the patient must observe to use the monitor safely and correctly.

A CAUTION





7.5 PRE-INSPECTION



Perform the pre-inspection every day to ensure safe and correct usage.

Introduction

Before using the monitor for the first time each day, perform the following pre-inspection.

Before switching the power on

- · Is there any external deformation or damage to the monitor?
- · Is the monitor wet?
- · Is the monitor in a stable location free of tilting, vibrations and impacts?
- Is there damage or abnormalities around where the arm is inserted (cuff area)?
- · Is the arm cuff cover attached?
- · Is the arm cuff cover overstretched?

Connection cable

Are the optional cables inserted firmly into the connectors of the monitor?

Power cable

Make sure that the electrical outlet is properly grounded and supplies the specified voltage and frequency (100 V-240 V \sim 50 Hz-60 Hz).

After switching the power on

If you observe smoke or a strange smell or noise - discontinue use and call A&D Customer Service.


Checking the time

Is the time set correctly? If the time is incorrect when recording data, the data will be incorrect.

Checking the display

After switching the power on, all LEDs switch on for several seconds. When the diastolic blood pressure display shows "0", the monitor is ready for blood pressure measurement.





When performing continuous measurements, wait 2 to 3 minutes between measurements for the patient to relax.

· To obtain accurate measurement results, ensure the patient sits with good posture and his/her back straight, and with his/

Blood pressure measurement results are affected by the posture and physical condition of the patient.

- 6. The measurement results are displayed.
- 7. The measurement results are printed on the printer paper. Remove the arm from the cuff.

her feet flat on the floor without crossing legs. Ensure the patient is relaxed and remains still.
Adjust the height of the chair such that the cuff is at the same height as the heart. If the cuff is not at the same height as

the heart, correct measurement is not possible.

Notes

blayed. 75 ted on the printer paper. Remove

9. Setting the clock


To set the date and time, use the select button on the back of the monitor.


Setting the date and time:

- While the monitor is in standby mode, press the SELECT button for 1 second until the Systolic blood pressure display starts to flash with the year value Select the correct year by pressing the ▲ button.
- 2. After reaching the desired year, press the **SELECT** button again.
- 3. Repeat the process for month, day, hour and minute.
- Once the desired date and time is selected, press the button to save the changes and return to standby mode.
 Note: If the COUNT button is pressed while configuring settings, changes are not saved and the monitor returns to standby mode.

Notes

- If no operation is performed for about 10 seconds, the specified settings are set. After Adult is displayed for 2 seconds, the monitor returns to the standby mode.
- Dates up to December 31, 2050 are supported.






10.1 INSTALLING THE PRINTER PAPER





2. Install the printer paper in the way shown in the illustration below.



3. With the end of the paper at the top and protruding out, secure the printer paper by closing the printer cover until you hear a click. If the cover is not completely closed, a paper jam may occur.



- If the high-speed printing mode is used, approximately 700 prints are possible from one printer paper roll. With 3-line printing mode, 600 prints are possible. When the end of the printer paper roll becomes pink, replace the paper.
- Use thermal paper only.
- If the following error codes are displayed in the systolic display section, a printer error has occurred. Perform the required countermeasure.

Error code	Error/countermeasure			
PE	No printer paper. Install a new printer paper roll.			
Po	The printer cover is open. Firmly close the printer cover.			
Pc	A printer cutter error. Open the printer cover, check the printer paper, and then close the printer cover.			

• When no printer error is displayed and the monitor is in standby mode, holding down the ▲ button for 2 seconds will cut the paper.

	 There will be no printing if the direction of the printer paper is incorrect. Use genuine A&D printer paper. If genuine A&D paper is not used, the print may be too light or paper jams may occur.
	• On the last 60 cm (24") of printer paper, there are pink end marks (pink lines on both sides). If these end marks appear, replace the printer paper.
Notes	Thermal printer paper is used. Note that discoloration or fading may occur.
Notes	 Items that will be discolored: Felt-tip pens and adhesive agents including starch and organic solvents. Items that can cause fading: Highlight pens, tape, transparent storage cases, desk pads, sunlight and ultraviolet. To ensure you don't lose your measurement due to discoloration or fading, make a copy of measurement results.

10. Printer

10.2 SELECTING THE PRINT FORMAT

Users can format the information on the printout. The printing area is divided into 4 sections: print header, measurement value, graph and bitmap. Each section has printing items available for selection. See below for the options and see "10. CHANGING FUNCTIONS" for more details on how to change the settings.

Print Header

Function Setting	Function Detail	Options		
		oFF	ID and name are not printed	
ENR		1	Only name printed	
- uo	ID and name printing	0	Only ID printed	
		П	Both name and ID printed	
	IHB printed	on	Will print (C)) if there is an IHB detected during measurement.	
FOS		o88	Will not print ${}^{(\!$	
		orr	measurement.	
626	Measurement start	ΕIJ	DD, month, YYYY month, DD, YYY	
rco	date format	US	Month, DD, YYY	
620	Measurement start	25	24 hour clock	
	time format	12	12 hour clock	
	Height and Weight format	off	Height and weight are not printed	
F 16		1	Printer mode printing	
		2	Integrated mode printing	

Measurement Value

Function Setting	Function Detail		Options
		1	High-speed printing
E ; ;	Measurement Value Printing	2	Normal 3-line printing
1 1		3	Big font printing
		4	Table printing
FD9	Mean arterial pressure	on	MAP not printed
-up	(MAP) printing	oFF	MAP printed

Graph

Function Setting	Function Detail		Options
C (D	One she aristic a	oFF	Graph not printed
	Graph printing	1	Pulse fluctuation graph printed

Bitmap (Logo)

Function Setting	Function Detail	Options	
	Bitmap printing	_o ff	oNo printing
E IS		1	Standard pattern printing
		2	User pattern printing

ICT

Function Setting	Function Detail	Options	
	ICT Printing	oFF	No printing
		1	Barcode printing (CODE39)
F29		2	QR code printing with ID
		Э	Barcode printing (CODE39 with check digit (modulus43))
		Ч	QR code printing V2 with ID

Print example 1: Initial settings				
Name		(()))		
17 Oct.,	2015	22:18	F26: Date format [1] (EU format)	
SYS	13	U mmHg	F27: Time format [24] (24 hour)	
DIA	9	6 mm+g ⁻	F11: Measurement value printing [2] (Normal 3-line printing)	
PUL		/min.	(Romano into printing)	


Printing example 2

ID: 1234567890123456 Name	F08: ID printing [3] F05: IHB [on] (no IHB detected)
17 Oct., 2015 22:18 SYS DIA PUL 130 06 71	F26: Date format [1] (EU format) F27: Time format [24] (24 hour)
mmHg mmHg /min. MAP	F11: Measurement value printing [1] (High-speed printing)
	F09: MAP printing [on]

Printing example 3

Name			F05: IHB [on] (no IHB detected)
Nov 5,	2015	3:37 PN	F26: Date format [1] (US format)
SYS	DIA	ΡUΙ	F27: Time format [12] (12 hour)
130	96	71	F11: Measurement value printing [1]
nmHį	g m	mHg /mir	(High-speed printing)
			F09: MAP printing [off]

11. Adding Logo or Other Patterns



	•	The bitmap must be in black and white.
Notes	•	For operable SD card standard, the device operation is checked with SD and SDHC. Some SD cards cannot be recognized with the device, if your card cannot be recognized, please use other SD card.
	•	For a file system, the device operation is checked with FAT16 and FAT32.


11. Adding Logo or Other Patterns

11.2 SENDING BITMAPS


- 1. Switch off the power of the monitor.



2. With the COUNT, ▲ and SELECT buttons pressed, switch the power on. The monitor enters the bitmap transfer mode.



3. Insert the SD card containing the bitmap file (Logo.bmp) saved in "15.1. Size of original bitmap patterns" into the SD socket. Press the socket. Press th




After transfer, restart the power, and then set the function **F15** to **2**. The bitmap is printed with the blood pressure value after blood pressure measurement.

12. Changing Functions

The multi-functional monitor can be configured for various applications by changing function settings. To change function settings, use the buttons located on the rear panel of the monitor while the monitor is in standby mode.

12.1 PROCEDURE TO CHANGE FUNCTION SETTINGS

- In power off mode, hold both the ▲ and SELECT buttons down and switch the power on.
 F01 is displayed in the systolic display section and the monitor enters the function changing mode.
- 2. Each time the SELECT button is pressed, the setting item changes to F02, F03...
- 3. Each item can be changed using the ▲ button.
- 4. After completing the settings, switch the power off and then on again.



All available functions:

Setting Items	Details	Default	Diastolic Display Section	Function		
F0 (Not used	_				
F02	Display time	20	oFF, S, 10, 20, 999	Measurement result display time (seconds)		
F03	Applied pressure	Rut	Rut, 160, 180, 200	Applied pressure setting (mmHg)		
FOH	Not used	_				
F05	IHB	on	oFF, on	IHB-mark printing on/off		
F05	Not used	—				
			0FF	Printing off		
FON	Print quality/light or		1	Light printing (high speed)		
	dark	0	2	Standard printing		
			3	Dark high-quality printing (low speed)		
			0FF	ID : No Name: No		
coo	ID and name printing		1	ID : No Name: Yes		
F08	ID and name printing	0	2	ID : Yes Name: No		
			Э	ID : Yes Name: Yes		
F09	Mean arterial pressure (MAP) printing	on	oFF, on	Mean arterial blood pressure (MAP) printing on/off		
F 10	Not used	—				
			1	High-speed printing		
F	Measurement value	0	2	Normal 3-line printing		
, , ,	printing		З	Big font printing		
			Ч	Table printing		
		0	0FF	Graph printing off		
F 12	Graph printing		1	Pulse fluctuation graph printing		
			2	AHA graph printing		
F 13	Comment printing	0	0FF	Comment printing off		
1 1	Comment printing		on	Blood pressure classification comment printing		
F 14	Not used	—				
		0	oFF	Bitmap printing off		
F 15	Bitmap printing		1	Standard pattern printing		
			2	User pattern printing		
	Lisisht and weight		0FF	Height and weight values printing OFF		
F 15	Height and weight values printing		1	Printer mode printing		
	values printing	0	2	Integrated mode printing		
F 17	Not used	—				
F 18	Beep sound	on	oFF, on	Beep sound on/off		
F 19	Not used	—				
			0FF	No connection		
	External input/output	0	1	Mini-DIN: Blood pressure result input/output (STD/RI/ RB/BP/RA) D-Sub: Blood pressure result input/output (STD/RI/ RB/BP/RA)		
F20	External input/output protocol		2	Mini-DIN: A&D weight scale D-Sub:Blood pressure result input/output (STD/RI/RB BP/RA)		
			Э	Mini-DIN: Blood pressure result input/output (STD/RI/ RB/BP/RA) D-Sub: ID reader		

12. Changing Functions

Setting items	Details	Default	Diastolic display section	Function
			Ч	Mini-DIN: Blood pressure result input/output (STD/RI/RB/BP/RA) D-Sub: Ux compatibility
6 20	External input/output		5	Mini-DIN: Blood pressure result input/output (STD/RI/RB/BP/RA) D-Sub: RVX compatibility
F20	protocol		6	Mini-DIN: Blood pressure result input/output (STD/RI/RB/BP/RA) D-Sub: A&D weight scale
			η	Mini-DIN: Blood pressure result input/output (STD/RI/RB/BP/RA) D-Sub: RVY compatibility
			120	1200 bps
F2 I	Transmission speed	0	240	2400 bps
	(Mini-DIN		480	4800 bps
			960	9600 bps
			120	1200 bps
C 7 7	Transmission speed	0	240	2400 bps
F22	(D-Sub)		480	4800 bps
			960	9600 bps
		0	1	Stop bit: 1
F23	Stop bit (Mini-DIN)		2	Stop bit: 2
		0	1	Stop bit: 1
F24	Stop bit (D-Sub)		2	Stop bit: 2
		0	1	RB (no ID, immediately after measurement) + STD
			2	RI (with ID, immediately after measurement) + STD
F25	Blood pressure result output		З	BP (with ID, immediately after measurement) only
	output		Ч	STD (command response) only
			5	RA (with ID, immediately after measurement)
F26	Date format	US	EU	DD month., YYYY
1 6 0	Date Ionnat	00	US	month. DD, YYYY
F27	Time format	12	24	24 hour
			12	12 hour (AM/PM)
F28	Not used	_		
		0	0FF	ICT printing OFF
			1	Bar code printing (CODE39)
F29	ICT printing		2	QR code printing, including ID
			Э	Bar code printing (CODE39 , with check digit (modulus43))
			Ч	QR code printing V2, including ID
F3 I	<i>Bluetooth</i> [®] connection timing	0	1	Connection at the end of measurement
	uning		2	Connection at the start of measurement
F3S	Airplane mode	0	1	Airplane mode OFF
'			2	Airplane mode ON

	 F16 setting is valid only if F20 setting is 2 or 6. Initial values are determined depending on each destination.
Notes	 F35 setting is valid only when TM2657-04 is installed. To reset all settings to factory default settings, hold the button for 5 seconds when any of the "FXX" numbers are displayed.

Note: The following options are changes by using the ▲ button to change the setting. The setting will appear in the diastolic display section.

12.2 DISPLAY TIME

					-	
F02	Time	the	display	s showina.	after	measurement results.

DIA LED	Display time setting	Default
oFF	No display of results (All values are displayed as "")	
5	5 seconds]
10	10 seconds	20
- 20	20 seconds	
999	Remains displayed]

12.3 APPLIED PRESSURE

F03: Applied pressure setting. Note: If automatic applied pressure (Aut) is set, pulsation is observed while pressure is applied and the applied pressure value is automatically determined.

DIA LED	Applied pressure setting	Default
Rut	Automatic applied pressure	
150	160 mmHg	Rut
180	180 mmHg	1 HUE
200	200 mmHg	

12.4 IHB

F05: IHB setting.

DIA LED	IHB setting	Default
0FF	IHB off	
on	IHB on	00

When IHB is on - Printing example

When IHB is detected				When IHE	s is not de	tected
Name		«Q»		Name		
17 Oct.,	2015 2	22:18	IHB	17 Oct.,	2015	22:18

12.5 PRINT QUALITY

F07: Print quality setting.

DIA LED	Print quality setting	Default
0FF	Printing off	
1	Light printing (high speed)	
2	Standard printing	
3	Dark high-quality printing (low speed)	


12.6 ID AND NAME PRINTING

To input an ID, set the function **F20** to **3**, and connect an ID reader. The ID data is maintained until the blood pressure is measured correctly and is cleared immediately after the result is displayed or printed.

F08: ID printing setting.

DIA LED	ID printing setting	Default
0FF	ID : No / Name : No	
1	ID : No / Name : Yes	
2	ID : Yes / Name : No] '
3	ID : Yes / Name : Yes	

When ID and name printing is on - Printing example




12. Changing Functions

12.7 MEAN ARTERIAL PRESSURE (MAP) PRINTING

F09: Mean arterial pressure (MAP) printing setting.

DIA LED	Mean arterial pressure printing	Default
oFF	Mean arterial pressure (MAP) printing off	oFF
on	Mean arterial pressure (MAP) printing on	000

When mean arterial pressure (MAP) printing is on - Printing examples




12.8 MEASUREMENT VALUE PRINTING

F11: Measurement value printing setting.

DIA LED	Measurement value printing mode	Default
1	High-speed printing	
2	Normal 3-line printing	-
3	Big font printing	_ C
Ч	Table printing	

When Mean arterial pressure (MAP) printing is off - Printing examples





Note: In the table printing mode, paper is not cut automatically. To cut paper, hold the ▲ button for 2 seconds while the monitor is in the standby mode.


12.9 GRAPH PRINTING

F12: Graph printing settings.

DIA LED	Graph printing	Default
0FF	Graph printing off	
1	Pulse fluctuation graph printing	oFF
2	AHA graph printing	







12.10 COMMENT PRINTING

F13: Blood pressure classification comment printing settings (pre-programmed comments based on your blood pressure value).

DIA LED	Comment Printing	Default
oFF	Comment printing off	OFF
00	Blood pressure classification comment printing	

Printing example - Blood pressure classification comment printing



12.11 BITMAP/LOGO PRINTING

F15: Bitmap printing setting.

DIA LED	Bitmap printing	Default
0FF	Bitmap printing off	
1	Standard pattern printing	oFF
2	User pattern printing	

For details about bitmap registration and user pattern (e.g. logo) printing, see "15. SENDING BITMAP PATTERNS".

Bitmaps up to 384 x 640 pixels can be printed.

Printing example – Standard pattern printing



12. Changing Functions

12.12 BEEP SOUND

F18: Sound setting. When a measurement starts/ends there will be an audible sound (beep).

DIA LED	Buzzer	Default
oFF	Beep sound off	
00	Beep sound on	00

12.13 EXTERNAL INPUT/OUTPUT PROTOCOL

F30: Protocol settings for connections.

External input/output unit <TM-2657-01>

DIA LED		External input/output unit (option) protocol	Default
oFF	No connectio	n	
1	Mini-DIN: O	Blood pressure result input/output (STD/RI/RB/BP/RA)	
1	D-Sub: 📼	Blood pressure result output (STD/RI/RB/BP/RA)	
2	Mini-DIN: O	A&D height and weight scale	
	D-Sub:	Blood pressure result input/output (STD/RI/RB/BP/RA)	
3	Mini-DIN: O	Blood pressure result input/output (STD/RI/RB/BP/RA)	
2	D-Sub:	ID reader	
ч	Mini-DIN: O	Blood pressure result input/output (STD/RI/RB/BP/RA)	1
	D-Sub: 📼	Ux compatibility	
5	Mini-DIN: O	Blood pressure result input/output (STD/RI/RB/BP/RA)	
	D-Sub: 📼	RVX compatibility	
6	Mini-DIN: O	Blood pressure result input/output (STD/RI/RB/BP/RA)	
0	D-Sub: 📼	A&D weight scale	
7	Mini-DIN: O	Blood pressure result input/output (STD/RI/RB/BP/RA)	
1	D-Sub: 📼	RVY compatibility	

External input/output unit <TM-2657-04>

DIA LED	External input/output unit (option) protocol	Default
oFF	No connection	
	Bluetooth® Low Energy: Blood pressure result input/output (RX)	
1	D-Sub: Blood pressure result output (STD/RI/RB/BP/RA)	
2	Bluetooth® Low Energy: Blood pressure result input/output (RX)	
C	D-Sub: Blood pressure result input/output (STD/RI/RB/BP/RA)	
3	Bluetooth® Low Energy: Blood pressure result input/output (RX)	
1	D-Sub: ID reader	
ч	Bluetooth® Low Energy: Blood pressure result input/output (RX)	1
1	D-Sub: I Ux compatibility	
-	Bluetooth® Low Energy: Blood pressure result input/output (RX)	
5	D-Sub: EVX compatibility	
6	Bluetooth® Low Energy: Blood pressure result input/output (RX)	
0	D-Sub: Come A&D weight scale	
7	Bluetooth® Low Energy: Blood pressure result input/output (RX)	
1	D-Sub: E RVY compatibility	

External input/output unit <TM-2657-05>

DIA LED	External input/output unit (option) protocol	Default
oFF	No connection	
1	D-Sub: Blood pressure result output (STD/RI/RB/BP/RA)	
2	D-Sub: Blood pressure result input/output (STD/RI/RB/BP/RA)	
3	D-Sub: ID reader	
닉	D-Sub: E Ux compatibility	
5	D-Sub: E	
6	D-Sub: E A&D weight scale	
7	D-Sub: E RVY compatibility	

For details on communication commands (STD/RI/RB/BP/RA), contact your local A&D dealer. For details on connecting ID readers, weight scales, or computers, contact your local A&D dealer.

12.14 TRANSMISSION SPEED (MINI-DIN)

F21: Mini-DIN ⁽ⁱⁱⁱ⁾ transmission speed setting.

DIA LED	Bitmap printing	Default
150	1200 bps	
240	2400 bps	240
480	4800 bps	0
960	9600 bps	

12.15 TRANSMISSION SPEED (D-SUB)

F22: D-Sub event transmission speed setting.

DIA LED	Transmission speed (D-Sub)	Default
150	1200 bps	
240	2400 bps	240
480	4800 bps	טרס ן
960	9600 bps	

12.16 STOP BIT (MINI-DIN)

F23: Stop bit (Mini-DIN ())setting.

DIA LED	Stop bit (Mini-DIN)	Default
1	Stop bit 1	1
2	Stop bit 2	i

12.17 STOP BIT (D-SUB)

F24: Stop bit (D-Sub eme)setting.

DIA LED	Stop bit (D-Sub)	Default
1	Stop bit 1	
2	Stop bit 2	1

12.18 BLOOD PRESSURE RESULT OUTPUT

F25: Blood pressure result output setting.

DIA LED	Blood pressure result output	
1	RB (no ID, immediately after measurement) + STD	
2	RI (with ID, immediately after measurement) + STD	
3	BP (with ID, immediately after measurement) only	
Ч	STD (command response) only	
5	RA (with ID, immediately after measurement)	1

12.19 DATE FORMAT

F26: Printing date format setting.

DIA LED	Date Format	Default
EU	DD month., YYYY	1.05
US	month DD, YYYY	00

12.20 TIME FORMAT

F27: Time format setting.

DIA LED	Time Format	Default
24	24 hour	
12	12 hour (AM/PM)	05

12.21 ICT PRINTING

F29: ICT printing setting.

DIA LED	ICT Printing		Default
oFF	ICT printing	OFF	
1	Bar code printing	CODE39	
2	QR code printing	including ID	0FF
3	Bar code printing	CODE39, with check digit (modulus43)	
닉	QR code printing V2	including ID	

The following information is included in code printing.

Bar code printing (CODE 39) Systolic blood pressure value, Mean blood pressure value, Diastolic blood pressure value, Pulse rate	Name Jan 29, 2015 13:37 SYS 119 mmHg DIA 68 mmHg
QR code printing including ID	PUL (min. +119094068080*
YYYY/MM/DD/HH/MM, ID (16 digits), Systolic blood pressure value, Mean blood pressure value, Diastolic blood pressure value, Pulse rate	29 Jan., 2015 11:44 sys 115 mmHg DIA 689 mmHg PUL 9/min.
Bar code printing (CODE39), with check digit (modulus43) Systolic blood pressure value, Diastolic blood pressure value, Pulse rate	Name Jan 29, 2015 13:32 SYS 118 mmHg DIA 67 mmHg PUL 80 /min.
QR code printing V2 including ID YYYY/MM/DD/HH/MM, ID (16digits), Systolic blood pressure value, Mean blood pressure value, Diastolic blood pressure value, Pulse rate, Height value, Weight value	Name Jan 29, 2015 13:33 SYS 1200 mmHg DIA 800/min, PUL 800/min,

NOTE: For details on ICT printing, contact your local A&D dealer.

12.22 Bluetooth® CONNECTION TIMING

F31: Bluetooth[®] connection timing setting.

DIA LED	ICT Printing	Default
1	Connect at the end of measurement – Connect with the host device after each	
1	measurement and start Bluetooth® transmission.	
	Connect at the start of measurement – Connect with the host device at the start	1
	of each measurement and start <i>Bluetooth®</i> transmission.	

12.23 Bluetooth® AIRPLANE MODE

F35: Airplane mode setting.

DIA LED	Airplane mode	Default
1	Airplane mode OFF	
2	Airplane mode ON	1

13. Transmission Specifications

The monitor can connect to the optional external input/output unit. Various settings for each channel are available from functions F20 to F25.

0	• The personal computer and medical equipment connected to the device must be located out of reach of the patient.		
	The personal computer or ID reader used must conform to EN60601-1.		

13.1 EXTERNAL INPUT/OUTPUT UNIT

Unit	Function	
TM-2657-01	Mini-DIN 8pin female, D-Sub 9pin male	
TM-2657-04	Bluetooth [®] Low Energy, D-Sub 9pin male	
TM-2657-05	Bluetooth [®] , D-Sub 9pin male	

Note: For details on EXTERNAL INPUT/OUTPUT UNIT (TM-2657-01, TM-2657-04, TM-2657-05), contact your local A&D dealer.

13.2 MINI-DIN 8 PIN FEMALE (EXTERNAL INPUT/OUTPUT UNIT : ONLY TM-2657-01)

Transmission specifications

Main standard	Complies with EIA RS-232C
Transmission format	Stop-start system (Full duplex)
Cignal anod	1200, 2400, 4800 and 9600 bps
Signal speed	(can be changed using F21)
Transmission format	Can be changed using F20
Data bit length	8 bits, 7 bits
Parity	None
Stop bit 1 bit, 2 bits (can be changed using F23)	
Code ASCII	

Pin Assignment



Pin No.	Signal Name	Description	
1	TXD	D Transmit data	
2	RXD	Receive data	
3	RTS Request to send		
4	– No connection		
5	CTS Clear to send		
6	GND Signal ground		
7	 – No connection 		
8	_	- No connection	

Note

Do not connect to Pins No. 4, 7, or 8. They are used for the blood pressure monitor.



13. Transmission Specifications

Cable specifications for computer connection

TM-2657WP- Mini-DIN 8 pin female

Content	Pin No.	Content	Pin No.
TXD	1	— —	1
RXD	2		2
RTS	3		3
—	4	DTR	4
CTS	5		5
GND	6		6
—	7	RTS	7
_	8	CTS	8
		—	9


Personal computer D-Sub 9 pin male

13.3 D-SUB 9-PIN MALE – (EXTERNAL INPUT/OUTPUT UNIT)

Transmission specifications

Output standards	Complies with EIA RS-232C
Transmission format	Stop-start system (Full duplex)
Signal speed	1200, 2400, 4800 and 9600 bps (can be changed using F22)
Transmission format	Can be changed using F20
Data bit length	8 bits
Parity	None
Stop bit	1 bit, 2 bits (can be changed using F24)
Code	ASCII

Pin Assignment



Pin No.	Signal Name	Description
1	_	-
2	RXD	Receive data
3	TXD	Transmit data
4	DTR	Data terminal ready
5	GND	Signal ground
6	DSR	Data set ready
7	RTS	Request to send
8	CTS	Clear to send
9	_	-

Note	The protocol depends on the equipment connected.
------	--

13. Transmission Specifications

Cable connection between the device and a personal computer

TM-2657WP – Sul D-sub connector	b 9 pin male		Personal computer D-sub connector	ter or ID Reader D-S	Sub 9 pin
Content	Pin No.		Content	Pin No.	
TXD	1		—	1	
RXD	2		RXD	2	
RTS	3		TXD	3	
—	4	Г	DTR	4	
CTS	5		GND	5	
GND	6	L	DSR	6	
—	7		RTS	7	
_	8		CTS	8	
				9	

Personal computer or ID Reader D-Sub 9 pin male

13.6 Bluetooth® (EXTERNAL INPUT/OUTPUT UNIT : TM-2657-04 & TM-2657-05)

In order to use the Bluetooth® transmission function of the TM-2657WP series safely and correctly, carefully read the following precautions before using the monitor. The following content summarizes general matters regarding the safety of patients and operators, in addition to the safe handling of the monitor.

Before using the monitor

\otimes	Do not use in places where wireless communication is prohibited, such as on airplanes or in hospitals. This monitor may have an adverse effect on electronic devices or medical electrical equipment.	
	If implantable heart pacemaker or implantable cardioverter defibrillator are used, please contact about the influence of radio waves individually to medical electrical equipment manufacture.	
9	 For such as warning and caution about the handling of sphygmomanometer body, please follow the description of the instruction manual of sphygmomanometer. 	
0	• This monitor has built-in wireless equipment with construction design certification as wireless equipment of a low electric power data communicating system based on regulations of the Radio Act. Therefore, when the wireless function of this equipment is used, wireless station permission is not necessary.	

• Disassembly or modification of this monitor may be punished by a law because this monitor has construction design certification.

During use of the wireless equipment

• We cannot accept any responsibility for any losses incurred such as operating malfunctions or loss of data that may occur through the use of this monitor.
This monitor is not guaranteed to connect with all <i>Bluetooth®</i> compatible devices.
• In the event of radio wave interference from the monitor to the other wireless station, change the location of use or stop using immediately.

For good wireless communication

\bigcirc	Do not use in the vicinity of cell phones. This could cause malfunction.

	• Ensure wireless device is within view of the monitor. Wireless range is affected by building structure and obstructions. Specifically, reinforced concrete can cause wireless interference.
Notes	• When communicating using this equipment close to a wireless communication device that is communicating using radio wave of near 2.4GHz, there are cases when the processing speed decreases at the both side. Do not use this equipment at a location where a magnetic field, static electricity or wave interference around the microwave ovens appears. (Doing so may prevent the devices from communicating properly.)
	 If the monitor cannot normally transmit data near a radio or broadcast station, please change the location.

13. Transmission Specifications

13.7 Bluetooth[®] (EXTERNAL INPUT/OUTPUT UNIT : TM-2657-05)

Transmission specifications

Main standard Bluetooth® Ver.2.1 class1	
Supported profiles	HDP
Communication distance	Maximum of 100m (depends on usage)
Frequency band	2,402 - 2480 MHz
Modulation	GFSK/QPSK
Maximum RF output power	< 20 dBm
	Continua certified devices
	 Applications and devices that are compatible with SSP and A&D specifications. However, each device needs an application to receive data. For connection methods, refer to the manual for each device. Bluetooth[®] devices described the Bluetooth[®] logo mark.
Devices that can be connected	 Bluetooth[®] Continua certified devices described with the Continua logo mark. Continua Continua

This monitor may be changed for improvement without any prior notice.

Pairing

A *Bluetooth*[®] device needs to be paired with a different specific device in order to communicate with that device. When this monitor is paired with a receiver device, measurement data is transmitted automatically to the receiver device each time a measurement is made.

Follow the steps below to pair the monitor with a *Bluetooth*[®] compatible receiver device. Also refer to pairing in the manual of the receiver device. Please use a pairing wizard if it provided.

- 1 Follow the instructions in the manual of the receiver device to switch it to the state that a pairing is possible. When pairing this monitor, place it as close as possible to the receiver device to be paired with.
- 2 Hold down the **SELECT** button and turn on the power. Press the button after "do" is displayed in the systolic display section and "PAr" is displayed in diastolic display section.

The monitor will be searchable from the receiver device for about one minute after pressing the S button.

- 3 Follow the manual of the pairing receiver device, the monitor performs a search, select, and pair. If a PIN code is requested by the receiver device, enter "123456".
- 4 "End" is displayed in the pulse rate display section when the pairing is over successfully on the receiver device side, and the pairing is finished.
- 5 If the pairing is failed, "Err" is displayed in the pulse rate display section. Turn off the monitor and back on again, and then retry from the step (1).

	• Other than the operation of the above (2), the monitor will be searchable from the receiver device for about one minute after turning on the power. In this operation, "End/Err" are not displayed in the pulse rate display section when the pairing is over.
Notes	 When reset with the FAST STOP button, searching is impossible.
	• Be sure to turn off the power of <i>Bluetooth</i> [®] devices other than the monitor when pairing. Multiple devices cannot be paired at the same time.

Measurement data transmission

Transmission after pairing is performed automatically by the following procedure. Enable wireless communication on the receiving device.

- 1 Press the Subtton to start blood pressure measurement.
- 2 After measurement, the measurement data is transmitted automatically to the receiver device.

Notes	• When the function setting F20 of the Automatic Blood Pressure Monitor on which the monitor is installed is OFF, data transmission and reception are not performed. Ensure F20 is not set to OFF.
	If the receiver device cannot receive measurement data, try pairing again.
	 The communication distance between this monitor and the receiver device is dependent on the Bluetooth[®] output class of the receiver device.
	When the receiver device is a Class 1 Bluetooth® device: Less than 100 m
	When the receiver device is a Class 2 Bluetooth® device: Less than 10 m
	This distance depends on the conditions in the surrounding environment. Please check that the distance is acceptable for transmitting measurement data.

In cases when the receiver device cannot receive measurement data, the measurement data is temporarily stored in the monitor memory along with the measurement time. A total of 200 sets of measurement data can be automatically stored. When the amount of data exceeds 200 sets, the oldest data is deleted and the new data is stored.

The data stored in the memory is transmitted the next time a connection is successfully made to the receiver device, and when the reception is confirmed, it is removed automatically. The amount of data that can be stored temporarily may vary with the receiver device.

Bluetooth® utility mode

Configure *Bluetooth*[®] settings for this monitor in *Bluetooth*[®] utility mode. To change function settings, use the buttons located on the rear panel of the monitor while the monitor is in standby mode.

- 1 Hold down the **SELECT** buttons and turn on the power. "do" is displayed in the systolic display section and "PAr" is displayed in the diastolic display section, when the *Bluetooth*[®] utility mode has started.
- 2 Each time the **SELECT** button is pressed, the setting changes to "un" / "PAr" → "cLr" / "dAt" → "do" / "PAr" →...
- 3 Each item can be performed using the 🏀 button.

Pairing

See "Pairing" described on previous page.

Unpairing


Devices can be unpaired.

- 1. Enter the *Bluetooth*[®] utility mode. Press the button with "un" in the systolic display section and "PAr" in the diastolic display section.
- 2. When "End" is displayed in the pulse rate display section, cancellation of the pairing is completed, but when "Err" is displayed in it, retry from the step (1).

Data clear

- 1. Erase data temporarily stored in the Automatic Blood Pressure Monitor.
- 2. Enter the *Bluetooth*[®] utility mode. Press the 🛞 button with "cLr" in the systolic display section and "dAt" in the diastolic display section.
- 3. When "End" is displayed in the pulse rate display section, cancellation of the data clear is completed, but when "Err" is displayed in it, retry from the step (1).

Note Th	This function is valid only with the TM-2657-05.
---------	--



13. Transmission Specifications

Time

This monitor has a built-in clock. The measurement data includes the date and time that a measurement was taken. The time is designed to be synced with the time of a receiver device side. Refer to the specifications of the receiver device side.

Notes	• The clock in the monitor can be automatically set by the receiver device side function. After the pairing, the time of the monitor is automatically set to the time of the receiver device 2 minutes after power on if there are no operations, or at the start of first measurement.
	• When the setting function F20 is off, the above clock synchronization is not performed.

Contents of transmission

Transmission data: Systolic blood pressure, diastolic blood pressure, pulse rate, measurement time, ID For more information, please contact A&D Customer Service.

13.8 Bluetooth® LOW ENERGY (EXTERNAL INPUT/OUTPUT UNIT : TM-2657-04)

Transmission specifications

Main standard	Bluetooth [®] Ver.2.4
Supported profiles	BLP (Blood Pressure Profile)
Frequency band	2.4 GHz (2400 - 2483.5 MHz)
Modulation	GFSK
Maximum RF output power	< 20 dBm
Maximum RF output power	< 20 dBm
Devices that can	 Applications and devices that are compatible with <i>Bluetooth</i>[®]. However, each device needs an application to receive data. For connection methods, refer to the manual for each device.
be connected	 Bluetooth[®] devices described the Bluetooth[®] logo mark.
	Bluetooth [°]

Communication data: SYS, DIA, MAP, PUL, Measurement time, ID

Bluetooth® Low Energy utility mode

Configure *Bluetooth*[®] settings for TM-2657-04 in *Bluetooth*[®] Low Energy utility mode. To change function settings, use the buttons located on the rear panel of the monitor while the monitor is in standby mode.


- 1 Hold down the **SELECT** button and turn on the power.
- 2 "INT" is displayed in the systolic display section and "____" is displayed in the diastolic display section, when the initialize has started.
 "E98" is displayed in the systolic display section and "-4" is displayed in the diastolic display section, when the initialize has failed. Turn off the monitor and back on again, and then retry from the step (1).
- "do" is displayed in the systolic display section and "PAr" is displayed in the diastolic display section, when the *Bluetooth*[®] Low Energy utility mode has started.
- 4 Each item can be performed using the 🛞 button.

Pairing

A *Bluetooth*[®] device needs to be paired with a different specific device in order to communicate with that device. When this monitor is paired with a receiver device, measurement data is transmitted automatically to the receiver device each time a measurement is made. Follow the steps below to pair the monitor with a *Bluetooth*[®] compatible receiver device. Also refer to pairing in the manual of the receiver device. Please use a pairing wizard if it provided.

- 1 Follow the instructions in the manual of the receiver device to switch it to the state that a pairing is possible. When pairing this monitor, place it as close as possible to the receiver device to be paired with.
- 2 Hold down the **SELECT** button and turn on the power. Press the button after "do" is displayed in the systolic display section and "PAr" is displayed in diastolic display section.

The monitor will be searchable from the receiver device for about one minute after pressing the 👹 button.



- 3 Follow the manual of the pairing receiver device, the monitor performs a search, select, and pair.
- 4 "End" is displayed in the pulse rate display section when the pairing is over successfully on the receiver device side, and the pairing is finished.
- 5 If the pairing is failed, "Err" is displayed in the pulse rate display section. Turn off the monitor and back on again, and then retry from the step (1).

	This function is valid only with the TM2657-04.
Notes	• When the setting item F35 in the function settings is set to 2 (Airplane mode ON), pairing cannot
	achieve. Set to 1 (Airplane mode OFF) for a pairing.

Measurement data transmission

Transmission after pairing is performed automatically by the following procedure. Enable wireless communication on the receiving device.

- 1 Press the button to start blood pressure measurement.
- 2 After measurement, the measurement data is transmitted automatically to the receiver device.

	• When the function setting F20 of the Automatic Blood Pressure Monitor on which the monitor is installed is OFF, data transmission and reception are not performed. Ensure F20 is not set to OFF.
	If the receiver device cannot receive measurement data, try pairing again.
	The measurement data will be erased when the monitor is turned off.
Notes	 The communication distance between this monitor and the receiver device is dependent on the Bluetooth[®] output class of the receiver device.
	The communication distance when there are no obstacles is about 10 m.
	• This distance depends on the conditions in the surrounding environment. Please check that the distance is acceptable for transmitting measurement data.

Time

This monitor has a built-in clock. The measurement data includes the date and time that a measurement was taken. The clock in the monitor can be set by sending a command from the receiver device.

Cybersecurity features for each communication interface:

- BLE, BT: Authentication via pairing
- · RS232C, USB: No access control, no encryption (measurement data transmitted unencrypted)

14.1 INSPECTION AND SAFETY MANAGEMENT

Do not open the device. It uses delicate electronic components and an intricate air unit that could be damaged. Do not change the variable resistor setting of the power supply. The device may not turn on if the settings are changed. If you cannot fix the problem using the troubleshooting instructions, request service from your local dealer or from A&D customer service. A&D customer service will provide technical information, spare parts and units to authorized dealers. Technical inspection procedures which should be done at least every two years, can be performed either by the manufacturer or by an authorized repair service in accordance with the regulations governing manufacturing of medical products


Checking pressure accuracy

	0	• When using a rubber pump, do not apply a pressure of 280 mmHg or higher to the monitor or inspection equipment (UM-102, accurate mercury sphygmomanometer or aneroid gauge).	
	Q	 Perform the inspection only as described below or the setting values and function settings may be changed. 	
	0	After inspection, check that the air connector plug is inserted into the blood pressure monitor. If the air connector plug is not inserted, pressure cannot be applied and measurement is not possible. When inserting the plug, push in until you hear a click.	

Objective: Compare the pressure values of the inspection equipment and the blood pressure monitor to check for errors in the monitor.

Inspection equipment: Inspection equipment (UM-102, Accurate mercury sphygmomanometer or aneroid gauge)

Connection: Connect the inspection equipment to the blood pressure monitor as shown below. Remove the armrest of the blood pressure monitor and then remove the cover of pressure inspection area. Remove the air connector plug from the air socket of the blood pressure monitor. Connect the coupling connector to the connection hose, and connect it to the air socket.



1. Hold the **COUNT** button on the rear of the blood pressure monitor, and turn the **POWER** switch on.

2. "L BD" appears in the clock display_section.

With "L∃⊡" displayed, press the button.
 Pressure inspection mode starts and the current pressure is displayed.

Using the rubber pump, apply the pressures listed below. Compare and check the pressures of the blood pressure monitor and the inspection equipment.

No	Pressure Setting	Instrumental error A-B (standard)
1	0 mmHg	0 mmHg
2	50 mmHg	Within ±6 mmHa
3	200 mmHg	within to himng

A: Pressure displayed by the inspection equipment

B: Diastolic and systolic pressures displayed by the monitor

5. Confirm that the values are within standards. To exit the pressure inspection mode and return to the standby mode, switch the power off and switch the power on again.

Note	Use the coupling connector for exclusive use with the TM-2657WP series.
------	---

14.2 CLEANING

	Before cleaning, switch the power off and disconnect the power cable from the electrical outlet.		
	When cleaning the monitor, never splash it with or soak it in water.		
 The blood pressure monitor is not waterproof device. Do not splash water on it and avoid exposure to moisture. 			
 When disinfecting the monitor, never use an autoclave or gas sterilization (EOG, formaldehy concentration of ozone). 			
Never clean the monitor with solvents such as thinner or benzene.			
	 Clean the monitor about once a month in the following manner based on policies and procedures determined by the hospital. 		

When the main body or the arm cuff cover is dirty, wipe them fully by using gauze or cloth dampened with warm water and a neutral detergent avoiding excess water.

To prevent a risk due to infection, disinfect the main body and the arm cuff cover regularly. When disinfecting them, wipe them gently by using the gauze or dampened cloth with local antiseptic solution then wipe the moisture off the surface by using a dry soft cloth.

The antiseptic solution should be used as a water solution by following a rule for notes for its product at the dilution ratio. The following shows the example in which can be used as antiseptic solution.

Example of useable antiseptic solution (Ingredient name)

Component Name	Product Name
Benzalkonium chloride	Benzalkonium chloride 10% solution
Isopropanol	70% in 1-propanol
Ethanol	Ethanol for disinfection 76.9 to 81.4 vol%

Check that the arm cuff cover is not damaged. If it is damaged, replace it. For the replacement procedure, see "Replacing the arm cuff cover".

	The arm cuff cover and cables are consumables. If there are frequent measurement errors or
Note	measurement is not possible, these items must be replaced. Before ordering replacements, see "15.
	ACCESSORIES AND OPTIONS LIST".


Printer head

If the printer head has paper debris, or other foreign matter has collected, printing will not be performed correctly. To prevent this, follow the procedure below to clean the printer head.




CAUTION
 Before cleaning, switch the power off and wait until the printer head has cooled completely. The printer head gets very hot and may cause burns.


- Some printer parts have sharp edges. Take great care when handling them to avoid injury.
- 1. Switch the power off.
- 2. Press the **Open printer cover** button to open the printer cover.




3. Using a soft cotton swab or cotton cloth moistened with alcohol (ethyl or isopropyl), clean the heating element very gently.



- 4. Clean the printer paper compartment to remove dust, paper debris and other foreign matter. Debris in the paper output path may lower the printing quality.
- 5. Wait for the cleaned parts to completely dry and install the printer paper.



6. With the end of the paper at the top and protruding out, secure the printer paper by closing the printer cover until you hear a click. If the cover is not completely closed, a paper jam may occur.



	•	When cleaning the printer head, be careful of static electricity. Static electricity can damage the printer head.
Notes • Do not use abrasive substances, such as sandpaper, to clean the pr heating element.	Do not use abrasive substances, such as sandpaper, to clean the printer head. They will damage the heating element.	
	•	Make sure that the printer head is completely dry before installing the printer paper and switching the power on.

14.3 PERIODIC INSPECTION

To ensure correct use of the monitor, perform a periodic inspection. The main items of the periodic inspection are as follows.

Before switching the power on


Item	Description		
	Check for deformations and damage from drops.		
Exterior	Check parts for dirt, rust, scratches.		
Exterior	Check panels for dirt, scratches, damage.		
	Check for moisture.		
Operation parts	Check switches and buttons for damage, looseness.		
Display Check display for dirt, scratches.			
Measurement parts Check the cuff and arm cuff cover for damage.			
A	Check that the arm cuff cover is installed.		
Arm cuff cover	Please use the arm cuff cover to prevent any foreign matter from entering into this device.		
Printer Check that the printer paper is the specified type			
	Check that the power cable is inserted correctly into the connector.		
Dowor porto	Check the power cable for damage (exposed core wires, disconnection).		
Power parts	Check that the electrical outlet is properly grounded and supplies the specified voltage and		
	frequency (100-240 V~ 50-60 Hz).		

After switching the power on


Item	Description		
Exterior	Check for smoke or unusual smells.		
EXTERIO	Check for unusual noise.		
Operation parts	Press the 🛞 button and check for errors.		
	Press the FAST STOP button during inflation to check that pressurization stops.		
	Check the blood pressure, pulse and clock display sections for missing numbers or characters.		
Display	Check that no error codes are displayed.		
	Check that measurement values are near normal values.		
	Check that the paper availability and run out are detected.		
Printer	Check that the printer paper is fed correctly.		
Finiter	Check that test printing has no missing items.		
	Check that the paper is cut after printing.		
Pookup function	Check that the date and time are correct.		
Backup function	Check that the contents of set values are saved.		

14.4 REPLACING THE ARM CUFF COVER


- 1. Use a flathead screwdriver to loosen the screw.
- 2. Slide the front frame down, and then pull forward.




3. Loosen the screws (armrest securing screws) on the rear side and remove the screws.




- 4. Lift the armrest and pull back.
- 5. Slide the rear frame down, then pull out.




6. Pull the arm cuff cover out from the vinyl ring groove to remove.



 Insert the new arm cuff cover and push the vinyl ring into the groove (on the inner side of the frame) to attach.



- 8. Fit the new arm cuff cover over the front vinyl ring groove.
- Reversing the steps used to remove, reattach the rear and front frames, return the armrest to its original position, then replace the armrest securing screws (2) and front frame screw (1).



Note

The arm cuff cover is consumable. New covers must be purchased separately. (arm cuff cover: AX-134005759-S)



Using a correct arm cuff cover and exchanging it are important for safety and measurement accuracy at this device.

14.5 CHECKING THE NUMBER OF MEASUREMENTS

The monitor can count the number of times blood pressure measurement has been performed. This function is designed to check usage frequency and provide a reference for scheduled cleaning. The count value is stored even after the power is switched off.

Displaying the Number of Measurements

To display the number of measurements:

Hold the **COUNT** button for 1 second while the monitor is in the standby mode. The number of measurements is displayed for about 60 seconds in the systolic and diastolic display sections.

In the example display below, the number of measurements is 2,382. (The maximum count is 999,999.)



Display thousands digits and higher


 Display hundreds, tens, and ones digits


Appears when displaying number of measurements

To reset the number of measurements:

Hold the ▲ button for 4 seconds to display the reset confirmation display.

Press the 🛞 button to reset the count.





Printing the Count Graph

To print the count graph

Press the **COUNT** button. While the number of measurements is displayed, press the button to print the count graph.


Total Count: Number of measurements since shipping

Trip Count: Number of measurements since the last reset

(See "Displaying the number of measurements")

Weekly Count: A distribution of the number of measurements in the last week.

Monthly Count: A distribution of the number of measurements in the last month.



	If the function F07 is set to off, the count graph is not printed. (See "Print quality")
Notes	• After the count graph is printed, the number of measurements remains displayed for about 60 seconds.
	 If "Low Battery" is printed in the lower left of the print out after the count graph is printed out, please contact your local A&D dealer.

14.6 DISPOSING OF THE COMPONENT PARTS

Dispose of or recycle the monitor in an environmentally friendly manner according to local regulations.

Arm cuff cover

As there is a danger of infection, dispose of the arm cuff cover as medical waste.

Internal backup battery

The monitor is equipped with a lithium battery to back up settings and other data. Before disposing of the main unit, remove the lithium battery and dispose of it according to local regulations.

Product Name	Structure Name	Material
	Box	Cardboard
Package	Packing material	Cardboard
	Bag	Vinyl
	Case	ABS/ABS plastic
Inside main unit	Internal parts	General parts
	Chassis	Steel
	Battery on PCB	Lithium battery
	Case	ABS/ABS plastic
Printer unit	Internal parts	General parts
	Chassis	Steel
External input/output unit (Option)	Case	ABS/ABS plastic
External input/output unit (Option)	(Option)	General parts

14.7 TROUBLESHOOTING

Before requesting service, please review the following checklist and the error code list in the next section.

Problem	Check	Countermeasure	
Nothing is displayed when the power is switched on.	Is the power cable connected correctly?	Connect the power cable correctly.	
E00 is displayed.	Is there air remaining in the cuff?	Wait until the air is released completely from the cuff, and then switch the power on again.	
There is no pressure.	Is the arm cuff cover pulled too far over the frames?	· · ·	
	Is the patient's posture correct?	Ensure that the arm and heart are at the same height and that the patient is relaxed.	
Measurement is not possible.	Is the patient relaxed?	Ensure that the patient does not move their arm.	
(An error code is displayed.)	Is clothing too thick? If so, measurement is not possible.	Remove the clothing from the arm.	
	Does the patient have arrhythmia or a weak pulse?	Measurement may not be possible with patients with arrhythmia or a weak pulse.	
	The printer paper is not installed. $(\frac{1}{2} \frac{1}{2}$ is displayed)	See "9.1. Installing the printer paper" to install a new roll of printer paper.	
No printing	The printer cover is open. $(\stackrel{\square}{=} \Box$ is displayed)	See "Installing the printer paper" to close the printer cover.	
no printing	A printer cutter error. (^[2] ⊂ is displayed)	See "9.1. Installing the printer paper" to temporarily open the printer cover and then close it again.	
	Is the printer paper causing a jam?	See "9.1. Installing the printer paper", readjust the paper.	
The printing content was not as expected.	Is the printing method selection appropriate?	See Sections "10.4. IHB" to "10.10. Bitmap printing" to select the printing method.	
	Check the clock setting.	Refer to "8. SETTING THE CLOCK"	
Date and/or time are off.	Is the Low Battery printed on the lower left of the print out after the count graph is printed as shown in 12.5.2?	The lithium battery for back up settings and other data is dead. Contact your local A&D dealer.	
	Check the clock setting on the Bluetooth [®] receiver.	See the specifications of the receiver device.	

 \otimes

Do not touch the interior of the monitor.

14.8 ERROR CODES

When an error occurs, one of the following error codes is displayed in the systolic display section.

Printer error codes

Error code	Error/countermeasure	
PE	o printer paper. Install a new roll of printer paper.	
Po	The printer cover is open. Firmly close the printer cover.	
Pc	A printer cutter error. Open the printer cover, check the printer paper, and then close the printer cover.	

Error code details

Error code	Details	Check items	
Error related	to blood pressure measuremen		
E00	When the power is switched on, the pressure detection is unstable.	Check if there is air remaining in the cuff. Restart and then try blood pressure measurement again. If the problem continues, stop using the monitor immediately.	
803	An electrical error is detected in the blood pressure measurement section.	Restart and then try blood pressure measurement again. If the problem continues, stop using the monitor immediately.	
609	The safety monitor of the blood pressure measurement section detected an error.	A condition that may affect the safety of the patient was detected during measurement. External vibrations may have been applied to the air system of the cuff or inside the monitor or an obstruction may have been mistakenly detected. Check the patient condition and measurement environment and try blood pressure measurement again. If the problem	
E I I,E IS	Pressure is not applied at the start of the measurement.	continues, stop using the monitor immediately. There may be an air leak in the air system inside the monitor. If the problem continues, stop using the monitor immediately.	
E 12	Pressure cannot be applied within a certain period of time.	There may be a leak in the air system inside the monitor or the cuff was applied loosely. If the problem continues, stop using the monitor.	
E 13	Inflation speed is too fast.	There may be a bend or blockage in the air system inside the monitor. If the problem continues, stop using the monitor.	
851	The exhaust speed is too slow.	Air is not being correctly exhausted. There may be a bend or blockage in the air system inside the monitor. If the problem continues, stop using the monitor.	
523	The exhaust speed is too fast.	The patient may have moved or a strong external pressure was applied during measurement. If the problem continues, stop using the monitor.	
623	Excess pressure was detected.	The cuff pressure during measurement exceeded 300 mmHg. The patient may have moved or a strong external pressure was applied t the cuff.	
624	The time limit for one measurement was exceeded.	Watch for errors and try measurement again. For the safety of the patient, measurement was cancelled because the measurement time exceeded 180 seconds. Measurement may have been repeated. Check the patient for body movement and arrhythmia.	
E42	The pressure is insufficient.	Blood pressure measurement was not possible because the pressure was insufficient. During inflation, patient movement or an external vibration introduced noise into the cuff pulse and the set pressure was detected or the patient's blood pressure rose greatly during blood pressure measurement.	
		Confirm the following conditions: The cuff is not loose; no thick clothing on the arm; the patient remains sti and no external vibrations on the cuff. And try measurement again.	
643	Pulse cannot be detected.	The pulse signal received by the cuff is too low. The circulation of the patient may be poor or the patient is wearing thick clothing. Check the condition of the patient.	
E45	Diastolic blood pressure cannot be determined.		
846	Mean arterial blood pressure cannot be determined.	Check the patient for body movement and arrhythmia.	
848	Systolic blood pressure cannot be determined.		
E6 (Pulse cannot be determined.		
663	The blood pressure value is inappropriate.	Check the patient for body movement and arrhythmia.	
863 1	SYS value is 'out of range'.	SYS measurement range : 40-270 mmHg Check the patient for body movement and arrhythmia.	
663 2	DIA value is 'out of range'.	DIA measurement range : 20-200 mmHg Check the patient for body movement and arrhythmia.	
663 3	PUL value is 'out of range'.	PUL measurement range : 30-240 mmHg Check the patient for body movement and arrhythmia.	

Error code	Details	Check items
Other errors		
697	Restart the power.	Restart the power.
i to 4	A power voltage error was detected inside the monitor.	If the problem continues, stop using the monitor immediately.
697	Restart the power.	The function settings have been initialized. Check the settings. Restart the
S	A setting error was detected inside the monitor.	power. If the problem continues, stop using the monitor immediately.
697	Restart the power.	The counting function has been initialized. Restart the power. If the
6	A setting error was detected inside the monitor.	problem continues, stop using the unit for the time being.
697	Restart the power.	Restart the power.
8,9	A setting error was detected inside the monitor.	If the problem continues, stop using the monitor immediately.
898	Restart the power.	
	A memory error was detected inside the monitor.	
898	Restart the power.	
3	A USB error was detected inside the monitor.	Restart the power. If the problem continues, stop using the monitor immediately.
	Restart the power.	
698 Ч	Bluetooth [®] Low Energy error was detected inside the monitor.	
899	There may be a malfunction.	
1	A font error was detected.	
699	There may be a malfunction.	Restart the power.
- 2	A cuff error was detected.	If the problem continues, stop using the monitor immediately and request
699	There may be a malfunction.	repairs.
3	A blood pressure module error was detected.	

Displaying the error status

Press the **COUNT** button. The count is displayed. Press the **SELECT** button within 60 seconds. The past error codes (systolic display section), error sub codes (diastolic display section) and the number of occurrences (pulse display section) are displayed. Each time the **SELECT** button is pressed, past error codes are displayed in numerical order.

After 60 seconds of no operation, the monitor returns to standby mode.

15. Accessories and Options List

Product Name	Product Description	Catalog Number
Printer paper	5 rolls	AX-PP147-S
Arm cuff cover	5 pieces	AX-134005759-S
Disposable Kiosk Sleeves	Pack of 400	AS-134010367
Dedicated Stand		TM-ST520
Adjustable Height Stool		TM-STA001
Power cable	Cord set Type C	AX-KO243
Power cable	Cord set Type BF Fuse rating: T3AH250V	AX-KO242
Power cable	Cord set Type A	AX-KO115-EX
External input/output unit	RS 2ch	TM-2657-01
External input/output unit	RS+Bluetooth® Low Energy	TM-2657-04
External input/output unit	RS+Bluetooth®	TM-2657-05

16. Performance Specifications and Standards

16.1 PERFORMANCE SPECIFICATIONS

General

AC Power supply	100 V-240 V ~ 5 0 Hz-60 Hz		
Power consumption	50-80 VA		
Safety standard	IEC60601-1:2020		
EMD compliance	Complies with EMD standard IEC60601-1-2:2020.		
Type of protection	Type B 🛧 Applied part: Cuff		
Type of protection against electrical shock	Class I		

Blood pressure measurement

Measurement method	Oscillometric measurement		
Pressure display range	0-299 mmHg		
Pressure display accuracy	Pressure: ±3 mmHg		
	SYS 40-270 mmHg		
NIBP Measurement range	DIA 20-200 mmHg		
	Pulse rate 30-240 bpm		
Measurement intervals	1 mmHg		
NIBP Clinical test	ISO81060-2:2018		
Pulse rate accuracy	±5%		
Cuff	Winding mechanism operated by geared motor		
Applicable arm circumference	7.1-16.5" / 18-42 cm		
Inflation	Automatic inflation by air pump		
Deflation	Automatic deflation by mechanical exhaust		
Rapid deflation	Automatic rapid deflation by solenoid valve		

Environment specifications

Operating environment	Temperature: 50 to104°F / 10 to 40°C
Operating environment	Humidity: 15-85% RH (no condensation)
Storage environment	Temperature: -4 to 140°F / -20 to 60°C
Storage environment	Humidity: 10-95% RH (no condensation)
Atmospheric pressure range 70-106 kPa (both for operation and storage)	

Physical specifications

External dimensions	9.5 (W) x 12.8 (H) x 15.3" / 241 (W) x 324 (H) x 390 (D) mm	
Weight	Approx. 12.1 lb / 5.5 kg	

Functional specifications

16. Performance Specifications and Standards

Display method	3-digit display LED & LED lamp		
Printer	Thermal printing, paper width: 58 mm		
	5 years from installation		
Usable life	According to A&D data (tested for use under recommended environment, including maintenance inspection. Results may be different under other conditions.)		

Transmission specifications

Standard	USB 2.0 (Option)	Bluetooth [®] Low Energy (Option)

16.2 STANDARDS

The TM-2657WP Series Automatic Blood Pressure Monitor complies with the following standards:

- IEC 60601-1:2020 (Medical electrical equipment Part 1: General requirements for safety and essential performance);
- IEC 60601-1-2:2014 +A1:2020 (Medical electrical equipment Part 1-2: General requirements for basic safety and essential performance Collateral Standard: Electromagnetic disturbances Requirements and tests);
- EN ISO81060-2:2018(Non-invasive sphygmomanometers Part 2: Clinical investigation of intermittent automated measurement type)
- IEC 80601-2-30: 2018 (Medical electrical equipment –Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers).
- The TM-2657WP Series is not made with natural rubber latex.

17. Warranty

LIMITED WARRANTY

A&D Medical

For purchasers within the US or Canda only:

Product	Warranty Term	
Monitor	2 year	

For outside of US and Canada, please contact local distributor or dealer.

Limited Warranty:

A&D Medical ("A&D") warrants to the first purchaser ("You") that the A&D product You purchased (the "Product") will be free from defects in material, workmanship and design for the applicable Warranty Term stated above from the date You purchased the Product under normal use. This Limited Warranty is personal to You and is not transferable. If the Product is defective, then You return the Product to A&D in accordance with the procedure set forth below. A&D's warranty obligation is limited to the repair or replacement, at A&D's option, of the defective Product that has been returned by You within the warranty period. Such repair or replacement will be at no charge to You. The repaired or replacement Product is warranted here-under for the longer of the remainder of the original warranty period or 90 days from the date of shipment of the repaired or replacement Product.

To obtain a warranty service, please contact us in **US at 1-888-726-9966** or in **Canada at 1-800-461-0991** for return address, shipping and handling fee, and other instructions for processing warranty. Please ensure you have satisfactory proof of the date of Your purchase and a description of the defect. Returns will not be accepted unless a Return Material Authorization (RMA) Number has been issued from A&D Customer Service Representative.

This Limited Warranty does not cover, and A&D will not be liable for (i) any shipment damage, (ii) any damage or defect due to misuse, abuse, failure to use reasonable care, failure to follow written instructions enclosed with the Product, accident, subjecting the Product to any voltage other than the specified voltage, improper environmental conditions, or modification, alteration or repair by anyone other than A&D or persons authorized by A&D, or (iii) expendable or consumable components.

THIS LIMITED WARRANTY IS THE ONLY WARRANTY PROVIDED BY A&D; THERE ARE NO OTHER EXPRESS WARRANTIES. If A&D cannot reasonably repair or replace the Product, A&D will refund the amount You paid for the Product (not including taxes), less a reasonable charge for usage. To receive a refund you must have returned the Product and all associated materials to A&D. The above remedy of repair, replacement or refund is your only and exclusive remedy. IN NO EVENT SHALL A&D BE LIABLE FOR ANY DAMAGES, INCLUDING, WITHOUT LIMITATION, LOST PROFITS, LOST INFORMATION OR REPLACEMENT COSTS, ARISING OUT OF YOUR USE OF OR INABILITY TO USE THE PRODUCT, INCLUDING, WITHOUT LIMITATION, ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, EVEN IF A&D HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. Some states do not allow the exclusion of incidental or consequential damages, so that the above exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that may vary from state to state.

No distributor, dealer or other party is authorized to make any warranty on behalf of A&D or to modify this warranty, or to assume for A&D any liability with respect to its products.

Medical electrical equipment requires special precautions regarding EMD and must be installed and put into service according to the EMD information provided below.

Portable and mobile RF communication equipment (e.g. cell phones) can affect medical electrical equipment.

The use of accessories and cables other than those specified may result in increased emissions or decreased immunity of the unit.

EMISSION Limits -

Phenomenon		Compliance
Conducted and radiated RF EMISSION	CISPR 11	Group 1, Class B
Harmonic distortion	IEC 61000-3-2	Class A
Voltage fluctuations and flicker	IEC 61000-3-3	Compliance

IMMUNITY TEST LEVELS: Enclosure Port -

Phenomenon		IMMUNITY TEST LEVELS
Electrostatic discharge	IEC 61000-4-2	±8 kV contact
	ILC 01000-4-2	±2 kV, ±4 kV, ±8 kV, ±15 kV air
		10 V/m
Radiated RF EM fields	IEC 61000-4-3	80 MHz - 2.7 GHz
		80 % AM at 1 kHz
Proximity fields from RF wireless		See table (Test specifications for ENCLOSURE PORT
communications equipment	IEC 61000-4-3	IMMUNITY to RF wireless communications equipment)
Rated power frequency magnetic fields		30 A/m
	IEC 61000-4-8	50 Hz or 60 Hz
Proximity magnetic fields		See table (Test specifications for ENCLOSURE PORT
	IEC 61000-4-39	IMMUNITY to proximity magnetic fields

IMMUNITY TEST LEVELS: Input a.c. power Port -

Phenomenon		IMMUNITY TEST LEVELS	
Electrical fast transients / bursts	IEC 61000-4-4	±2 kV 100 kHz repetition frequency	
Surges, Line-to-line	IEC 61000-4-5	±0.5 kV, ±1 kV ±0.5 kV, ±1 kV, ±2kV(Line-to-Earth)	
Conducted disturbances induced	by RF fields IEC 61000-4-6	 3 V 0.15 MHz - 80 MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz 	
Voltage dips	IEC 61000-4-11	$ \begin{array}{c ccccccccccccccccccccccccccccccccccc$	
Voltage interruption	IEC 61000-4-11	I 0% U ₁ ; 250/300 cycle	
NOTE U_{τ} is the AC mains voltage prior to application of the test level.			

Appendix: EMD Information

IMMUNITY TEST LEVELS: Signal input/output Port -

Phenomenon		IMMUNITY TEST LEVELS	
Electrostatic discharge	IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	
Electrical fast transients / bursts IEC 61000-4-4		±1 kV	100 kHz repetition frequency
Conducted disturbances induced by RF fields IEC 61000-4-6		3 V 6 V in IS 80 %	0.15 MHz - 80 MHz M and amateur radio bands between 0.15 MHz and 80 MHz AM at 1 kHz

Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment -

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380 - 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 - 470	GMRS 460 FRS 460	FM ±5 kHz deviation 1 kHz sine	2	0.3	28
710	704 - 787	LTE Band 13,17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800 - 960	GSM 800/900		2	0.3	28
870		TETRA 800	Pulse modulation			
930		CDMA 850	18 Hz			
930		LTE Band 5				
1720	1700 - 1990	GSM 1800		2	0.3	28
1845		CDMA 1900	1			
1970		GSM 1900	Pulse modulation			
		DECT	217 Hz			
		LTE Band 1,3,4,25				
		UMTS				
		Bluetooth®				
2450	2400 - 2570	WLAN 802.11 b/g/n	Pulse modulation 217 Hz	2	0.3	28
		RFID 2450				
		LTE Band 7				
5240	5100 - 5800	WLAN 802.11 a/n	Pulse modulation	ion 0.2	0.3	9
5500			217 Hz			
5785						

Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment -

Test frequency	Modulation	IMMUNITY TEST LEVEL(A/m)	
30kHz	CW	8	
134.2kHz	Pulse modulation 2.1 kHz sine	65	
13.56MHz	Pulse modulation 50kHz	7.5	



Manufactured by (Fabriqué par):

A&D Company, Limited 1-243 Asahi, Kitamoto-shi, Saitama-ken 364-8585 Japan

Manufactured for (Fabriqué pour): A&D Engineering, Inc. 4622 Runway Boulevard Ann Arbor, MI 48108 USA 1-888-726-9966

Distributed by (Distribué par): Auto Control Médical – an A&D Company / une compagnie A&D 6695 Millcreek Drive, Unit 6, Mississauga, Ontario, L5N 5R8 Canada 1-800-461-0991

© 2024 A&D Company, Limited. All rights reserved.

No part of this publication may be reproduced, transmitted, transcribed or translated into any language in any form by any means without the express written consent of A&D Company, Limited.

The contents of this manual and the specifications of the instruments covered by this manual are subject to change without notice.

Windows is a registered trademark of the Microsoft Corporation.

The Bluetooth $^{\circ}$ word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by A&D is under license.

QR code is a registered trademark of DENSO WAVE Incorporated.

Other trademarks and trade names are those of their respective owners.