VOICE GUIDED

Blood Pressure Monitor
Premium Upper Arm
UA-1040TCNBL

Instruction Manual
ORIGINAL

Table of Contents

1. Introduction	E-3
Intended Use / FCC Compliance / Precautions	
About Blood Pressure / About Hypertension	
About Irregular Heartbeat / What is AFib? / Pressure Bar Indicator	
Tips for Taking Blood Pressure	
AHA Classification	
2. Getting to Know Your Monitor	E-11
Parts Identification / Display	
Symbols	
3. Initial Setup	E-16
Installing and Changing Batteries	
Connecting and Disconnecting AC Adapter / Connecting Air Hose	
Adjusting Clock / Adjusting Pressure	
TriCheck Mode / Adjusting Volume	
Selecting Language / Turning Instructions On / Off	
Pairing with a Mobile Device	
Selecting the Correct Cuff Size	
4. Taking Your Blood Pressure	E-24
Applying and Removing the Cuff	
Normal Measurement	
Using TriCheck	
Measurement with Set Pressure	
Measurement with Desired Systolic Pressure	
Reviewing and Deleting Memory	
5. Maintenance	E-31
6. Troubleshooting	E-32
7. Optional Accessories	E-34
8. Warranty	E-35
9. Technical Data	E-37

Congratulations on purchasing a state-of-the-art A&D blood pressure monitor. Designed for ease of use and accuracy, this monitor will facilitate your daily blood pressure regimen. We recommend that you read through

Intended Use

- The monitor is designed for use on adults only. Do not use on newborns or infants.
- Environment for use. The monitor is intended for use in the home healthcare environment.
- This monitor is designed to measure blood pressure and pulse rate of people for diagnosis.

FCC Compliance Information

This monitor complies with Part 15 of FCC rules. Operation is subject to the following two conditions:

1. This monitor may not cause harmful interference, and
2. This monitor must accept any interference received, including interference that may cause undesired operations.

Precautions

- This is a medical device for use by lay person. Please consult your healthcare provider with any questions or concerns you may have regarding your condition.
- If any serious incident occurs in relation to this device, report to its manufacturer and the competent authority in your country.
- This device complies with Part 15 of FCC rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operations.

Introduction

- This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions: (1) This device may not cause interference. (2) This device must accept any interference, including interference that may cause undesired operation of the device.
- Precision components are used in the construction of this monitor. Extremes in temperature, humidity, direct sunlight, shock or dust should be avoided.
- Clean the monitor and cuff with a dry, soft cloth or a cloth dampened with water and a neutral detergent. Never use alcohol, benzene, thinner or other harsh chemicals to clean the monitor or cuff.
- Avoid tightly folding the cuff or storing the hose tightly twisted for long periods, as such treatment may shorten the life of the components.
- Do not excessively bend or squeeze the air hose.
- The monitor and cuff are not water resistant. Prevent rain, sweat and water from soiling the monitor and cuff.
- Measurements may be distorted if the monitor is used close to televisions, microwave ovens, cellular telephones, X-ray or other devices with strong electrical fields.
- When using the monitor, confirm that the monitor is clean.
- Used equipment, parts and batteries are not treated as ordinary household waste, and must be disposed of according to the applicable local regulations.
- When the AC adapter is used, make sure that the AC adapter can be readily removed from the electrical outlet when necessary.
- Do not modify the monitor. It may cause accidents or damage to the monitor.
- To measure blood pressure, the arm must be squeezed by the cuff hard enough to temporarily stop blood flow through the artery. This may cause pain, numbness or a temporary red mark to the arm. This condition will appear especially when measurement is repeated successively. Any pain, numbness, or red marks will disappear with time.

- Do not apply the cuff on an arm in which another medical device is attached. The equipment may not function properly.
- People who have a severe circulatory deficit in the arm must consult a doctor before using the monitor, to avoid medical problems.
- Do not self-diagnose the measurement results and start treatment by yourself. Always consult your doctor for evaluation of the results and treatment.
- Do not apply the cuff on an arm with an unhealed wound.
- Do not apply the cuff on an arm receiving an intravenous drip or blood transfusion. It may cause injury or accidents.
- Do not use the monitor where flammable gases such as anesthetic gases are present. It may cause an explosion.
- Do not use the monitor in highly concentrated oxygen environments, such as a high-pressure oxygen chamber or an oxygen tent. It may cause a fire or explosion.
- Take care to avoid accidental strangulation of babies or infants with the hose and cable.
- Do not twist the air hose during measurement. This may cause injury due to continuous cuff pressure.
- Wireless communication devices, such as home networking devices, mobile phones, cordless phones and their base stations, walkie-talkies can affect this blood pressure monitor. Therefore, a minimum distance of 30 cm (12") should be kept from such devices.
- Measuring blood pressure too frequently may cause harm due to blood flow interference. Check that the operation of the monitor does not result in prolonged impairment of blood circulation, when using the monitor repeatedly.
- Clinical testing has not been conducted on newborn infants and pregnant woman.
- If you have had a mastectomy or lymph node clearance, please consult a doctor before using the monitor.
- Do not let children use the monitor by themselves and do not use the monitor in a place within the reach of infants. It may cause accidents or damage.

Introduction

- There are small parts that may cause a choking hazard if swallowed by mistake by infants.
- Do not touch the batteries, the DC jack and the patient at the same time. That may result in electrical shock.
- Unplug the AC adapter when not in use during the measurement.
- Use of accessories not detailed in this manual may compromise safety.
- Should the battery short-circuit, it may become hot and potentially cause burns.
- Allow the monitor to adapt to the surrounding environment before use (about one hour).
- Do not inflate without wrapping the cuff around the upper arm.
- When applying the cuff, ensure the arm protector flap is flat against arm to prevent injury to the skin.

About Blood Pressure

Blood pressure is the force exerted by blood against the walls of the arteries. Systolic pressure occurs when the heart contracts; diastolic pressure occurs when the heart expands. Blood pressure is measured in millimeters of mercury (mmHg).

Blood pressure is affected by many factors: age, weight, time of day, activity level, climate and altitude. Certain activities can alter one's blood pressure. Walking will likely raise blood pressure. Sleeping will likely decrease blood pressure. Not relaxing for several moments before measurements can influence readings as well. In addition to these factors, beverages containing caffeine or alcohol, certain medications, emotional stress and even tight-fitting clothes can make a difference in the readings.

About Hypertension

Hypertension (high blood pressure) is the diagnosis given when readings consistently rise above normal. Hypertension can lead to stroke, heart attack or other illness if left untreated. Referred to as a "silent killer" because it does not always produce symptoms that alert you to the problem, hypertension is treatable when diagnosed early.

In many individuals, hypertension can be controlled by altering lifestyle and minimizing stress, and by appropriate medication prescribed and monitored by your doctor. The American Heart Association recommends the following lifestyle suggestions to prevent or control hypertension:

- Do not smoke
- Exercise regularly
- Reduce salt and fat intake
- Have regular physical checkups
- Maintain proper weight

Blood pressure readings taken in a clinical setting might be elevated as a result of apprehension and anxiety. This is commonly called "white coat hypertension." Self-measurement at home supplements your doctor's readings and provides a more accurate, complete blood pressure history.

In addition, clinical studies have shown that the detection and treatment of hypertension is improved when patients both consult their physicians and monitor their own blood pressure at home.

Introduction

About Irregular Heartbeat

An irregular heartbeat is defined as a heartbeat that varies from the average of all heartbeats. When the monitor detects an irregular rhythm during the measurements, the IHB / AFib indicator will appear on the display with the measurement values.



Note: We recommend contacting your physician if you see this  IHB/AFib indicator frequently

What is AFib?

The heart contracts due to electrical signals occurring in heart and sends blood through the body. Atrial fibrillation (AFib) occurs when the electrical signal in the atrium becomes confused and leads to disturbances in the pulse interval. AFib can cause blood to stagnate in the heart, which can easily create clots of blood, a cause of stroke and heart attack.

Pressure Bar Indicator

The indicator monitors the progress of pressure during measurement.

Inflation in Progress	Inflation Complete	Deflation/ Measurement in Progress
		

This blood pressure monitor is designed to detect the pulse and to inflate the cuff to a systolic pressure level automatically.

Tips for Taking Blood Pressure

Before Your Measurement:

- For 30 minutes prior to taking your blood pressure
 - » Do not exercise
 - » Do not drink coffee, caffeinated soda or alcohol
 - » Do not smoke
- Sit quietly for 5-10 minutes
- Ensure you are using the right size cuff

During Your Measurement:

- Do not talk
- Sit with your back straight and supported
- Uncross your legs and place feet flat on floor
- Rest arm on a table so that the cuff is at heart level
- Measure two times a day, in the morning and evening

Notes for Accurate Measurement

- This monitor bases its measurements on the heartbeat. If you have a very weak or irregular heartbeat, the monitor may have difficulty determining your blood pressure.
- Should the monitor detect a condition that is abnormal, it will stop the measurement and display an error symbol. Refer to the section "Symbols" for the description of the symbols.
- If you have emotional stress, the measurement will reflect this stress as a higher (or lower) than normal blood pressure reading and the pulse reading will usually be faster than normal.
- An individual's blood pressure varies constantly, depending on what you are doing and what you have eaten. What you drink can have a very strong and rapid effect on your blood pressure.

Introduction

AHA Classification

Each segment of the bar indicator corresponds to the AHA blood pressure classification.

Example

Hypertension Stage 2	Hypertension Stage 1	Elevated
		

AHA Classification Indicator


- ← Hypertension Stage 2
- ← Hypertension Stage 1
- ← Elevated
- ← Normal

■: The indicator displays a segment, based on the current data, corresponding to the AHA classification.


Blood Pressure Classification	Systolic mm Hg		Diastolic mm Hg
Normal	< 120	and	< 80
Elevated	120 – 129	and	< 80
Stage 1 Hypertension (High Blood Pressure)	130 – 139	or	80 – 89
Stage 2 Hypertension (High Blood Pressure)	≥ 140	or	≥ 90
Severe Hypertension (If you don't have symptoms*, call your health care professional)	> 180	and/or	> 90
Hypertensive Emergency (Consult your doctor immediately)	> 180	and/or	> 120

*symptoms: chest pain, shortness of breath, back pain, numbness, weakness, change in vision, or difficulty speaking

Getting to Know Your Monitor

Display

Getting to Know Your Monitor

Symbols

Symbols that appear on the display

Symbols	Function/Meaning	Recommended Action
	Appears while measurement is in progress. It blinks when the pulse is detected.	Measurement is in progress. Remain as still as possible.
"♡"	Irregular Heartbeat symbol (I.H.B. / AFib) Appears when an irregular heartbeat is detected. It may light when a very slight vibration like shivering or shaking is detected.	_____
	TriCheck mode Takes three consecutive measurements automatically and displays the average values of the three measurements.	_____
	Appears when a body or arm movement is detected.	The reading may yield an incorrect value. Take another measurement. Remain still during measurement.
	Appears during measurement when the cuff is applied loosely.	The reading may yield an incorrect value. Apply the cuff correctly, and take another measurement.
	Previous measurements stored in memory	_____
	Average data	_____

Getting to Know Your Monitor

Symbols	Function/Meaning	Recommended Action
	The battery power indicator during measurement.	_____
	The battery is low when it blinks.	Replace all batteries with new ones when the mark blinks.
P _r	Pairing in progress.	_____
	The device is connecting to the Bluetooth® devices.	_____
SYS	Systolic blood pressure in mmHg	_____
DIA	Diastolic blood pressure in mmHg	_____
PUL/ min	Pulse per minute	_____
AM	Data taken in the AM	_____
PM	Data taken in the PM	_____
	Pressure settings Indicates the pressure value as set by the user.	_____

Getting to Know Your Monitor

Symbols that are printed on the device case

Symbols	Function/Meaning
START STOP	Standby and Turn the device on.
	Battery installation guide
	Direct current
	Serial number
	Global trade item number
	Lot number
	Refer to instruction manual/booklet
	Polarity of DC jack
IP	International Protection Symbol
	Keep dry
	Bluetooth® address


Symbols that are printed on the AC adapter

Symbols	Function/Meaning
	Alternating current
	Direct current
	Polarity of DC jack
	Polarity of DC jack
	Class II device
	For indoor use only
	UL/cUL certification device label
	Energy efficiency level VI certification device label
IP	International protection symbol

Getting to Know Your Monitor

Symbols that are printed on the cuff


Symbols	Function/Meaning	Recommended Action
●	Artery position mark	Set the ● mark on the artery of the upper arm or in line with the ring finger on the inside of the arm.
▲	Index	_____
REF	Catalog number	_____
LOT	Lot number	_____
	Type BF: applied part	_____
W	Proper fit range for the Wide Range cuff.	_____
L	Proper fit range for the Large cuff.	_____
M	Proper fit range for the Medium cuff.	_____
S	Proper fit range for the Small cuff.	_____
S	Range to use the Small cuff. Under range printed on the Medium cuff and Wide Range cuff.	Use the Small cuff
M	Range to use the Medium cuff. Under range printed on the Large cuff and over range printed on the Small cuff.	Use the Medium cuff
L	Range to use the Large cuff. Over range printed on the Medium cuff and Wide Range cuff.	Use the Large cuff

Initial Setup

Installing/Changing the Batteries

1. Remove the battery cover.
2. Remove the used batteries and insert new batteries into the battery compartment as shown, taking care that the polarities (+ and -) are correct. Use only R6P, LR6 or AA batteries.
3. Attach the battery cover.


CAUTION

- Insert the batteries as shown in the battery compartment. If installed incorrectly, the device will not work.
- When  (LOW BATTERY mark) blinks on the display and the device announces that the battery needs to be replaced, replace all batteries with new ones. Do not mix old and new batteries. It may shorten the battery life, or cause the device to malfunction. Replace the batteries two seconds or more after the device turns off. If  (LOW BATTERY mark) appears even after the batteries are replaced, make a blood pressure measurement. The device may then recognize the new batteries.
-  (LOW BATTERY mark) does not appear when the batteries are drained.
- The battery life varies with the ambient temperature and may be shorter at low temperatures.
- Use the specified batteries only. The batteries provided with the device are for testing monitor performance and may have a limited life.
- Remove the batteries if the device is not to be used for a long time. The batteries may leak and cause a malfunction.

Connecting the AC adapter to the monitor (purchased separately)

Insert the AC adapter plug into the DC jack.


Next, connect the AC adapter to an electrical outlet.

- Use the specified AC adapter. (Refer to Technical Data)
- When disconnecting the AC adapter from the electrical outlet, grasp and pull the AC adapter body out of the outlet.
- When disconnecting the AC adapter plug from the blood pressure monitor, grasp and pull the AC adapter plug out of the monitor.

Connecting the Air Hose


Insert the air connector plug into the air socket firmly.

Adjusting the Clock

The UA-1040TCNBLE blood pressure monitor has a built-in clock that gives you the time and date and can tell you when a blood pressure reading stored in memory was taken. We recommend that you adjust the clock prior to use. To adjust the clock, follow these simple steps:

1. When the monitor is turned off, press the set button  located to the right of the display.
2. Use the arrow   buttons to set the year.
3. After reaching the desired year, press the set button  again.
4. Repeat the process for the month, day and hour.
5. The next setting will be for pressure.

Initial Setup

Adjusting the Pressure

1. After setting the year, month, day and hour, the next mode is the pressure setting.
2. Default is  , which stands for automatic.
3. For even more rapid readings, you can adjust the default to about 30 mmHg or more above your expected systolic pressure.
4. Use the arrow ▲ ▼ buttons to set the pressure.
5. After reaching the desired pressure, press the set button  again.
6. The next setting will be for TriCheck Mode.

TriCheck Mode

The UA-1040TCNABLE blood pressure monitor has an optional feature to automatically take three consecutive readings one minute apart and then provide you with the average. This complies with the American Heart Association's recommended measuring method.

1. After setting the year, month, day, hour and pressure, the next mode is the  TriCheck Mode.
2. Default is off.
3. To turn on, use the arrow ▲ ▼ buttons.
4. To turn off again, use the arrow ▲ ▼ buttons.
5. After reaching the desired setting, press the set  button again.
6. This will end your settings session.
7. To adjust any values, push the set  button and begin again.

Adjusting the Speaker Volume

Note: Pressing  will toggle between volume setting → language setting → instruction setting → original state (standby) .

1. Press the  button to enter the volume setting mode.
2. There will be a confirmation tone and the the current volume will be displayed by the pressure bar (volume off, 1, 2, 3, 4, 5, 6).
If it is off , the LCD will display “ oFF ”.
3. Change the volume value with the ▲ ▼ buttons.

Selecting the Language

1. Press  to enter the language setting mode. L1, L2 or L3 will be displayed.
L1: English
L2: French
L3: Spanish
2. Change the language with the ▲ ▼ buttons.

Instructions On / Off

1. Press  to enter the instruction setting mode.
2. Turn the audio ON and OFF with the ▲ ▼ buttons.

Initial Setup

Pairing with a Mobile Device

Wireless Function

- In the unlikely event that this monitor causes radio wave interference to a different wireless station, change the location where this monitor is used or stop use immediately.
- The communication distance between this monitor and the mobile device is about 10 m. This distance is reduced by the conditions in the surrounding environment, so be sure to check that the distance is short enough for a connection to be made after measurement is complete.
- Be sure to use in a location where visibility between the two devices that you want to connect is good. The connection distance is reduced by the structure of buildings or other obstructions. In particular, connection may be impossible when devices are used on either side of reinforced concrete.
- Do not use *Bluetooth*® connection in the range of a wireless LAN or other wireless devices, near devices that emit radio waves such as microwaves, in locations where there are many obstructions, or in other locations where signal strength is weak. Doing so may result in frequent loss of connection, very slow communication speeds and errors.
- In this case, switch off the power supply to the device that is not being used or use the monitor in a different location.
- If the monitor does not connect normally when used near a wireless station or broadcast station, use the monitor in a different location.
- A&D Medical cannot accept liability for any damages incurred due to impaired operation or data loss, etc that occur through the use of this product.
- This product is not guaranteed to connect to all *Bluetooth*® compatible devices.

Bluetooth® Transmission

Bluetooth® devices carry the Bluetooth® logo mark.

To connect with your mobile device – download and install “A&D Heart Track” app.

Follow the instructions in the app to connect.


Cautions for Pairing

- Only one device can be paired with this monitor at one time. If the mobile device cannot receive measurement data, try pairing again.
- The monitor is capable of registering 10 mobile devices.
- In case an 11th mobile device is registered, the monitor will delete the oldest mobile device.

Pairing Procedure

Open the A&D Heart Track app, click on Add Bluetooth Devices and follow instructions.

1. Turn on *Bluetooth*® settings on your mobile device. Blinking → 
2. Press and hold the START/STOP button until “Pr” is displayed, and then release the button. The monitor will be in a state that can be found by the mobile device for about one minute.
3. If “Err 10” is displayed or pairing is failed, remove the batteries and try steps 1-2 again.
4. Accept the pairing request on the A&D Heart Track app.

Bluetooth
Pairing Request
“UA-1040TCNBLE”
would like to pair with your
iPhone.

Cancel

Pair

Initial Setup

Transmitting Temporarily Stored Data

In cases when the mobile device cannot receive measurement data, the measurement data is temporarily stored in the monitor memory. The data stored in the memory is transmitted the next time a connection is successfully made to the mobile device. A total of 90 sets of measurement data can be stored. When the amount of data exceeds 90, the oldest data is deleted and the new data is stored.

Time

This monitor has a built-in clock. The date and time that a measurement was taken is included in the measurement data. The built-in clock is automatically adjusted by syncing the clock of a mobile device. Sync of the time is done in the timing *Bluetooth*® icon lights up, in the pairing process.

Selecting the Correct Cuff Size

Using the correct cuff size is important for an accurate reading. If the cuff is not the proper size, the reading may yield an incorrect blood pressure value.

- The arm size is printed on each cuff.
- The index ▲ and proper fit range, on the cuff, tell you if you are applying the correct cuff.
- If the index ▲ points outside of the range, contact your local dealer to purchase a replacement cuff. The arm cuff is a consumable. If it becomes worn, purchase a new one.

Catalog Number	Recommended Cuff Size	Arm Size
UA-289A	Small Cuff	16-24 cm (6.3-9.4")
UA-290A	Medium Cuff	23-37 cm (9.0-14.6")
UA-291A	Large Cuff	31-45 cm (12.2-17.7")
UA-420A	Wide Range Cuff	22-42 cm (8.6-16.5")
UA-372A	Preformed Wide Range Cuff	22-42 cm (8.6-16.5")

Arm size: The circumference at the biceps.



Taking Your Blood Pressure

Applying the arm cuff

1. Insert arm (preferably the left arm) into cuff as shown in the figure to the right.
2. Wrap the cuff around the upper arm, about 1-2 cm (0.5") above the inside of the elbow, as shown in the figure to the right.
3. Place the cuff directly against the skin, as clothing may cause a faint pulse and result in a measurement error. Constriction of the upper arm, caused by rolling up a shirtsleeve, may prevent accurate readings.
4. Confirm the arm protector is flat against the arm as shown in the figure below.


Do not roll up shirtsleeve tightly

5. Confirm that the index ▲ points within the proper fit range.


Note: During measurement, it is normal for the cuff to feel very tight. (Do not be alarmed).

Taking Your Blood Pressure

Removing the Arm Cuff

1. Wait for the cuff to deflate.
2. Unfasten the Velcro.

3. Remove the cuff.

Taking Your Blood Pressure

Normal Measurement

1. Press the START/STOP button. All of the display segments are displayed. Zero (0) is displayed blinking briefly. The display changes, as indicated in the figure at the right, as the measurement begins. The cuff starts to inflate. It is normal for the cuff to feel very tight. A pressure bar indicator is displayed, on the left edge of the display, during the inflation.

Note: If you wish to stop inflation at any time, press the START/STOP button again.


2. When inflation is complete, deflation starts automatically and ♥ blinks, indicating that the measurement is in progress. Once the pulse is detected, the mark blinks with each pulse beat.

Note: If an appropriate pressure is not obtained, the device starts to inflate again automatically. To avoid re-inflation, see "Measurement with Set Pressure".

3. The display measurement is complete, the systolic and diastolic pressure readings and pulse rate are displayed. The device announces the measurement results. The cuff exhausts the remaining air and deflates completely.

4. Press the START/STOP button to turn the device off. After one minute of non-operation, the device will turn off automatically.


5. The reading is then stored into memory. See "About Memory" for more information.

Taking Your Blood Pressure

Measurement using TriCheck

1. Press the START/STOP button.
All of the display segments are displayed. Zero (0) is displayed blinking briefly and the first measurement starts.
2. When the measurement is complete, the systolic and diastolic pressure readings and pulse rate are displayed, then a one-minute countdown starts for the second measurement.
3. After one minute, the second measurement starts.
4. When the measurement is complete, the average readings of the first and second measurements are displayed, then a one-minute countdown starts for the third measurement.
5. After one minute, the third measurement starts.
6. When the measurement is complete, the average readings of the three measurements are displayed and stored in memory.

Notes: During the measurement, the TriCheck  3 symbol is displayed.

To cancel the measurement, press the START/STOP button. No data is stored in memory.

When the cuff fit error symbol appears during the first measurement, cancel the measurement, reapply the cuff correctly and start a new measurement.


After measurement, the average data of the three measurements is stored in memory. No data will be stored, when canceling the measurements, or before completion of the third measurement.

Taking Your Blood Pressure

Measurement with Set Pressure

The UA-1040TCNBLE is designed to detect the pulse and to inflate the cuff to a systolic pressure level automatically. If re-inflation occurs repeatedly, use the following methods. During the blood pressure measurement, re-inflation may occur. A fixed pressure value can be set to avoid re-inflation.

1. At step 5 of the clock adjustment procedure (see page E-17), press the set button to go to the pressure setting mode. The current setting blinks.
2. Press the ▲ or ▼ button to select a pressure value about 30mmHg or more above your expected systolic pressure from the following.

AUTO	: Automatic pressurization (default value)
180	: Pressure value of 180mmHg (fixed)
210	: Pressure value of 210mmHg (fixed)
240	: Pressure value of 240mmHg (fixed)

3. Press the set button to go to the TriCheck setting mode.


Press the **START/STOP** button to turn the device off. After three minutes of non-operation, the device will turn off automatically. The next measurement will be performed with the new pressure value.

Taking Your Blood Pressure

Measurement with desired Systolic Pressure

Use this method when re-inflation occurs repeatedly even if the pressure value is set to 240 in the procedure above or when the results are not displayed even if the pressure decreases to 20mmHg or less.

1. Place the cuff on the arm (preferably the left arm).
2. Press the START/STOP button.
3. While the zero blinks, press and hold the START/STOP button until a number about 30 to 40 mmHg higher than your expected systolic pressure appears.
4. When the desired number is reached, release the START/STOP button to start measurement. Continue to measure your blood pressure as described on the E-26.

Taking Your Blood Pressure

Note: This device stores the last 90 measurements in memory. The device announces the memory data values as they are displayed.

1. When the monitor is turned off, press the ▼ or ▲ button. The average of all measurements is displayed. (If no data, "0" is displayed. Press the ▼, ▲ or START/STOP button to turn device off.)
2. Each time the ▼ button (or the ▲ button to display the data in the reverse order) is pressed, the memory data is displayed as follows.

Average of all AM (morning) measurements taken between 4:00 am and 9:59 am. (In the example, 10 measurements. If no data, "--" is displayed.)

Average of all PM (evening) measurements taken between 6:00 PM and 1:59 AM. (In the example, 9 measurements. If no data, "--" is displayed.)


Most recent data (See No. 35). Three seconds after the data number display, the measurement data is displayed.

Last data (No.1). Three seconds after the data number display, the measurement data is displayed.

3. After the last data is displayed, press the ▼ button to return the average display of all measurements.
4. Press the START/STOP button to turn the device off. After one minute of non-operation, the device will turn off automatically.

Taking Your Blood Pressure

Deleting all Data Stored in Memory

Press both the ▲ and ▼ buttons at the same time. The **(M)** mark and the battery indicator appear. Press and hold both the ▲ and ▼ buttons until the illuminated **(M)** starts blinking to delete all data stored in memory.

Maintenance

Do not open the device. It uses delicate electrical components and an intricate air unit that could be damaged. If you cannot fix the problem using the troubleshooting instructions, contact the authorized dealer in your area or our customer service department. A&D customer service can provide technical assistance and spare parts.

Troubleshooting

Problem	Probable Cause	Corrective Action
Nothing appears on the display, even when the power is turned on.	Battery is drained.	Replace all batteries with new ones.
	Battery terminals are not in the correct position.	Reinstall the batteries with negative and positive terminals matching those indicated on the battery compartment.
The cuff does not inflate.	Battery voltage is too low.  (LOW BATTERY mark) blinks. If the batteries are drained completely, the mark does not appear.	Replace all batteries with new ones.
The device does not measure. Readings are too high or too low.	The cuff is not applied properly.	Apply the cuff correctly.
	Readings are too high or too low.	Make sure you remain very still and quiet during measurement.
	The cuff position is not correct.	Sit comfortably and still. Place your arm on a table with your palm facing upward and the cuff at the same level as your heart.
	_____	If you have a very weak or irregular heart beat, the device may have difficulty in determining your blood pressure.
Other	The value is different from that measured at a clinic or doctor's office.	Refer to "Tips for Taking Blood Pressure".
	_____	Remove the batteries. Place them back properly and take another measurement.

Note: If the actions described above do not solve the problem, contact the dealer. Do not attempt to open or repair this product, as any attempt to do so will make your warranty invalid.

Troubleshooting

Symbols	Function/Meaning	Recommended Action
Err	Unstable blood pressure due to movement during measurement.	Take another measurement. Remain very still during measurement.
	The systolic and diastolic values are within 10mmHg of each other.	Apply the cuff correctly, and take another measurement.
	The pressure value did not increase during the inflation.	
Err CUFF	The cuff is not applied correctly.	
E	PUL DISPLAY ERROR The pulse is not detected correctly.	
Err E	Blood pressure monitor internal error	Remove the batteries and press the START/STOP button, and then install the batteries again. If the error still appears, contact the dealer.
Err 9		
Err 10	Communication has failed.	Press the START/STOP button or remove and re-insert the batteries. Try communication again. If the error still appears, contact the dealer.

Optional Accessories

Accessory AC adapter

The adapter is to connect the device to a power source at home. Please contact your local A&D dealer for purchasing. The AC adapter is required to be inspected or replaced periodically.

Accessories (sold separately)

Part Number	Description	Specification
UA-289A	Small Cuff	16-24 cm (6.3-9.4")
UA-290A	Medium Cuff	23-37 cm (9.0-14.6")
UA-291A	Large Cuff	31-45 cm (12.2-17.7")
UA-420A	Wide Range Cuff	22-42 cm (8.6-16.5")
UA-372A	Preformed Wide Range Cuff	22-42 cm (8.6-16.5")
<i>Arm size: The circumference at the biceps.</i>		
TB-233	AC Adapter	Input: 120V ~ 60Hz 0.15A Output: 6V  500mA            IP21

Arm size: The circumference at the biceps.

Note: Specifications are subject to change without prior notice. IP classification is the degrees of protection provided by enclosures in accordance with IEC 60529. This device is protected against solid foreign objects of 12 mm diameter and greater such as a finger. This device is not protected against water.

LIMITED WARRANTY

A&D Medical

For purchasers within the US and Canada only:

Product	Consumer Warranty Term
Monitor	5 year
Cuff	2 year

For purchases made outside the United States or Canada, please contact the local distributor or retailer.

Limited Warranty:

A&D Medical ("A&D") warrants to the first purchaser ("You") that the A&D product You purchased (the "Product") will be free from defects in material, workmanship and design for the applicable Warranty Term stated above from the date You purchased the Product under normal use. This Limited Warranty is personal to You and is not transferable. If the Product is defective, then You return the Product to A&D in accordance with the procedure set forth below. A&D's warranty obligation is limited to the repair or replacement, at A&D's option, of the defective Product that has been returned by You within the warranty period. Such repair or replacement will be at no charge to You. The repaired or replacement Product is warranted here-under for the longer of the remainder of the original warranty period or 90 days from the date of shipment of the repaired or replacement Product.

To obtain a warranty service, please contact us in **US and Canada at 1-844-682-2322** for return address, shipping and handling fee, and other instructions for processing warranty. Please ensure you have satisfactory proof of the date of Your purchase and a description of the defect. Returns will not be accepted unless a Return Material Authorization (RMA) Number has been issued from A&D Customer Service Representative.

Warranty

This Limited Warranty does not cover, and A&D will not be liable for (i) any shipment damage, (ii) any damage or defect due to misuse, abuse, failure to use reasonable care, failure to follow written instructions enclosed with the Product, accident, subjecting the Product conditions, or modification, alteration or repair by anyone other than A&D or persons authorized by A&D, or (iii) expendable or consumable components.

THIS LIMITED WARRANTY IS THE ONLY WARRANTY PROVIDED BY A&D; THERE ARE NO OTHER EXPRESS WARRANTIES. If A&D cannot reasonably repair or replace the Product, A&D will refund the amount You paid for the Product (not including taxes), less a reasonable charge for usage. To receive a refund you must have returned the Product and all associated materials to A&D.

The above remedy of repair, replacement or refund is your only and exclusive remedy. IN NO EVENT SHALL A&D BE LIABLE FOR ANY DAMAGES, INCLUDING, WITHOUT LIMITATION, LOST PROFITS, LOST INFORMATION OR REPLACEMENT COSTS, ARISING OUT OF YOUR USE OF OR INABILITY TO USE THE PRODUCT, INCLUDING, WITHOUT LIMITATION, ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, EVEN IF A&D HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. Some states do not allow the exclusion of incidental or consequential damages, so that the above exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that may vary from state to state.

No distributor, dealer or other party is authorized to make any warranty on behalf of A&D or to modify this warranty, or to assume for A&D any liability with respect to its products.

Technical Data

Model	UA-1040TCNBLE
Type	Oscillometric
Memory	90 readings
Measurement range	Pressure: 0 - 299mmHg Systolic pressure: 60 - 279mmHg Diastolic pressure: 40 - 200mmHg Pulse: 40 - 180 beats per minute
Static accuracy	Pressure: ± 3 mmHg Pulse: $\pm 5\%$
Rating range	DC 6V 3W
Power source	4 x 1.5 volt batteries (AA, LR6 or R6P) or 120 volt AC adapter (TB-233)
Classification	Internally powered ME equipment (Supplied by batteries) /Class II (Supplied by adapter) Continuous operation mode.
Clinical test	According to ISO81060-2 : 2020
EMD	IEC 60601-1-2: 2014+A1:2020
Wireless Communication	Bluetooth:Ver.5.1LE BLP
Frequency band	2402 MHz to 2480 MHz
Maximum RF output power	< 10 dBm
Modulation	GFSK
Supported Data	Systolic Pressure, Diastolic Pressure, Pulse Rate

Technical Data

Communication distance	About 10 m (This distance is reduced by the conditions in the surrounding environment)
Paired device	10 devices
Operating Environment	50°F to 104°F (10°C to 40°C) 15%RH to 85%RH / 800hPa to 1060hPa
Transport/Storage Environment	-4°F to 140°F (-20°C to 60°C) 10%RH to 95%RH / 700hPa to 1060hPa
Dimensions	5.5" (W) x 2.5" (H) x 4.6" (D) 140mm (W) x 64mm (H) x 118mm (D)
Weight	11.3 oz (320g) without batteries
Ingress protection	Device: IP20
Number of measurements	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries)With pressure value 180 mmHg, roomtemperature 23 °C
Applied part.	Cuff Type BF 
Useful life	Device: 5 years(when used six times a day) Cuff: 2 years (when used six times a day) AC adapter: 5 years (when used six times a day)

Note: Specifications are subject to change without prior notice. IP classification is the degrees of protection provided by enclosures in accordance with IEC 60529. This device is protected against solid foreign objects of 12 mm diameter and greater such as a finger. This device is not protected against water.

EMD Technical Data

Battery-operated or AC Adapter-operated Blood Pressure Monitor

Medical Electrical Equipment needs special precautions regarding EMD and needs to be installed and put into service according to the EMD information provided in the following.

Portable and mobile RF communication equipment (e.g. cell phones) can affect Medical Electrical Equipment.

The use of accessories and cables other than those specified may result in increased emissions or decreased immunity of the unit.

Table 1 - EMISSION Limits -

Phenomenon	Compliance
Conducted and radiated RF EMISSION CISPR 11	Group 1, Class B
Harmonic distortion IEC 61000-3-2	Class A
Voltage fluctuations and flicker IEC 61000-3-3	Compliance

Table 2 - IMMUNITY TEST LEVELS : Enclosure Port -

Phenomenon	IMMUNITY TEST LEVELS
Electrostatic discharge IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Radiated RF EM fields IEC 61000-4-3	10 V/m 80 MHz - 2.7 GHz 80 % AM at 1 kHz
Proximity fields from RF wireless communications equipment IEC 61000-4-3	See table 4
Rated power frequency magnetic fields IEC 61000-4-8	30 A/m 50 Hz or 60 Hz
proximity magnetic fields IEC 61000-4-39	See table 5

Technical Data

Table 3 - IMMUNITY TEST LEVELS : Input a.c. power Port -

Phenomenon	IMMUNITY TEST LEVELS
Electrical fast transients / bursts IEC 61000-4-4	± 2 kV 100 kHz repetition frequency
Surges Line-to-line IEC 61000-4-5	± 0.5 kV, ± 1 kV
Conducted disturbances induced by RF fields IEC 61000-4-6	3 V 0.15 MHz - 80 MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz
Voltage dips IEC 61000-4-11	0 % U_T ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°
	0 % U_T ; 1 cycle and 70 % U_T ; 25/30 cycle Single phase: at 0°
Voltage interruption IEC 61000-4-11	0% U_T ; 250/300 cycle
NOTE U_T is the AC mains voltage prior to application of the test level.	

Table 4 - Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment -

Test frequency (MHz)	Band (MHz)	Service	Modulation	IMMUNITY TEST LEVEL (V/m)
385	380 - 390	TETRA 400	Pulse modulation 18 Hz	27
450	430 - 470	GMRS 460 FRS 460	FM ± 5 kHz deviation 1 kHz sine	28
710	704 - 787	LTE Band 13, 17	Pulse modulation 217 Hz	9
745				
780				
810	800 - 960	GSM 800/900 TETRA 800 iDEN 820 CDMA 850 LTE Band 5	Pulse modulation 18 Hz	28
870				
930				

Technical Data

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)
1720	1700 - 1990	GSM 1800 CDMA 1900 GSM 1900 DECT LTE Band 1, 3, 4, 25 UMTS	Pulse modulation 217 Hz	28
1845				
1970				
2450	2400 - 2570	Bluetooth WLAN 802.11 b/g/n RFID 2450 LTE Band 7	Pulse modulation 217 Hz	28
5240	5100 - 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9
5500				
5785				

Table 5 - Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields-

Test frequency	Modulation	IMMUNITY TEST LEVEL(A/m)
30 kHz	CW	8
134.2 kHz	Pulse modulation 2.1 kHz	65
13.56 MHz	Pulse modulation 50 kHz	7.5

Technical Data

FCC CAUTION

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This transmitter must not be co-located or operated in conjunction with any other antenna or transmitter.

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment and meets the FCC radio frequency (RF) Exposure Guidelines in Supplement C to OET65. This equipment has very low levels of RF energy that it deemed to comply without maximum permissive exposure evaluation (MPE). But it is desirable that it should be installed and operated keeping the radiator at least 20 cm or more away from person's body (excluding extremities: hands, wrists, feet and ankles).

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

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© 2024 All trademarks are property of their rightful owners
A&D Company, Limited
A&D Engineering, Inc.