MULTI USER

Blood Pressure Monitor

Premium Upper Arm

UA-770CNBLE

Instruction Manual

ORIGINAL

Table of Contents

1. Introduction	E-3
Intended Use / FCC Compliance / Precautions	
About Blood Pressure / About Hypertension	
About Irregular Heartbeat	
What is AFib? / Pressure Bar Indicator	
Tips for Taking Blood Pressure	
AHA Classification	
2. Getting to Know Your Monitor	E-10
Parts Identification / Display	
Symbols	
3. Initial Setup	E-14
Installing and Changing Batteries	
Connecting Air Hose / Connecting and Disconnecting AC Adapter	
Adjusting Clock	
Pairing with a Mobile Device	
Selecting the Correct Cuff Size	
4. Taking Your Blood Pressure	E-24
Applying and Removing the Cuff	
Normal Measurement	
Measurement with Desired Systolic Pressure	
Reviewing and Deleting Memory	
5. Maintenance	E-29
6. Troubleshooting	E-30
7. Optional Accessories	E-32
8. Warranty	E-33
9. Technical Data	E-35

Congratulations on purchasing a state-of-the-art A&D blood pressure monitor. Designed for ease of use and accuracy, this monitor will facilitate your daily blood pressure regimen. We recommend that you read through this manual carefully before using the monitor for the first time.

Intended Use

- The monitor is designed for use on adults only. Do not use on newborns or infants.
- Environment for use: the monitor is for use in the home healthcare environment.
- This monitor is designed to measure blood pressure and pulse rate of people for diagnosis.

FCC Compliance Information

This monitor complies with Part 15 of FCC rules. Operation is subject to the following two conditions:

1. This monitor may not cause harmful interference, and
2. This monitor must accept any interference received, including interference that may cause undesired operations.

Precautions

- This is a monitor for use by lay person. Please consult your healthcare provider with any questions or concerns you may have regarding your condition.
- If any serious incident occurs in relation to this monitor, report to its manufacturer and the competent authority in your country.
- This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions: (1) This device may not cause interference. (2) This device must accept any interference, including interference that may cause undesired operation of the device.
- Precision components are used in the construction of this monitor. Extremes in temperature, humidity, direct sunlight, shock or dust should be avoided.

Introduction

- Clean the monitor and cuff with a dry, soft cloth or a cloth dampened with water and a neutral detergent. Never use alcohol, benzene, thinner or other harsh chemicals to clean the monitor or cuff.
- Avoid tightly folding the cuff or storing the hose tightly twisted for long periods, as such treatment may shorten the life of the components.
- The monitor and cuff are not water resistant. Prevent rain, sweat and water from soiling the monitor and cuff.
- Measurements may be distorted if the monitor is used close to televisions, microwave ovens, cellular telephones, X-ray or other devices with strong electrical fields.
- When using the monitor, confirm that the monitor is clean.
- Used equipment, parts and batteries are not treated as ordinary household waste, and must be disposed of according to the applicable local regulations.
- When the AC adapter is used, make sure that the AC adapter can be readily removed from the electrical outlet when necessary.
- Do not modify the monitor. It may cause accidents or damage to the monitor.
- To measure blood pressure, the arm must be squeezed by the cuff hard enough to temporarily stop blood flow through the artery. This may cause pain, numbness or a temporary red mark to the arm. This condition will appear especially when measurement is repeated successively. Any pain, numbness, or red marks will disappear with time.
- People who have a severe circulatory deficit in the arm must consult a doctor before using the monitor, to avoid medical problems.
- Do not self-diagnose the measurement results and start treatment by yourself. Always consult your doctor for evaluation of the results and treatment.
- Do not apply the cuff on an arm with an unhealed wound.
- Do not apply the cuff on an arm receiving an intravenous drip or blood transfusion. It may cause injury or accidents.
- Do not use the monitor where flammable gases such as anesthetic gases are present. It may cause an explosion.
- Do not use the monitor in highly concentrated oxygen environments, such as a high-pressure oxygen chamber or an oxygen tent. It may cause a fire or explosion.

- Do not excessively bend or squeeze the air hose.
- Do not twist the air hose during measurement. This may cause injury due to continuous cuff pressure.
- Measuring blood pressure too frequently may cause harm due to blood flow interference. Check that operation of the monitor does not result in prolonged impairment of blood circulation, when using the monitor repeatedly.
- If you have had a mastectomy or lymph node clearance, please consult a doctor before using the monitor.
- Take care to avoid accidental strangulation with babies or infants by the hose and cable.
- Wireless communication devices, such as home networking devices, mobile phones, cordless phones and their base stations, walkie-talkies can affect this blood pressure monitor. Therefore, a minimum distance of 30 cm (12") should be kept from such devices.
- Do not let children use the monitor by themselves and do not use the monitor in a place within the reach of infants. It may cause accidents or damage.
- There are small parts that may cause a choking hazard if swallowed by mistake by infants.
- Unplug the AC adapter when not in use during the measurement.
- Use of accessories not detailed in this manual may compromise safety.
- Should the battery short-circuit, it may become hot and potentially cause burns.
- Allow the monitor to adapt to the surrounding environment before use (about one hour).
- Clinical testing has not been conducted on newborn infants and pregnant woman.
- Do not touch the batteries, the DC jack, and the patient at the same time. That may result in electrical shock.
- Do not inflate without wrapping the cuff around the upper arm.
- Do not apply the cuff on an arm with another medical device is attached. The equipment may not function properly.
- When applying the cuff, ensure the arm protector flap is flat against arm to prevent injury to the skin.

Introduction

About Blood Pressure

Blood pressure is the force exerted by blood against the walls of the arteries. Systolic pressure occurs when the heart contracts; diastolic pressure occurs when the heart expands. Blood pressure is measured in millimeters of mercury (mmHg).

Blood pressure is affected by many factors: age, weight, time of day, activity level, climate and altitude. Certain activities can alter one's blood pressure. Walking will likely raise blood pressure. Sleeping will likely decrease blood pressure. Not relaxing for several moments before measurements can influence readings as well. In addition to these factors, beverages containing caffeine or alcohol, certain medications, emotional stress and even tight-fitting clothes can make a difference in the readings.

About Hypertension

Hypertension (high blood pressure) is the diagnosis given when readings consistently rise above normal. Hypertension can lead to stroke, heart attack or other illness if left untreated. Referred to as a "silent killer" because it does not always produce symptoms that alert you to the problem, hypertension is treatable when diagnosed early.

In many individuals, hypertension can be controlled by altering lifestyle and minimizing stress, and by appropriate medication prescribed and monitored by your doctor. The American Heart Association recommends the following lifestyle suggestions to prevent or control hypertension:

- Do not smoke
- Exercise regularly
- Reduce salt and fat intake
- Have regular physical checkups
- Maintain proper weight

Blood pressure readings taken in a clinical setting might be elevated as a result of apprehension and anxiety. This is commonly called "white coat hypertension." Self-measurement at home supplements your doctor's readings and provides a more accurate, complete blood pressure history.

In addition, clinical studies have shown that the detection and treatment of hypertension is improved when patients both consult their physicians and monitor their own blood pressure at home.

About Irregular Heartbeat

An irregular heartbeat is defined as a heartbeat that varies from the average of all heartbeats. When the monitor detects an irregular rhythm during the measurements, the IHB / AFib indicator will appear on the display with the measurement values.


Note: We recommend contacting your physician if you see this  IHB/AFib indicator frequently

What is AFib?

The heart contracts due to electrical signals occurring in heart and sends blood through the body. Atrial fibrillation (AFib) occurs when the electrical signal in the atrium becomes confused and leads to disturbances in the pulse interval. AFib can cause blood to stagnate in the heart, which can easily create clots of blood, a cause of stroke and heart attack.

Pressure Bar Indicator

The indicator monitors the progress of pressure during measurement

Introduction

This blood pressure monitor is designed to detect the pulse and to inflate the cuff to a systolic pressure level automatically.

Tips for Taking Blood Pressure

Before Your Measurement:

- For 30 minutes prior to taking your blood pressure
 - » Do not exercise
 - » Do not drink coffee, caffeinated soda or alcohol
 - » Do not smoke
- Sit quietly for 5-10 minutes
- Ensure you are using the right size cuff

During Your Measurement:

- Do not talk
- Sit with your back straight and supported
- Uncross your legs and place feet flat on floor
- Rest arm on a table so that the cuff is at heart level
- Measure two times a day, in the morning and evening
- It is normal for the cuff to feel very tight. (Do not be alarmed).

Notes for Accurate Measurement

- This monitor bases its measurements on the heartbeat. If you have a very weak or irregular heartbeat, the monitor may have difficulty determining your blood pressure.
- Should the monitor detect a condition that is abnormal, it will stop the measurement and display an error symbol. Refer to the section "Symbols" for the description of the symbols.
- If you have emotional stress, the measurement will reflect this stress as a higher (or lower) than normal blood pressure reading and the pulse reading will usually be faster than normal.
- An individual's blood pressure varies constantly, depending on what you are doing and what you have eaten. What you drink can have a very strong and rapid effect on your blood pressure.

AHA Classification

Each segment of the bar indicator corresponds to the AHA blood pressure classification.

Example

Hypertension Stage 2	Hypertension Stage 1	Elevated
		

AHA Classification Indicator



-  ← Hypertension Stage 2
-  ← Hypertension Stage 1
-  ← Elevated
-  ← Normal

: The indicator displays a segment, based on the current data, corresponding to the AHA classification.


Blood Pressure Classification	Systolic mm Hg		Diastolic mm Hg
Normal	< 120	and	< 80
Elevated	120 – 129	and	< 80
Stage 1 Hypertension (High Blood Pressure)	130 – 139	or	80 – 89
Stage 2 Hypertension (High Blood Pressure)	≥ 140	or	≥ 90
Severe Hypertension (If you don't have symptoms*, call your health care professional)	> 180	and/or	> 90
Hypertensive Emergency (Consult your doctor immediately)	> 180	and/or	> 120

*symptoms: chest pain, shortness of breath, back pain, numbness, weakness, change in vision, or difficulty speaking

Getting to Know Your Monitor

Display

Getting to Know Your Monitor

Symbols

Symbols that appear on the display


Symbols	Function/Meaning	Recommended Action
	Appears while measurement is in progress. It blinks when the pulse is detected.	Measurement is in progress. Remain as still as possible.
	Irregular Heartbeat symbol (I.H.B. / AFib) Appears when an irregular heartbeat is detected. It may light when a very slight vibration like shivering or shaking is detected.	_____
	Appears when a body or arm movement is detected.	The reading may yield an incorrect value. Take another measurement. Remain still during measurement.
	Appears during measurement when the cuff is attached loosely.	The reading may yield an incorrect value. Apply the cuff correctly, and take another measurement.
	User	_____
	Previous measurements stored in memory.	_____
	Average data	_____

Getting to Know Your Monitor

Symbols	Function/Meaning	Recommended Action
	FULL BATTERY The battery power indicator during measurement.	_____
	LOW BATTERY The battery is low when it blinks.	Replace all batteries with new ones when the indicator blinks.
	Pairing in progress.	_____
	The device is connecting to the <i>Bluetooth</i> ® devices.	_____
AM	Data taken between 12:01am and 12:00pm	_____
PM	Data taken between 12:01pm and 12:00am	_____

Getting to Know Your Monitor

Symbols that are printed on the device case

Symbols	Function/Meaning
	Standby and Turn the device on.
SYS	Systolic blood pressure in mmHg
DIA	Diastolic blood pressure in mmHg
PUL/min	Pulse per minute
	Battery installation guide
	Direct current
	Serial number
	Lot number
	Global trade item number
	Bluetooth® address
	Refer to instruction manual/booklet
	Polarity of DC jack
IP	International Protection Symbol
	Keep dry

Getting to Know Your Monitor


Symbols that are printed on the AC adapter

Symbols	Function/Meaning
	Alternating current
	Direct current
	Polarity of DC jack
	Polarity of DC jack
	Class II device
	For indoor use only
	UL/cUL certification device label
	Energy efficiency level VI certification device label
IP	International protection symbol

Getting to Know Your Monitor

Symbols that are printed on the cuff


Symbols	Function/Meaning	Recommended Action
●	Artery position mark	Set the ● mark on the artery of the upper arm or in line with the ring finger on the inside of the arm.
▲	Index	_____
REF	Catalog number	_____
LOT	Lot number	_____
W	Proper fit range for the Wide Range cuff.	_____
L	Proper fit range for the Large cuff.	_____
M	Proper fit range for the Medium cuff.	_____
S	Proper fit range for the Small cuff.	_____
S	Range to use the Small cuff. Under range printed on the Medium cuff and Wide Range cuff.	Use the Small cuff
M	Range to use the Medium cuff. Under range printed on the Large cuff and over range printed on the Small cuff.	Use the Medium cuff
L	Range to use the Large cuff. Over range printed on the Medium cuff and Wide Range cuff.	Use the Large cuff
人	Type BF: applied part	

Initial Setup

Installing/Changing the Batteries

1. Remove the battery cover.
2. Remove the used batteries and insert new batteries into the battery compartment as shown, taking care that the polarities (+ and -) are correct. Use only R6P, LR6 or AA batteries.
3. Attach the battery cover.


CAUTION

- Insert the batteries as shown in the battery compartment. If installed incorrectly, the device will not work.
- When  (LOW BATTERY mark) blinks on the display, replace all batteries with new ones. Do not mix old and new batteries. It may shorten the battery life, or cause the device to malfunction.
- Replace the batteries two seconds or more after the device turns off.
- If  (LOW BATTERY mark) appears even after the batteries are replaced, make a blood pressure measurement. The device may then recognize the new batteries.
-  (LOW BATTERY mark) does not appear when the batteries are drained.
- The battery life varies with the ambient temperature and may be shorter at low temperatures.
- Use the specified batteries only.
- Remove the batteries if the device is not to be used for a long time. The batteries may leak and cause a malfunction.

Connecting the Air Hose

Insert the air connector plug into the air socket firmly.


Connecting the AC Adapter

(purchased separately)

Insert the AC adapter plug into the DC jack.

Next, connect the AC adapter to an electrical outlet.

- Use the specified AC adapter.
(Refer to Technical Data)
- When disconnecting the AC adapter from the electrical outlet, grasp and pull the AC adapter body out of the outlet.
- When disconnecting the AC adapter plug from the blood pressure monitor, grasp and pull the AC adapter plug out of the monitor.

Initial Setup


Adjusting the Built-in Clock

Adjust the clock prior to use.

1. When the monitor is turned off, press and hold the  button until the year starts blinking.
2. Select the year using the  or  button. Press the  button to set the current year and move to month/day selection. The date can be set anywhere between the years 2024 and 2079.
3. Select the month using the  or  button. Press the  button to set the current month and move to day selection.
4. Select the day using the  or  button. Press the  button to set the current day and move to hour/minute selection.
5. Select the hour using the  or  button. Press the  button to set the current hour and move to minute selection.
6. Select the minute using the  or  button. Press the START/STOP or  button to turn the device off.

- Holding down the arrow buttons will change the value continuously.
- Pressing the START/STOP button will turn the device off.

Note: After three minutes of non-operation, the device will turn off automatically. When the clock has not been set, the clock display indicates dashes.

Pairing with a Mobile Device

Wireless Function

- In the unlikely event that this monitor causes radio wave interference to a different wireless station, change the location where this monitor is used or stop use immediately.
- The communication distance between this monitor and the mobile device is about 10 m. This distance is reduced by the conditions in the surrounding environment, so be sure to check that the distance is short enough for a connection to be made after measurement is complete.
- Be sure to use in a location where visibility between the two devices that you want to connect is good. The connection distance is reduced by the structure of buildings or other obstructions. In particular, connection may be impossible when devices are used on either side of reinforced concrete.
- Do not use *Bluetooth*[®] connection in the range of a wireless LAN or other wireless devices, near devices that emit radio waves such as microwaves, in locations where there are many obstructions, or in other locations where signal strength is weak. Doing so may result in frequent loss of connection, very slow communication speeds and errors.
- In this case, switch off the power supply to the device that is not being used or use the monitor in a different location.
- If the monitor does not connect normally when used near a wireless station or broadcast station, use the monitor in a different location.
- A&D Medical cannot accept liability for any damages incurred due to impaired operation or data loss, etc that occur through the use of this product.
- This product is not guaranteed to connect to all *Bluetooth*[®] compatible devices.

Initial Setup

Bluetooth® Transmission Bluetooth®

Bluetooth® devices carry the Bluetooth® logo mark.

To connect with your mobile device – download and install A&D Heart Track app.


Follow the instructions in the app to connect.

Cautions for Pairing

- Only one device can be paired with this monitor at one time. If the mobile device cannot receive measurement data, try pairing again.
- The monitor is capable of registering 10 mobile devices.
- In case an 11th mobile device is registered, the monitor will delete the oldest mobile device.

Pairing Procedure

1. Turn on *Bluetooth®* settings on your mobile device.
2. Press and hold the START/STOP button until “Pr” is displayed, and then release the button. The monitor will be in a state that can be found by the mobile device for about one minute.
3. If “Err 10” is displayed or pairing is failed, remove the batteries and try steps 1-2 again.
4. Accept the pairing request on the A&D Heart Track app.

Bluetooth
Pairing Request
“UA-770CNBLE”
would like to pair with your
iPhone.

Cancel

Pair

Transmitting Temporarily Stored Data

In cases when the mobile device cannot receive measurement data, the measurement data is temporarily stored in the monitor memory. The data stored in the memory is transmitted the next time a connection is successfully made to the mobile device. A total of 60 sets of measurement data can be stored per user. When the amount of data exceeds 60, the oldest data is deleted and the new data is stored.

Time

This monitor has a built-in clock. The date and time that a measurement was taken is included in the measurement data. The built-in clock is automatically adjusted by syncing the clock of a mobile device. Sync of the time is done in the timing *Bluetooth*[®] icon lights up, in the pairing process.

Initial Setup

Selecting the Correct Cuff Size

Using the correct cuff size is important for an accurate reading. If the cuff is not the proper size, the reading may yield an incorrect blood pressure value.

- The arm size is printed on each cuff.
- The index ▲ and proper fit range, on the cuff, tell you if you are applying the correct cuff. (Refer to “Symbols that are printed on the cuff”)
- If the index ▲ points outside of the range, contact your local dealer to purchase a replacement cuff. The arm cuff is a consumable. If it becomes worn, purchase a new one.


Cuff Model	Cuff Size	Arm Size
UA-420A	Wide Range Cuff	22-42 cm (8.6-16.5")
UA-291A	Large Cuff	31-45 cm (12.2-17.7")
UA-290A	Medium Cuff	23-37 cm (9.0-14.6")
UA-289A	Small Cuff	16-24 cm (6.3-9.4")
UA-372A	Preformed Wide Range Cuff	22-42 cm (8.6-16.5")

Arm size: The circumference at the biceps.



Taking Your Blood Pressure

Applying the arm cuff

1. Insert arm (preferably the left arm) into cuff as shown in the figure to the right.
2. Wrap the cuff around the upper arm, about 1-2 cm (0.5") above the inside of the elbow, as shown in the figure to the right.
3. Place the cuff directly against the skin, as clothing may cause a faint pulse and result in a measurement error. Constriction of the upper arm, caused by rolling up a shirtsleeve, may prevent accurate readings.
4. Confirm the arm protector is flat against the arm as shown in the figure.


Do not roll up shirtsleeve tightly

5. Confirm that the index ▲ points within the proper fit range.


Note: During measurement, it is normal for the cuff to feel very tight. (Do not be alarmed).

Taking Your Blood Pressure

Removing the arm cuff

1. Wait for the cuff to deflate
2. Unfasten the Velcro

3. Remove the cuff

Taking your Blood Pressure

Normal Measurement

1. Press the  button to select a user from 1 to 4.

2. Place the cuff on the arm (preferably the left arm).

Sit quietly during measurement.

3. Press the START/STOP button.
All of the display segments are displayed.
Zero (0) is displayed blinking briefly.

The display changes, as indicated in the figure at the right, as the measurement begins. The cuff starts to inflate. It is normal for the cuff to feel very tight. A pressure bar indicator is displayed, on the left edge of the display, during the inflation.


Note: If you wish to stop inflation at any time, press the START/STOP button again.

4. When inflation is complete, deflation starts automatically and  (heart mark) blinks, indicating that the measurement is in progress. Once the pulse is detected, the mark blinks with each pulse beat.

Note: If an appropriate pressure is not obtained, the device starts to inflate again automatically.

To avoid re-inflation, see "Measurement with the Desired Systolic Pressure".

5. When the measurement is complete, the systolic and diastolic pressure readings and pulse rate are displayed. The cuff exhausts the remaining air and deflates completely.
6. Press the START/STOP button to turn the device off. After one minute of non-operation, the device will turn off automatically.


Taking your Blood Pressure

Measurement with the Desired Systolic Pressure

If re-inflation occurs repeatedly, or if your systolic pressure is expected to exceed 230 mmHg, use this procedure.

1. Place the cuff on the arm at heart level (preferably the left arm).
2. Press the START/STOP button.
3. While the zero blinks, press and hold the START/STOP button until a number about 30 to 40 mmHg higher than your expected systolic pressure appears.
4. When the desired number is reached, release the START/STOP button to start measurement. Then continue to measure your blood pressure as described on section "Normal Measurement".

Taking Your Blood Pressure

Reviewing Memory

Note: This device stores the last 60 measurements for each of the four user in memory.

1. When the monitor is turned off, press the ▲ or ▼ button.

The average of last 3 memory and the number of data are displayed.
(If no data, "0" is displayed. Press the ▲, ▼ or START/STOP button to turn the device off.)

2. Each time the ▲ button (or the ▼ button to display the data in the reverse order) is pressed, the memory data is displayed as follows.

Most recent data (No.n, in the example, No.35)


Three seconds after the data number display, the measurement data is displayed.

Last data (No.1)

Three seconds after the data number display, the measurement data is displayed.


3. After the last data is displayed, press the ▼ button to return the average display of the last 3 measurements.
4. Press the START/STOP button to turn the device off. After one minute of non-operation, the device will turn off automatically.

Taking Your Blood Pressure


Changing user for memory display

Press the  button while a memory is displayed. The user is changed, and the average of measurement values for that user is displayed.

Deleting All Data Stored in Memory

In standby, press both the ▲ and ▼ buttons. The  mark, battery indicator and user mark will appear. When you would like to delete the memory data of the currently displayed user, press and hold both the ▲ and ▼ buttons until the illuminated  mark starts blinking.

Do not open the device. It uses delicate electrical components and an intricate air unit that could be damaged. If you cannot fix the problem using the troubleshooting instructions, contact the authorized dealer in your area or our customer service department. A&D customer service will provide technical information, spare parts and units to authorized dealer.

Troubleshooting

Problem	Probable cause	Corrective action
Nothing appears on the display, even when the power is turned on.	Battery is drained.	Replace all batteries with new ones.
	Battery terminals are not in the correct position.	Reinstall the batteries with negative and positive terminals matching those indicated on the battery compartment.
The cuff does not inflate.	Battery voltage is too low.  (LOW BATTERY mark) blinks. If the batteries are drained completely, the mark does not appear.	Replace all batteries with new ones.
The device does not measure. Readings are too high or too low.	The cuff is not applied properly.	Apply the cuff correctly.
	Readings are too high or too low.	Make sure you remain very still and quiet during measurement.
	The cuff position is not correct.	Sit comfortably and still. Place your arm on a table with your palm facing upward and the cuff at the same level as your heart.
	_____	If you have a very weak or irregular heart beat, the device may have difficulty in determining your blood pressure.
Other	The value is different from that measured at a clinic or doctor's office.	Refer to "Tips for Taking Blood Pressure".
	_____	Remove the batteries. Place them back properly and take another measurement.

Note: If the actions described above do not solve the problem, contact the dealer. Do not attempt to open or repair this product, as any attempt to do so will make your warranty invalid.

Troubleshooting


Symbols	Function/Meaning	Recommended Action
Err	Unstable blood pressure due to movement during measurement.	Take another measurement. Remain very still during measurement.
	The systolic and diastolic values are within 10 mmHg of each other.	Apply the cuff correctly, and try the measurement again.
	The pressure value did not increase during the inflation.	
Err CUFF	The cuff is not applied correctly.	
E	Pulse display error. The pulse is not detected correctly.	
Err E	Blood pressure monitor internal error	Remove the batteries and press the START/STOP button, and then install the batteries again. If the error still appears, contact the dealer.
Err 9		
Err 10	Communication has failed.	Press the START/STOP button or remove and re-insert the batteries. Try communication again. If the error still appears, contact the dealer.

Optional Accessories

Accessory AC adapter

The adapter is to connect the device to a power source at home. Please contact your local A&D dealer for purchasing. The AC adapter is required to be inspected or replaced periodically.

Accessories (sold separately)

Part Number	Description	Specification
UA-420A	Wide Range Cuff	22-42 cm (8.6-16.5")
UA-291A	Large Cuff	31-45 cm (12.2-17.7")
UA-290A	Medium Cuff	23-37 cm (9.0-14.6")
UA-289A	Small Cuff	16-24 cm (6.3-9.4")
UA-372A	Preformed Wide Range Cuff	22-42 cm (8.6-16.5")
<i>Arm size: The circumference at the biceps.</i>		
TB-233	AC Adapter	Input: 120V ~ 60Hz 0.15A Output: 6V  500mA   IP21

Arm size: The circumference at the biceps.

LIMITED WARRANTY

A&D Medical

For purchasers within the US and Canada only:

Product	Consumer Warranty Term
Monitor	5 year
Cuff	2 year

For purchases made outside the United States or Canada, please contact the local distributor or retailer.

Limited Warranty:

A&D Medical ("A&D") warrants to the first purchaser ("You") that the A&D product You purchased (the "Product") will be free from defects in material, workmanship and design for the applicable Warranty Term stated above from the date You purchased the Product under normal use. This Limited Warranty is personal to You and is not transferable. If the Product is defective, then You return the Product to A&D in accordance with the procedure set forth below. A&D's warranty obligation is limited to the repair or replacement, at A&D's option, of the defective Product that has been returned by You within the warranty period. Such repair or replacement will be at no charge to You. The repaired or replacement Product is warranted here-under for the longer of the remainder of the original warranty period or 90 days from the date of shipment of the repaired or replacement Product.

To obtain a warranty service, please contact us in **US and Canada at 1-844-682-2322** for return address, shipping and handling fee, and other instructions for processing warranty. Please ensure you have satisfactory proof of the date of Your purchase and a description of the defect. Returns will not be accepted unless a Return Material Authorization (RMA) Number has been issued from A&D Customer Service Representative.

This Limited Warranty does not cover, and A&D will not be liable for (i) any shipment damage, (ii) any damage or defect due to misuse,

Warranty

abuse, failure to use reasonable care, failure to follow written instructions enclosed with the Product, accident, subjecting the Product conditions, or modification, alteration or repair by anyone other than A&D or persons authorized by A&D, or (iii) expendable or consumable components.

THIS LIMITED WARRANTY IS THE ONLY WARRANTY PROVIDED BY A&D; THERE ARE NO OTHER EXPRESS WARRANTIES. If A&D cannot reasonably repair or replace the Product, A&D will refund the amount You paid for the Product (not including taxes), less a reasonable charge for usage. To receive a refund you must have returned the Product and all associated materials to A&D.

The above remedy of repair, replacement or refund is your only and exclusive remedy. IN NO EVENT SHALL A&D BE LIABLE FOR ANY DAMAGES, INCLUDING, WITHOUT LIMITATION, LOST PROFITS, LOST INFORMATION OR REPLACEMENT COSTS, ARISING OUT OF YOUR USE OF OR INABILITY TO USE THE PRODUCT, INCLUDING, WITHOUT LIMITATION, ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, EVEN IF A&D HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. Some states do not allow the exclusion of incidental or consequential damages, so that the above exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that may vary from state to state.

No distributor, dealer or other party is authorized to make any warranty on behalf of A&D or to modify this warranty, or to assume for A&D any liability with respect to its products.

Technical Data

Type	UA-770CNBLE
Measurement method	Oscillometric measurement
Memory	Last 60 measurements per user
Measurement range	Pressure: 0-299 mmHg Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg Pulse: 40-180 beats / minute
Static accuracy	Pressure: ± 3 mmHg Pulse: $\pm 5\%$
Power supply	4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TB-233)
Number of measurements	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries) With pressure value 180 mmHg,
Classification	Internally powered ME equipment (Supplied by batteries) / Class II (Supplied by adapter) Continuous operation mode
Clinical test	ISO81060-2 : 2020
EMD	IEC 60601-1-2: 2014 + A1: 2020
Wireless Communication	Bluetooth: Ver.5.1LE BLP
Frequency band	2402 MHz to 2480 MHz
Maximum RF output power	< 10 dBm
Modulation	GFSK

Technical Data

Supported Data	Systolic Pressure, Diastolic Pressure, Pulse Rate
Communication distance	About 10 m (This distance is reduced by the conditions in the surrounding environment)
Paired device	10 devices
Operating conditions	+50°F to +104°F (10°C to 40°C) / 15% to 85% RH / 800 to 1060 hPa
Transport/Storage conditions	-4°F to +140°F (-20°C to 60°C) / 10% to 95% RH / 700 to 1060 hPa
Dimensions	140 mm (W) x 64 mm (H) x 118 mm (D) 5.5" (W) x 2.5" (H) x 4.6" (D)
Weight	10.4 oz. (295 g), excluding the batteries
Ingress protection	Device: IP20
Applied part	Cuff Type BF 
Useful life	Device: 5 years (when used six times a day) Cuff: 2 years (when used six times a day) Endurance 10,000 times AC Adapter: 5 years (when used six times a day)

Note: Specifications are subject to change without prior notice. IP classification is the degrees of protection provided by enclosures in accordance with IEC 60529. This device is protected against solid foreign objects of 12 mm diameter and greater such as a fingers. This device is not protected against water.

EMD Technical Data

Battery-operated or AC Adapter-operated Blood Pressure Monitor

Medical Electrical Equipment needs special precautions regarding EMD and needs to be installed and put into service according to the EMD information provided in the following.

Portable and mobile RF communication equipment (e.g. cell phones) can affect Medical Electrical Equipment.

The use of accessories and cables other than those specified may result in increased emissions or decreased immunity of the unit.

Table 1 - EMISSION Limits -

Phenomenon	Compliance
Conducted and radiated RF EMISSION CISPR 11	Group 1, Class B
Harmonic distortion IEC 61000-3-2	Class A
Voltage fluctuations and flicker IEC 61000-3-3	Compliance

Table 2 - IMMUNITY TEST LEVELS : Enclosure Port -

Phenomenon	IMMUNITY TEST LEVELS
Electrostatic discharge IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Radiated RF EM fields IEC 61000-4-3	10 V/m; 80 MHz - 2.7 GHz 80 % AM at 1 kHz
Proximity fields from RF wireless communications equipment IEC 61000-4-3	See table 4
Rated power frequency magnetic fields IEC 61000-4-8	30 A/m 50 Hz or 60 Hz
Proximity magnetic fields IEC 61000-4-39	See table 5

Table 3 - IMMUNITY TEST LEVELS : Input a.c. power Port -

Phenomenon	IMMUNITY TEST LEVELS
Electrical fast transients / bursts IEC 61000-4-4	± 2 kV 100 kHz repetition frequency

Technical Data

Table 3 -Continued

Phenomenon	IMMUNITY TEST LEVELS
Surges Line-to-line IEC 61000-4-5	± 0.5 kV, ± 1 kV
Conducted disturbances induced by RF fields IEC 61000-4-6	3 V 0.15 MHz - 80 MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz
Voltage dips IEC 61000-4-11	0 % U_T ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°
	0 % U_T ; 1 cycle and 70 % U_T ; 25/30 cycle Single phase: at 0°
Voltage interruption IEC 61000-4-11	0% U_T ; 250/300 cycle
NOTE U_T is the AC mains voltage prior to application of the test level.	

Table 4 - Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment -

Test frequency (MHz)	Band (MHz)	Service	Modulation	IMMUNITY TEST LEVEL (V/m)
385	380 - 390	TETRA 400	Pulse modulation 18 Hz	27
450	430 - 470	GMRS 460 FRS 460	FM ± 5 kHz deviation 1 kHz sine	28
710	704 - 787	LTE Band 13, 17	Pulse modulation 217 Hz	9
745				
780				
810	800 - 960	GSM 800/900 TETRA 800 iDEN 820 CDMA 850 LTE Band 5	Pulse modulation 18 Hz	28
870				
930				

Technical Data

Table 4 -Continued

Test frequency (MHz)	Band (MHz)	Service	Modulation	IMMUNITY TEST LEVEL (V/m)
1720	1700 - 1990	GSM 1800 CDMA 1900 GSM 1900 DECT LTE Band 1, 3, 4, 25 UMTS	Pulse modulation 217 Hz	28
1845				
1970				
2450	2400 - 2570	Bluetooth WLAN 802.11 b/g/n RFID 2450 LTE Band 7	Pulse modulation 217 Hz	28
5240	5100 - 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9
5500				
5785				

Table 5 - Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields-

Test frequency	Modulation	IMMUNITY TEST LEVEL(V/m)
30 kHz	CW	8
134.2 kHz	Pulse modulation 2.1 kHz	65
13.56 MHz	Pulse modulation 50 kHz	7.5

Technical Data

FCC CAUTION

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This transmitter must not be co-located or operated in conjunction with any other antenna or transmitter.

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment and meets the FCC radio frequency (RF) Exposure Guidelines in Supplement C to OET65. This equipment has very low levels of RF energy that it deemed to comply without maximum permissive exposure evaluation (MPE). But it is desirable that it should be installed and operated keeping the radiator at least 20 cm or more away from person's body (excluding extremities: hands, wrists, feet and ankles).

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

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A&D Company, Limited
A&D Engineering, Inc.