








# DIGITAL THERMOMETER



## KEY FEATURES:

-  Flexible tip
-  10 Second Measurement Time
-  Automatic Last Recall
-  One Button Operation
-  Audible Beep
-  Easy to read LCD
-  2 year warranty

## Accurate temperature readings at the touch of a button.

This easy to use thermometer provides 10 second measurements, a flexible tip for added comfort, and the ability to switch between measurements in Fahrenheit and Centigrade. Measurements can be taken orally, rectally or under the arm.

Features include temperature readout, last measurement in memory, an audible beep when a measurement is complete, and automatic shut-off. It includes a protective carrying case for safe storage, and a bilingual instruction guide with illustrations in English and Spanish.

## Technical Specifications

<b>Unit Dimensions (WxDxH) (in)</b>	4.9" x 0.9" x 0.5" (12.4 x 2.4 x 1.2 cm)
<b>Weight (lb)</b>	0.42 oz (12 grams) including battery
<b>Measurement Range</b>	89.6°F-109.2°F (32.0 to 42.9°C)
<b>Measurement Accuracy</b>	± 0.2°F (± 0.1°C)
<b>Power Supply</b>	1.5V button size battery
<b>Storage Range</b>	14°F to 131°F (-10°C to 55°C) when stored ≤ 95% RH
<b>Reading Duration</b>	10 seconds
<b>Number of Measurements (Battery)</b>	Approximately 200 hours of continuous operation

## Packaging and Shipping Specifications

	<b>Dimensions (LxWxH) (in)</b>	<b>Weight</b>	<b>Individual Item Cube (cubic feet)</b>	<b>UPC / GTIN</b>
<b>Individual Unit</b>	8.1 x 4.1 x 0.9	1.94 oz		0 93764 60312 9
<b>Case (20 units)</b>	9 x 9 x 4.5	48 oz	0.21	1 00 93764 60312 6
<b>Case (200 units)</b>	24 x 19 x 10	30 lb	2.70	2 00 93764 60312 3

## What is In The Box

- 1 Thermometer (pre-installed battery)
- 1 Instruction manual
- 1 Case

