








INSTANT READ DIGITAL EAR THERMOMETER



KEY FEATURES:

-  One Second Reading
-  10 Memory Recall
-  Audible Beep
-  Automatic Shutoff
-  One button, fully automatic operation
-  Easy to Read LCD
-  2 Year Warranty

Lightweight, Compact Ear Thermometer

A lightweight, compact ear thermometer that provides quick and accurate temperature measurements at the push of a button.

Quick and Easy Measurement

- Easy measurements in 1 second
- Measurements are taken through the ear – without the need for a protective lens cover
- Illuminated LCD display for easy reading
- 10 memory recall

Trusted Accuracy

- Scans the ear 512 times in 1 second – providing the utmost in accuracy
- An audible beep alerts when the measurement is complete
- Measurements in Fahrenheit or Celsius

User Comfort

- Ergonomically designed wide grip
- Unique shape allows for a gentle and safe ear temperature reading

Convenience

- Easy to understand fever indicator
- Suitable for adults and children of all ages
- Easy to clean
- Measurement can be taken on a sleeping child
- Auto off to save battery

Technical Specification

Measurement range	50°F – 122°F (10°C – 50°C)
Measurement accuracy	±0.4°F (0.2°C): From 96.8°F to 102.2°F (36°C – 39°C) ±0.5°F (0.3°C): Less than 96.8°F (36°C) or greater than 102.2°F (39°C)
Dimensions	2.2" x 1.75" x 6.2" (55.3 x 44.4 x 158.5 mm)
Weight	2.3 oz. (65.5 g)
Display	0.1°F or °C
Reading duration	1.0 Second
Memory function	10 memory recall
Power supply	One 3 Volt CR-2032 Lithium battery
Number of measurements (Battery)	Approximately 5,000 measurements
Operating conditions	60.8°F – 95°F (16°C – 35°C) with relative humidity up to 95%
Storage environment	-4°F – 122°F (-20°C – 50°C) with relative humidity up to 95%

Packaging and Shipping Specifications

	Dimensions (in) (WxDxH)	Weight (lb)	Individual Item Cube (cubic feet)	UPC
Selling Unit (Clamshell) UT-302	5.9 x 2.0 x 9.1	0.40	0.06	0-93764-60224-5
Unit Case (10 units) UT-302	6.2 x 22 x 9.7	5.0	0.77	

What is in the Box

- Digital ear thermometer
- CR-2032 lithium battery
- Trilingual instruction manual
- Trilingual quick start card

