

BASIC BLOOD PRESSURE MONITOR



CLINICALLY
VALIDATED

KEY FEATURES:

-  One Button Operation
-  Irregular Heartbeat Detection
-  Blood Pressure Rating Indicator
-  15 Reading Memory
-  5 Year Warranty
-  Cuff
9.0-14.6" (23-37 cm)

A classic blood pressure monitor with the benefits of speed and accuracy at a touch of a button.

Quick and Easy Measurement

- Easy one button operation for quick and accurate measurement
- Stores 15 measurements in memory with recall
- Displays average of all measurements

Trusted Accuracy

- Clinically validated for accuracy
- Indicates an irregular heartbeat (IHB) during a measurement
- Pressure Rating Indicator classifies blood pressure in accordance with medically approved standards

User Comfort

- Gentle Inflation enables a more comfortable measurement

All Our Blood Pressure Monitors Come With:

- 5 year warranty
- Cuffs are not made with natural rubber latex
- Battery life indicator
- A&D Connect companion app, download for free to manually log and electronically track measurements (available on iOS and Android)

Technical Specification

Measurement Method	Oscillometric measurement
Measurement Range	Pressure: 0-299 mmHg Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg Pulse: 40-180 beats /min
Measurement Accuracy	Pressure: ±3 mmHg Pulse: ±5%
Power Supply	4 x 1.5V batteries (R6P, LR6 or AA)
Number of Measurements	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries)
Monitor Dimensions	3.8 x 2.7 x 5.1"
Monitor Weight	8.1 oz without batteries
Display Size	1.7 x 1.9"
Clinical Test	ANSI / AAMI SP-10 1992
Operating Conditions	+50°F to +104°F (10°C to 40°C) / 15% to 85% RH / 800 to 1060 hPa
Storage / Transport Conditions	-4°F to +140°F (-20°C to 60°C) / 10% to 95% RH / 700 to 1060 hPa

Packaging and Shipping Specifications

	Dimensions (in) (HxLxD)	Weight (lb)	Case Cube (cubic feet)	UPC
Individual Unit Box	3.4 x 4.7 x 6.5	0.83	0.03	093764603174
Case (10 units)	18.5 x 10 x 7.3	11	0.78	

What is In The Box

- Monitor
- Instruction Manual
- Quick Start Card

