









PREMIUM UPPER ARM BLOOD PRESSURE MONITOR



KEY FEATURES:

-  4 User Compatible
-  Snapshot Averaging
-  Irregular Heartbeat Detection
-  Pressure Rating Indication
-  60 Memory Recall per user
-  Date and Time Stamps
-  Body Movement Sensor
-  5 year warranty

Top rated, multi-user upper arm blood pressure monitor that is clinically validated for accuracy.

Trusted Accuracy

- Clinically validated for accuracy
- Top-rated by several independent publications
- Irregular Heartbeat (IHB) alerts the user to irregular heartbeat patterns
- Pressure Rating Indicator classifies blood pressure in accordance with the medically approved standards

User Friendly

- Simple operation for up to 4 different users with a unique 60 reading memory for each user
- Snapshot Averaging automatically calculates the average for all readings to determine the average baseline blood pressure

- AC adapter available for operation without batteries (included in UA-767FAC model)
- Easy to apply wide-range cuff that comfortably fits most arms (8.6-16.5")
- Gentle inflation enables a more comfortable measurement

The A&D Promise

- Professional grade performance you can trust
- 5 year warranty on monitor, 2 year warranty on cuff
- Cuffs are not made with natural rubber latex
- Battery life indicator
- A&D Connect companion app, download for free to manually log and electronically track measurements (available on iOS and Android)

Technical Specifications

Measurement Method	Oscillometric
Memory	60 measurements for four users
Cuff Size	8.6-16.5" / 22-42cm
Measurement Range	Pressure: 0-299 mmHg Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg Pulse: 40-180 beats / minute
Measurement Accuracy	Pressure: ±3 mmHg Pulse: ±5%
Power Supply	4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TB-233)
Number of Measurements	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries)
Monitor Dimensions	5.5 x 2.3 x 4.1"
Monitor Weight	9.0 oz without batteries
Display Size	1.7 x 1.9"
Clinical Test	ISO81060-2: 2013
Operating Conditions	+50°F to +104°F (10°C to 40°C) / 15% to 85% RH / 800 to 1060 hPa
Storage/Transport Conditions	-4°F to +140°F (-20°C to 60°C) / 10% to 95% RH / 700 to 1060 hPa

Packaging and Shipping Specifications

	Dimensions (WxDxH) (in)	Weight (lb)	Individual Item Cube (cubic feet)	UPC / GTIN
UA-767F Individual Unit (Box)	5.1 x 3.9 x 7.7	1.57	0.05	0-93764-60336-5
UA-767F Case (10 units)	20 x 11 x 9	16.54	1.08	
UA-767FAC Individual Unit (Box)	5.1 x 3.9 x 7.7	1.67	0.05	0-93764-60337-2
UA-767FAC Case (10 units)	20 x 11 x 9	17.6	1.08	

What is In The Box

- 1 Upper arm blood pressure monitor
- 1 Blood pressure cuff (8.6-16.5")
- 1 Instruction manual
- 1 Quick start guide
- 1 Carrying case
- 1 AC adapter (UA-767FAC)
- 4 AA batteries (UA-767F)

