

THE NEXT GENERATION OF WRIST BLOOD PRESSURE MONITOR



ULTRACONNECT couples industry leading technology with the ultimate in ease of use

A simple path to better health. ULTRACONNECT Wireless Wrist Blood Pressure Monitors have a modern, portable design – making it easy to store results and track your progress. Using the A&D Connect App, you can wirelessly share readings with doctors and care providers. Like all of the trusted A&D blood pressure monitors, the ULTRACONNECT monitor is clinically validated for accuracy.

KEY FEATURES:



COMPACT DESIGN

All-in-One device makes it the smallest and lightest product on the market



MEMORY

Stores up to 100 readings on the device. Stores unlimited readings on the A&D Connect app



SMART CONNECT

The A&D Connect App automatically senses the presence of the monitor – providing the ultimate ease of use



CORRECT FIT & MOVEMENT SENSORS

Ensures that the cuff is positioned correctly giving accurate readings every time



OFFLINE MODE

Blood pressure measurements can be taken without the phone – allowing increased flexibility



IRREGULAR HEARTBEAT DETECTION

Alerts the user to an irregular heartbeat that may require further consultation



MULTI USER CAPABILITY

Allows up to 5 users and a guest. Seamless use by the entire family



WARRANTY

5 Year Warranty

Technical Specification

Measurement method	Oscillometric measurement
Measurement range	Pressure: 0 - 299mmHg; Systolic pressure: 60 - 279mmHg; Diastolic pressure: 40 - 200mmHg; Pulse: 40 - 180 beats /min
Measurement accuracy	Pressure: ± 3 mmHg; Pulse: $\pm 5\%$
Monitor dimensions	2.1" x 2.9" x 0.8" (53 x 74 x 20mm)
Weight	0.20 lb including cuff
Cuff size	5.3" x 8.5" (13.5 x 21.6cm)
Power supply	3.7V Li-ion 325mAh
Battery life	100 measurements
Clinical test	According to ISO81060-2
Wireless communication	ZY (MURATA Manufacturing Co.Ltd.) Bluetooth Ver. 4.1LE BLP
Operating condition	+10 to +40 °C / 15 to 85% RH / 800 to 1060 hPa
Transport/Storage condition	-20 to +60 °C / 10 to 95% RH / 700 to 1060 hPa
Ingress protection	Device: IP20

Packaging and Shipping Specifications

	Dimensions (in) (HxLxD)	Weight (lb)	Case Cube (cubic feet)	UPC
Individual Unit Box	3.75 x 2.85 x 4.6	0.50	0.028	093764603709
Case (20 units)	15 x 9 x 9.5	10.7	0.70	

Unlock the potential



A&D Connect is compatible with Bluetooth Smart mobile devices from Apple and many Android device manufacturers.

APPLE:

iOS version 6 or newer
iPhone 4s or newer, iPad 3rd & 4th generation or newer, iPad Mini

ANDROID:

OS version 4.4 or 5.1, 6.0
Requires Bluetooth Smart Ready™ device (Bluetooth 4.0).

For a list of supported devices, visit
www.wellnessconnected.jp/english/compatibility

What is in the box

- Monitor
- Quick Start Card
- Instruction Manual
- Charging Cable

