

DIGITAL WRIST BLOOD PRESSURE MONITOR



KEY FEATURES:

-  One Button Operation
-  Irregular Heartbeat Detection
-  Blood Pressure Rating Indicator
-  90 Reading Memory
-  5 Year Warranty
-  Wrist Cuff
5.3-8.5" (13.5-21.5 cm)

A classic wrist blood pressure monitor with the benefits of speed and accuracy at a touch of a button.

Quick and Easy Measurement

- Easy one button operation for quick and accurate measurement
- Stores 90 measurements in memory with recall
- Displays average of all measurements

Trusted Accuracy

- Indicates an irregular heartbeat (IHB) during a measurement
- Pressure Rating Indicator classifies blood pressure in accordance with medically approved standards

User Comfort

- Easy to use wrist cuff - 5.3-8.5" (13.5-21.5 cm)
- Gentle Inflation enables a more comfortable measurement

Convenience

- Compact size for storage and travel

All Our Blood Pressure Monitors Come With:

- 5 year warranty
- Cuffs are not made with natural rubber latex
- Battery life indicator
- A&D Connect companion app, download for free to manually log and electronically track measurements (available on iOS and Android)

Technical Specification

Measurement Method	Oscillometric measurement
Measurement Range	Pressure: 20-280 mmHg Pulse: 40-180 beats /min
Measurement Accuracy	Pressure: ±3 mmHg Pulse: ±5%
Power Supply	Two (2) 1.5V alkaline batteries (LR03 or AAA)
Battery Life	Approximately 400 measurements
Operating Conditions	+10 to +40 °C / 15 to 85% RH / 800 to 1060 hPa
Transport/Storage Conditions	-20 to +60 °C / 10 to 95% RH / 700 to 1060 hPa
Monitor Dimensions	Approx. 2.4 [W] x 2.6 [H] x 1.2 [D] in
Warranty	5 years

Packaging and Shipping Specifications

	Dimensions (in) (HxLxD)	Weight (lb)	Case Cube (cubic feet)	UPC
Individual Unit Box	4.6 x 3.5 x 6.1	0.71	0.06	093764602139
Case (10 units)	18 x 10 x 7	8.40	0.73	

What is In The Box

- Monitor
- Instruction Manual
- Quick Start Card
- Carrying Case
- 2 AAA Batteries

