









# PREMIUM WRIST BLOOD PRESSURE MONITOR



## KEY FEATURES:

-  Multi-user
-  Correct Positioning Guide
-  Irregular Heartbeat Detection
-  Pressure Rating Indication
-  60 Memory Recall per user
-  Date and Time Stamps
-  Averaging for AM, PM and All Readings
-  5 year warranty

## Advanced Multi-user Blood Pressure Monitor

The UB-543 wrist blood pressure monitor provides a superior set of features all in a simple, easy to use design. The large display and one button operation allow for quick and accurate measurements.

### Advanced Features

**Clinical Validation:** This blood pressure monitor is clinically validated using the European Society of Hypertension International Protocol, a gold standard in clinical validation that ensures consistently accurate blood pressure readings.

**Correct Positioning Guide:** The guide indicates when the wrist is in the correct position, ensuring accurate readings taken at heart level.

**Irregular Heartbeat Detection:** Identifies the presence of an irregular heartbeat caused by arrhythmia or other conditions without impacting the accuracy of the blood pressure reading.

**Automatic Averaging:** This feature automatically calculates the average for all readings, AM readings, and PM readings to determine the average baseline blood pressure. This will allow users to more easily identify spikes and abnormal readings.

## Technical Specifications

<b>Dimensions</b>	2.2" [W] x 0.7" [D] x 3.5" [H]
<b>Weight</b>	3.6 oz, without batteries
<b>Display</b>	2.4" x 1.6" (60.6 x 41.0 mm)
<b>Memory</b>	60 readings per user
<b>Pressurization</b>	Automatic, using micropump
<b>Measurement Range</b>	<b>Pressure:</b> 0-299 mmHg <b>Pulse:</b> 40-180 beats /min
<b>Accuracy</b>	<b>Pressure:</b> ±3 mmHg <b>Pulse:</b> ±5%
<b>Deflation</b>	Constant air release valve
<b>Type</b>	Oscillometric
<b>Operating Conditions</b>	50°F to 104°F (+10°C to +40°C) / 15% RH to 85% RH
<b>Storage Conditions</b>	-4°F to 140°F (-20°C to +60°C) / 10% RH to 95% RH
<b>Cuff Size</b>	Wrist 5.3"-8.5" (13.5 - 21.5 cm)
<b>Power Source</b>	2 x 1.5v alkaline batteries (LR03 or AAA)

## Packaging and Shipping Specifications

	<b>Dimensions (WxDxH) (in)</b>	<b>Weight (lb)</b>	<b>Individual Item Cube (cubic feet)</b>	<b>UPC</b>
<b>Individual Unit (Box)</b>	4.5" x 3.4" x 6"	0.64 lb	0.05 cubic feet	093764603358
<b>Case (10 units)</b>	18" x 10" x 7"	7.7 lb	0.73 cubic feet	

## What is In The Box

- Wrist blood pressure monitor
- Instruction manual
- Quick start guide
- Carrying case
- 2 AAA batteries

