

A&D Medical ... *Clearly a Better Value*

PRECISION SCALE UC-324



The UC-324 Precision Scale features an easy-to-read LCD display in a low profile, glass top form factor. The scale is highly accurate, and is designed with ease of operation in mind. There are no buttons to press - just tap the scale with your foot to turn it on, step on and view your reading on the large display. In addition, an auditory “beep” will indicate every action. The non-slip rubber feet ensure increased safety while stepping on to the scale to take a measurement. The UC-324 helps you manage and reach your personal weight goals to keep you moving towards better health. The scale will fit into your daily weight tracking record.

- One-Touch Activation
- Fast and easy to get your results and keep on moving
- Hi Capacity - 330lb (150kg)
- Precision Accuracy - 0.2lb (0.1kg) resolution
- Auditory Alert - Indicates the readiness to step on and the measurement is done
- Weigh units – Easy to change lb – kg mode

For questions about A&D products, call us toll free at 1-888-726-9966 and a specially trained representative will assist you, or visit us at www.andmedical.com.



- One Touch Activation
- Capacity of 330lb (150kg)
- Precision Accuracy
- Auditory alert


UC-324 Wireless Precision Scale

Contents Include:

- Wireless Precision Scale
- Instruction Manual
- 4 AA Batteries

Technical Specifications

Measuring Method	Electric Resistance
Digital Display	Liquid Crystal Display
Maximum Capacity (MAX)	330 lbs / 150 kg
Minimum Display (d)	0.2 lbs / 0.1 kg
Power Source	Four Type AA Batteries (1.5 V)
Dimensions	13.8 (D) x 13.8 (W) x 1.5 (H) in 350 (D) x 350 (W) x 39 (H) mm
Weight	4.8 lbs (2.2 kg), including the batteries



Packaging and Shipping Specifications

Selling Unit (Box)	W x D x H (inches)	Wt (lb.)	Individual Item Cube (cubic feet)
UC-324	16.4 x 3.3 x 16.3	6.64	0.51
Unit Case (5 units)	W x D x H (inches)	Wt (lb.)	Case Cube (cubic feet)
UC-324	17.8 x 17.6 x 17.5	36.0	3.17

Ask your A&D Medical Sales Representative about the entire Wellness Connected family of products.

Ask your A&D Medical Sales Representative about other A&D Medical & Wellness Connected products.

- Blood Pressure Monitors
- Digital Thermometers
- Personal Scales
- Stethoscopes
- Activity Monitors

For questions about A&D Medical products, call us toll-free at 1-888-726-9966 or visit one of our websites at:
www.andmedical.com
www.WellnessConnected.com
www.LifeSourceOnline.com

A&D
A&D Medical
1756 Automation Parkway
San Jose, CA 95131
1-888-726-9966
www.andmedical.com

LIFE SOURCE
A&D Medical

