












*Using the correct size cuff is important for an accurate reading*

To determine what cuff size is appropriate for you, measure the circumference of your upper arm at the midpoint between the shoulder and elbow



	Small UA-279 6.3-9.4"	Medium UA-280 9.4-14.2"	Large UA-281 14.2-17.7"	Small UA-289 6.3-9.4"	Medium UA-290 9.0-14.6"	Large UA-291 12.2-17.7"	Pre-Formed UA-381 9-17"	Smooth Fit UA-392 9.0-14.6"	X-Large UA-282 16.5-23.6"	AccuFit Plus (Wide Range) UA-420 8.6-16.5"
<b>UA-611</b> 					✓	✓				✓
<b>UA-651</b> 					✓	✓				✓
<b>UA-651BLE</b> 					✓	✓				✓
<b>UA-767F</b> 					✓	✓				✓
<b>UA-767FAC</b> 					✓	✓				✓
<b>UA-705V/VL</b> 		✓	✓							
<b>UA-767PSAC</b> 	✓									
<b>UA-787EJ</b> 	✓	✓	✓				✓			
<b>UA-789AC</b> 		✓	✓						✓	
<b>UA-1030T</b> 				✓	✓	✓		✓		