


Easy on-PC

PC spirometer



The Easy on-PC is a fully integrated spirometer that leverages our patented TrueFlow technology to deliver unprecedented accuracy and reliable results.

Simple. Certain. Proven.



Features

- Real-time animated incentives
- Proven EMR integration with top EMR/EHR systems
- Long term trending of patient results
- Automated quality feedback and interpretation
- Intuitive software interface
- Integrated provocation test with customizable protocols
- Instant quality feedback and interpretation
- Wide range of selectable predictable values including the new GLI predicted set
- Calibration free technology

Specifications

Test Types: FVC, FVL, SVC, MVV, Challenge and Pre/Post testing

Languages: English, French, German, Spanish, Brazilian-Portuguese, Dutch, Russian, Vietnamese, Turkish

EMR connectivity: EasyOne Connect

Standards: Meets ATS/ERS, NIOSH/OSHA and SSA Disability

Hard disk capacity: Installation/system 1 GB Data up to 4 GB

RAM: 2 GB

Operating system: Windows 7, Windows 8 and 8.1(32 and 64 bit), Windows 10 (32 and 64 Bit)

Part Numbers

2700-3: Easy on-PC Spirometer

2700-3K: Easy on-PC Spirometry Kit

2050-1: Spirette box of 50

2050-5: Spirette box of 200

TrueFlow

The original ultrasonic flow measurement is highly accurate in all flow ranges, independent of gas composition, pressure, temperature and humidity.

ndd TrueFlow is a hygienic and resistance-free solution that does not require calibration during its lifetime.

EasyOneConnect

EasyOne Connect software offers one database and one platform for all EasyOne point-of-care-solutions. Data management has never been easier.

Connect with Us



ndd Medical Technologies



nddmedical



nddMedical



nddmed