

Fast. Easy. Accurate.

ANEROID BLOOD PRESSURE KITS UA-101 & UA-201



Precision craftsmanship and easy to use blood pressure monitors for consumers (UA-101, left) and professionals (UA-201, right)


The UA-101 and UA-201 Aneroid Blood Pressure Kits provide consumers and health care professionals with reliable monitoring products at home or in a clinical setting. These aneroid blood pressure kits can assist in the management of hypertension, improve patient compliance to treatment of high blood pressure, and be used as tools in a preventative health management program.

- Nylon cuff for superior durability and appearance
- Latex bulb with standard air release valve
- Includes stethoscope (Model UA-101 only)
- Two year limited warranty for manometer accuracy

For questions about A&D Medical or LifeSource products, call us toll-free at 1-800-726-9966 and a specially trained representative will assist you, or visit us at www.andmedical.com.



- Professional Accuracy
- Includes Zippered Carrying Case
- Calibration Screwdriver Included



UA-101



UA-201

A&D Medical Aneroid Blood Pressure Kits

A&D Medical offers two aneroid blood pressure kits. The Home Blood Pressure Kit with Attached Stethoscope (Model UA-101) is designed to give consumers a way to track their readings at home. The Professional Sphygmomanometer (Model UA-201) is designed for health care professionals to measure a patient's blood pressure.

PART #	UPC NUMBER	CUFF SIZE
UA-101	0-93764-60308-2	ADULT 10.0" - 14.0" (25 - 36cm)
UA-201	0-93764-60309-9	ADULT 10.0" - 16.0" (25 - 41cm)

Home Blood Pressure Kit with Attached Stethoscope Model UA-101



- No-pin-stop manometer with easy to calibrate feature
- Comes with attached stethoscope and D-ring cuff
- Latex 2-tube bladder
- Includes bilingual instruction manual
- Zippered storage case and calibration screwdriver included

Professional Sphygmomanometer Model UA-201



- Arrow indicator on cuff ensures proper cuff placement
- Comes with extended adult size cuff
- Additional cuff sizes are available (see below)
- Zippered storage case and calibration screwdriver included

Accessories

PART #	DESCRIPTION	PART #	DESCRIPTION	PART #	DESCRIPTION
11-001	Inflator bulb	11-006	Adult cuff w/bladder	11-012	Large cuff w/bladder
11-002	Exhaust valve	11-007	Adult bladder	11-013	Large bladder
11-003	End valve	11-008	Infant cuff w/bladder	11-014	Thigh cuff w/bladder
11-004	Gauge	11-009	Infant bladder	11-015	Thigh bladder
11-005	Zippered case	11-010	Child cuff		
UA-700CH	Cuff hose	11-011	Child bladder		

Packaging and Shipping Specifications

	UA-101	UA-201
Selling Unit (box)		
L x W x H (inches)	7.3 x 4.1 x 3.0	7.3 x 4.1 x 3.0
Weight (lbs.)	1.08	0.86
Item Cube (cubic feet)	0.05	0.05
Unit Case (20 units)		
L x W x H (inches)	15.4 x 8.9 x 15.4	15.4 x 8.9 x 15.4
Weight (lbs.)	24.3	19.8
Case Cube (cubic feet)	1.21	1.21



1756 Automation Parkway
San Jose, CA 95131
1-888-726-9966
www.andmedical.com

LifeSource is a registered trademark of A&D Medical. ©2012 A&D Medical. All rights reserved.



Ask your A&D Medical Sales Representative about other A&D Medical products.

- Blood Pressure Monitors
- Digital Thermometers
- Personal Scales
- Stethoscopes
- Activity Monitors

For questions about A&D Medical products, call us toll-free at 1-888-726-9966 or visit one of our websites at: www.andmedical.com

LIT-UA-101-201-DS 08/2012