

ULTRACONNECT BLOOD PRESSURE MONITORS



Upper Arm
UA-1200BLE

Wrist
UB-1100BLE



KEY FEATURES:

-  Wireless
-  One button operation
-  Irregular Heartbeat Detection
-  Smart Connect
-  Unlimited memory on app
-  Multi-User Capability
-  Portable All-In-One Design
-  Fit & Movement Sensors

The next generation of blood pressure monitoring, providing the ultimate user experience and accuracy.

Trusted Accuracy

- Clinically validated for accuracy
- Indicates an irregular heartbeat (IHB) during a measurement
- Body movement sensor guides user to remain still during measurement
- Cuff fit error prevents inaccurate measurements by alerting user when the cuff is too loose

User Friendly

- Easy to use for up to 5 users and a guest
- Stores 100 measurements in monitor memory, plus unlimited readings on the Heart Track app
- Free companion A&D Heart Track app automatically stores readings, tracks progress with simple graphs and allows you to share results with others

- Send your measurements wirelessly via Bluetooth to your supported phone/tablet
- Smart Connect automatically detects and syncs measurements between the monitor and the Heart Track app
- Blood pressure measurements can be taken without your phone or tablet – allowing increased flexibility
- Compact, all-in-one device is small, light, and perfect for travel in both upper arm and wrist configurations

The A&D Promise

- Professional grade performance you can trust
- 5 year warranty
- Cuffs are not made with natural rubber latex



Technical Specifications

Measurement Method	Oscillometric
Memory	Unlimited on App
Cuff Size	8.6-16.5" (22-42 cm) (UA-1200BLE)
	5.3-8.5" (13.4-22 cm) (UB-1100BLE)
Measurement Range	Pressure: 0-299 mmHg Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg Pulse: 40-180 beats / minute
Measurement Accuracy	Pressure: ±3 mmHg Pulse: ±5%
Power Supply	3.7V Li-ion 325mAh
Battery Life	100 measurements
Monitor Dimensions	5 x 1 x 1.7" (13 x 3 x 4.3 cm) (UA-1200BLE)
	2.1 x 0.8 x 2.9" (5.3 x 2 x 7.4 cm) (UB-1100BLE)
Monitor Weight	0.45 lb (204 g) (including cuff (UA-1200BLE)
	0.20 lb (91 g) including cuff (UB-1100BLE)
Clinical Test	According to ISO81060-2
Operating Conditions	+50°F to +104°F (10°C to 40°C) / 15% to 85% RH / 800 to 1060 hPa
Storage/Transport Conditions	-4°F to +140°F (-20°C to 60°C) / 10% to 95% RH / 700 to 1060 hPa
Ingress Protection	Device: IP20

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
UA-1200BLE Individual Unit Box	5.3 x 2.8 x 8.1 in (13.4 x 7.1 x 21 cm)	1.52 lb (0.69 kg)	0.09 ft ³ (0.0025 m ³)	0-93764-60369-3 / 1-00-93764-60369-0
UA-1200BLE Case	10.7 x 8.2 x 13.9 in (27 x 21 x 35 cm)	11.69 lb (5.3 kg)	0.70 ft ³ (0.02 m ³)	
UB-1100BLE Individual Unit Box	3.6 x 2.9 x 4.6 in (9.5 x 7.2 x 12 cm)	0.50 lb (0.23 kg)	0.028 ft ³ (0.0008 m ³)	0-93764-60370-9 / 1-00-93764-60370-6
UB-1100BLE Case	15 x 9 x 9.5 in (38 x 23 x 24 cm)	10.7 lb (4.9 kg)	0.70 ft ³ (0.02 m ³)	

What is In The Box

- Blood Pressure Monitor
- Quick Start Guide
- Instruction Manual
- USB Charging Cable