




BLOOD PRESSURE MONITOR WITH EXTRA LARGE CUFF




KEY FEATURES:


 Extra Large Cuff
16.5-23.6" (42-60 cm)

 Averages all
measurements

 Reading Reminder

 Irregular Heartbeat
Detection

 60 Reading Memory

 5 Year Warranty



CLINICALLY
VALIDATED

The only home blood pressure monitor designed to provide accurate readings for people with larger arms.

Trusted Accuracy

- Clinically validated for accuracy
- Indicates an Irregular Heartbeat during a measurement
- Pressure Rating Indicator classifies blood pressure in accordance with medically approved standards

User Friendly

- Stores 60 measurements in memory
- Displays average of all measurements
- Clock display with time and date stamp for logging and trending measurements
- 3 programmable reminders to ensure user takes measurements at the same time each day

- Extra large cuff to accommodate people with larger arms (16.5-23.6")
- Fast and accurate readings with the high speed motor

The A&D Promise

- 5 year warranty on the monitor/ 2 years on the cuff
- Cuffs are not made with natural rubber latex
- Battery life indicator
- Download the A&D Heart Track App for Free to manually log and electronically track measurements (available on iOS and Android)

Technical Specification

Measurement Method	Oscillometric
Memory	60 measurements
Cuff Size	16.5-23.6" (42-60 cm)
Measurement Range	Pressure: 0-280 mmHg Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg Pulse: 40-200 beats /min
Measurement Accuracy	Pressure: ±3 mmHg Pulse: ±5 %
Power Supply	AC adapter or 4 x 1.5V batteries (R6P, LR6 or AA) (not included)
Battery Life	Approximately 120 measurements
Monitor Dimensions	6.4 x 4.4 x 2.4" (16.3 x 11.2 x 6.2 cm)
Monitor Weight	12.3 oz (350 g)
Display Size	2.2 x 2.8" (5.5 x 7 cm)
Clinical Test	According to ANSI/AAMI SP-10 1992
Operating Conditions	50°F to 104°F (+10°C to +40°C) / 15% to 85% RH / 800 to 1060 hPa
Storage/Transport Conditions	-4°F to 140°F (-20°C to +60°C) / 10% to 95% RH / 700 to 1060 hPa

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
UA-789AC Individual Unit (Box)	8.25 x 4.25 x 5.5" (21 x 11 x 14 cm)	2.2 lb (1.0 kg)	0.18 ft³ (0.0051 m³)	0-093764-60062-3
UA-789AC Case (10 units)	28.4 x 15.8 x 11.1" (72 x 40 x 28 cm)	36.8 lb (16.7 kg)	2.86 ft³ (0.081 m³)	1-00-93764-60062-0

What is In The Box

- 1 Blood pressure monitor
- 1 Blood pressure cuff (16.5-23.6")
- 1 Instruction manual
- 1 Quick start guide
- 1 Carrying case
- 1 AC adapter