

PREMIUM WRIST BLOOD PRESSURE MONITOR



KEY FEATURES:



2 User Compatible



Correct Positioning Guide



Irregular Heartbeat Detection



Pressure Rating Indication



60 Memory recall per user



Date and Time Stamps



Averaging for AM, PM and All Readings



5 year warranty

Advanced, multi-user blood pressure monitor that provides trusted accuracy and reliable measurements every time.

Trusted Accuracy

- Clinically validated for accuracy
- Indicates an irregular heartbeat (IHB) during a measurement
- Correct Positioning Guide helps ensure accurate readings
- Pressure Rating Indicator classifies blood pressure in accordance with medically approved standards

User Friendly

- 2 different users with a unique 30 reading memory for each user
- Easy, one button operation for quick and accurate measurement
- Stores 60 measurements in memory with recall, for two users, including average of all measurements

- Clock display with time and date stamp for logging and trending measurements
- Snapshot averaging calculates the average for all readings, AM readings, and PM readings to determine the average baseline blood pressure.
- Easy to use wrist cuff (fits 5.3-8.5" wrists)

A&D Promise

- Professional grade performance you can trust
- 5 year warranty on monitor
- Cuffs are not made with natural rubber latex
- Battery life indicator
- Download the A&D Heart Track App for Free to manually log and electronically track measurements (available on iOS and Android)

Technical Specifications

Measurement Method	Oscillometric
Memory	60 readings per user
Wrist Circumference	Wrist 5.3-8.5" (13.5-21.5 cm)
Measurement Range	Pressure: 0-299 mmHg Systolic pressure: 30-280 mmHg Diastolic pressure: 20-200 mmHg Pulse: 40-180 beats /min
Measurement Accuracy	Pressure: ± 3 mmHg Pulse: $\pm 5\%$
Power Supply	2 x 1.5v alkaline batteries (LR03 or AAA)
Battery Life	Approximately 180 measurements
Monitor Dimensions	2.2 x 0.7 x 3.5" (5.6 x 1.8 x 8.8 cm)
Unit Weight	3.6 oz (102 g) without batteries
Display Size	1.6 x 2.4" (4.1 x 6.1 cm)
Clinical Test	According to ISO81060-2:2013
Operating Conditions	50°F to 104°F (+10°C to +40°C) / 15% to 85% RH / 800 to 1060 hPa
Storage/Transport Conditions	-4°F to 140°F (-20°C to +60°C) / 10% to 95% RH / 700 to 1060 hPa

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
Individual Unit (Box)	4.5 x 3.4 x 6" 11.4 x 8.6 x 15.2 cm	0.64 lb (0.29 kg)	0.05 ft ³ (0.0014 m ³)	0-93764-60335-8
Case (10 units)	18 x 10 x 7" 46 x 25 x 18 cm	7.7 lb (3.49 kg)	0.73 ft ³ (0.0207 m ³)	1-00-93764-60335-5

What is In The Box

- 1 Wrist blood pressure monitor
- 1 Instruction manual
- 1 Quick start guide
- 1 Carrying case
- 2 AAA batteries