

# INSTANT READ DIGITAL EAR THERMOMETER



## KEY FEATURES:



Easy Operation



10 Memory Recall



Big Illuminated Display



Measurement in Celcius or Fahrenheit



No Lens Cover Required



Fever Indicator



2 Year Warranty

## Lightweight, Compact Ear Thermometer

A lightweight, compact ear thermometer that provides quick and accurate temperature measurements at the push of a button.

### Quick and Easy Measurement

- Easy measurements in 1 second
- Measurements are taken through the ear – without the need for a protective lens cover
- Illuminated LCD display for easy reading
- 10 memory recall

### Trusted Accuracy

- Scans the ear 512 times in 1 second – providing the utmost in accuracy
- An audible beep alerts when the measurement is complete
- Measurements in Fahrenheit or Celsius

### User Comfort

- Ergonomically designed wide grip
- Unique shape allows for a gentle and safe ear temperature reading

### Convenience

- Easy to understand fever indicator
- Suitable for adults and children of all ages
- Easy to clean
- Measurement can be taken on a sleeping child
- Auto off to save battery

## Technical Specification

<b>Measurement range</b>	50°F – 122°F (10°C – 50°C)
<b>Measurement accuracy</b>	±0.4°F (0.2°C): From 96.8°F to 102.2°F (36°C – 39°C) ±0.5°F (0.3°C): Less than 96.8°F (36°C) or greater than 102.2°F (39°C)
<b>Dimensions</b>	2.2" x 1.75" x 6.2" (55.3 x 44.4 x 158.5 mm)
<b>Weight</b>	2.3 oz. (65.5 g)
<b>Display</b>	0.1°F or °C
<b>Reading duration</b>	1.0 Second
<b>Memory function</b>	10 memory recall
<b>Power supply</b>	One 3 Volt CR-2032 Lithium battery
<b>Number of measurements (Battery)</b>	Approximately 5,000 measurements
<b>Operating conditions</b>	60.8°F – 95°F (16°C – 35°C) with relative humidity up to 95%
<b>Storage environment</b>	-4°F – 122°F (-20°C – 50°C) with relative humidity up to 95%

## Packaging and Shipping Specifications

	<b>Dimensions (in) (WxDxH)</b>	<b>Weight (lb)</b>	<b>Individual Item Cube (cubic feet)</b>	<b>UPC</b>
<b>Selling Unit (Clamshell)</b> UT-302	5.9 x 2.0 x 9.1	0.40	0.06	0-93764-60224-5
<b>Unit Case (10 units)</b> UT-302	6.2 x 22 x 9.7	5.0	0.77	

## What is in the Box

- Digital ear thermometer
- CR-2032 lithium battery
- Bilingual instruction manual
- Bilingual quick start card

