

PROFESSIONAL OFFICE BLOOD PRESSURE MONITOR

KEY FEATURES:



Large LCD Display



Dual Measurement
(Auscultatory/Oscillometric)



Wide Range of Cuffs



Adjustable Inflation



Portable



3 Year Warranty



CLINICALLY
VALIDATED

The A&D Medical Professional Office Blood Pressure Monitor is clinical grade and simple to use, designed to deliver clinicians the highest level of precision and accuracy required to measure blood pressure. This economical blood pressure monitor is perfect for pharmacies, nursing homes, clinics, hospitals, or other medical office environments.

Dual Measurement Mode

Auscultatory and Oscillometric modes offer healthcare professionals the flexibility to switch between manual and automatic.

Portable Design

Carrying grip designed with convenient storage for cuffs. Optional rolling stand available for maximum portability.

Clinical-Grade Materials

Monitor can be cleaned using hospital grade sanitizers including ethanol, chlorhexidine gluconate, and benzalkonium chloride.

Comfortable Inflation

Customize the right amount of inflation with manual controls; 220 mmHg, 250 mmHg, 280 mmHg or use in AUTO mode.

Technical Specifications

Measurement Method	Oscillometric and Auscultatory
Measurement Range	Pressure: 0 to 299 mmHg Systolic Pressure: 60 to 279 mmHg Diastolic Pressure: 40 to 200 mmHg Pulse: 40 to 200 beats/minute
Measurement Accuracy	Pressure: ±3 mmHg Pulse: ±5%
Temperature Accuracy	±2.5°C (36.5°F) (+5°C to +40°C)
Temperature Unit	°C or °F
Battery Life	Approx. 300 Measurements
Clinical Test	According to ISO81060-2 2013 IEC 60601-1-2:2007
Memory	99 Measurements
Power Supply	Built-in rechargeable 3.6V Nickel-Metal Hydride Battery battery or AC adapter (TB-268)
Dimensions (WxHxD)	Approx. 4.72 x 7.87 x 5.51" (120 x 200 x 140 mm)
Unit Weight	Approx. 1.43 lbs (650 g) excluding batteries

Optional Accessories

Name	Order Code	Description (Arm circumference)
Small Cuff	CUF-KS-SA	6.3-9.4" (16-24 cm)
Medium Cuff	CUF-KS A	8.7-12.6" (22-32 cm)
Large Cuff	CUF-KS LA	12.2-17.7" (31-45 cm)
Extra Large Cuff	CUF-KS LL	16.1-19.7" (41-50 cm)
AC Adapter	TB-268	
Rolling Stand with Basket	UM-211-ST	
	UM-21x-ST-OR27	
Rechargeable Battery Pack	UM-211-20	
UM-211 with all 4 Cuffs Sizes	UM-211Kit	

CUF-KS-LL



CUF-KS-LA



CUF-KS-A



CUF-KS-SA



UM-211



UM-211-ST

Packaging and Shipping Specifications

Individual Box Dimensions	8.6 x 6.6 x 8.6" (22 x 17 x 22 cm)
Individual Box Weight	3.2 lbs (1.5 kg)
Master Case Dimensions	28 x 10 x 10" (71 x 25 x 25)
Master Case Weight	14 lbs (6.4 kg)
Units Per Master Case	4
UPC	0-9376460393-8

What is In The Box

- 1 Blood Pressure Monitor
- 1 Blood Pressure Cuff
- 1 Nickel-Metal Hydride Battery
- 1 AC Adapter