

WIRELESS UPPER ARM BLOOD PRESSURE MONITOR



KEY FEATURES:



Wireless



One button operation



Irregular Heartbeat
Detection



Unlimited memory on app



5 year warranty



CLINICALLY
VALIDATED

A wireless blood pressure monitor that automatically records and tracks measurements via mobile app to enable trending and sharing.

Trusted Accuracy

- Clinically validated for accuracy
- Irregular Heartbeat (IHB) alerts the user to irregular heartbeat patterns

User Friendly

- Easy, one button operation for quick and accurate measurement
- Send your measurements wirelessly via Bluetooth to your supported phone/tablet
- Free companion A&D Heart Track app automatically stores readings, tracks progress with simple graphs and allows you to share results with others
- 30 reading memory in off-line mode, unlimited on app
- Easy to apply wide range cuff that comfortably fits most arms (8.6-16.5" / 22-42 cm)
- Gentle inflation enables a more comfortable measurement

The A&D Promise

- Professional grade performance you can trust
- Simple integration to telehealth systems with full service technical support
- 5 year warranty on monitor, 2 year warranty on cuff
- Cuffs are not made with natural rubber latex
- Battery life indicator



Heart
Track
POWERED BY MMHG



Technical Specifications

| | |
|-------------------------------------|---|
| Measurement Method | Oscillometric |
| Memory | Unlimited on App |
| Cuff Size | 8.6-16.5" (22-42 cm) |
| Measurement Range | Pressure: 0-299 mmHg Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg Pulse: 40-180 beats / minute |
| Measurement Accuracy | Pressure: ±3 mmHg Pulse: ±5% |
| Power Supply | 4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TB-233) |
| Battery Life | Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries) |
| Monitor Dimensions | 3.8 x 5.1 x 2.7" (9.6 x 13 x 6.8 cm) |
| Monitor Weight | 8.8 oz (249 g) without batteries |
| Display Size | 1.4 x 2.1" (3.6 x 5.3 cm) |
| Clinical Test | ISO81060-2: 2013 |
| Operating Conditions | 50°F to 104°F (+10°C to +40°C) / 15% to 85% RH / 800 to 1060 hPa |
| Storage/Transport Conditions | -4°F to 140°F (-20°C to +60°C) / 10% to 95% RH / 700 to 1060 hPa |

Packaging and Shipping Specifications

| | Dimensions (WxLxH) | Weight | Individual Item Cube | UPC / GTIN |
|--|---------------------------------------|-------------------|-----------------------------|--------------------|
| UA-651BLE Individual Unit (Box) | 5.1 x 4 x 7.7" (13 x 10 x 19 cm) | 1.52 lb (0.7 kg) | 0.09 ft³ (0.0024 m³) | 0-93764-60332-7 |
| UA-651BLE Case (10 units) | 21.7 x 11 x 8.5" (55 x 28 x 22 cm) | 15.87 lb (7.2 kg) | 1.18 ft³ (0.033 m³) | 1-00-93764-60332-4 |

What is In The Box

- 1 Upper arm blood pressure monitor
- 1 Blood pressure cuff (8.6-16.5")
- 1 Instruction manual
- 1 Quick start guide
- 1 Carrying case
- 4 AA batteries