

TABLE-TOP BLOOD PRESSURE KIOSK

KEY FEATURES:



One Button Operation



Small Footprint



Result Printout



Irregular Heartbeat Indicator



Antibacterial Cuff



2 Year Warranty



*Stand and stool sold separately



CLINICALLY
VALIDATED

The A&D Medical Professional blood pressure kiosk has a small footprint, perfect for small pharmacies, mobile clinics or workplace screening. This durable monitor is designed to be used multiple times a day with simple directions to take proper measurements.

Accurate Readings

Clinically validated with rigorous protocols performed by independent researchers and published in peer reviewed journals.

Innovative Technology

The unique cuff inflation system stabilizes the arm and customizes the cuff circumference to the patient. This provides fast, accurate and reliable results.

Easy and Fast Operation

Simply press a button to operate the machine. An integrated instruction panel allows patients to easily take their blood pressure without the aid of a healthcare professional.

Integration

Seamlessly collect data via the optional Bluetooth or RS232C communication modules.

Flexible Applications

This monitor can be used as a health screening station, alone on a tabletop, or as part of system. This makes it the ideal solution for use in clinics and retail pharmacies.

Customized Result Printout

A printer is included in all models. The printout can be customized with a store or hospital name and the user can keep it as a log or provide it to the nurse or doctor. Measurements can be uploaded to the Heart Track Blood Pressure Management System via a QR code on the result printout.

Sanitary Solution

An anti-bacterial cuff cloth provides a sanitary solution. The cuff cloth is easily replaced and disposable covers are available.

Technical Specifications

Measurement Method	Oscillometric Measurement
Cuff Size	Fits Arms 7.1-13.8" (18-35 cm)
Measurement Range	Systolic pressure: 40-270 mmHg Diastolic pressure: 20-200 mmHg Pulse: 30-240 beats/min
Measurement Accuracy	Pressure: ±3 mmHg or 2% Pulse: ±5%
Display	3 Digit Display LED & LED Lamp
Time/Date	12 Hour Format, Year, Month, and Day
Power Supply	100-240V, 50-60 HZ
AC Power Consumption	50-80 VA
Clinical Test	EN1060-4:2004
Wired Communication	Optional RS-232C (TM-2657P-RS)
Wireless Connection	Optional <i>Bluetooth</i> ® (TM-2657P-BT)
Safety Standard	IEC60601-1:2012
EMD Compliance	Complies with Electromagnetic Disturbances Standard IEC60601-1-2:2014
Operating Environment Temperature	50°F to 104°F (10 to 40°C), Humidity: 15-85% RH (no condensation)
Storage Environment Temperature	-4°F to 140°F (-20 to 60°C), Humidity: 10-95% RH (no condensation)

Optional Accessories

Name	Order Code
Arm Cuff Cover (Pack of 5)	AX-134005759-S
Dedicated Stand	TM-ST520
Adjustable Height Stool	TM-STA001
Printer Paper (Pack of 5)	AX:PP147-S
Disposable Kiosk Sleeves (Pack of 400)	AS-134010367

Included in the Box

- 1 Kiosk Monitor
- 1 Power Cable
- 1 Printer Paper Roll
- 1 Instruction Manual
- 1 Patient Instruction Panel

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	UPC
Individual Unit Box	13.3 x 22.8 x 18.1" (34 x 60 x 46 cm)	19.55 lb (8.9 kg)	0-93764-60350-1
Kiosk Dimensions	9.5 x 18.8 x 12.8" (24 x 48 x 33 cm)	12.4 lb (5.6 kg)	N/A
Stand Dimensions	4.63 x 25.38 x 22.38" (12 x 65 x 57 cm)	18.9 lb (8.6 kg)	N/A