

UC-352BLE

PREMIUM WIRELESS WEIGHT SCALE



KEY FEATURES:



A premium wireless weight scale with capacity up to 450 pounds - automatically records and tracks measurements via mobile app to enable trending and sharing.

Trusted Accuracy

- · 4 precision sensors to ensure accurate readings
- Motion technology for accurate measurements even with slight
 movement
- · Measurements are recorded in pounds or kilograms
- High capacity of 450 lbs (200 kg)

User Friendly

- · Step on start no buttons to push
- · LCD display for easy reading
- Stores 90 measurements, plus unlimited readings on the A&D Heart Track app
- Send your measurements wirelessly via Bluetooth to your supported phone/tablet

 Free companion A&D Heart Track app automatically stores readings, tracks progress with simple graphs and allows you to share results with others

Convenience

- · Date and time stamp sync with app to assist in trending
- · Easy to share with family or healthcare providers
- Auto off to save battery life





Technical Specifications

Maximum capacity	450 lb (200 kg)			
Internal Memory	90 measurements			
Minimum display	0.2 lb units from 0 to 220.0 lb			
	0.5 lb units from 220.0 to 450 lb			
	0.1 kg units from 0 to 100.0 kg			
	0.2 kg units from 100.0 to 200 kg			
Wireless communication	VZ (MURATA Manufacturing Co.Ltd.) Bluetooth Ver.4.0LE			
Setting Items	Units: Ib or kg			
Power source	Four AAA size alkaline batteries			
Battery life	Approx. 1 year (two measurements per day when alkaline batteries are used)			
Dimensions	11.8 x 11.8 x 1.2" (300 x 300 x 31 mm)			
Weight	4.7 lb (2.1 kg) including batteries			

Packaging and Shipping Specifications

	Dimensions (WxDxH)	Weight	Individual Item Cube	UPC / GTIN
UC-352BLE Individual Unit (Box)	13.6 x 1.9 x 13.6" (34.5 x 4.8 x 34.5 cm)	5.3 lb (2.4 kg)	0.21 ft³ (0.0059 m³)	0-93764-60333-4
UC-352BLE Case (5 units)	14 x 10.2 x 14.6" (36 x 26 x 37 cm)	27.56 lb (12.5 kg)	1.2 ft³ (0.034 m³)	1-00-93764-60333-1

What is In The Box

- 1 Weight Scale
- 1 Instruction manual
- 1 Quick start guide
- 4 AAA batteries



4622 Runway Boulevard Ann Arbor, MI 48108 1-888-726-9965 andmedical.com