

ESSENTIAL UPPER ARM BLOOD PRESSURE MONITOR



CLINICALLY
VALIDATED

KEY FEATURES:



One Button Operation



Irregular Heartbeat
Detection



Blood Pressure Rating
Indicator



60 Reading Memory



5 Year Warranty



Cuff
8.6-16.5" (22-42 cm)

A classic blood pressure monitor with the benefits of speed and accuracy at a touch of a button.

Trusted Accuracy

- Clinically validated for accuracy
- Indicates an irregular heartbeat (IHB) during a measurement
- Pressure Rating Indicator classifies blood pressure in accordance with medically approved standards

User Friendly

- Easy one button operation for quick and accurate measurement
- Stores 60 measurements in memory with recall
- Displays average of all measurements
- Gentle Inflation enables a more comfortable measurement

The A&D Promise

- 5 year warranty on the monitor / 2 years on the cuff
- Cuffs are not made with natural rubber latex
- Battery life indicator
- Download the A&D Heart Track App for Free to manually log and electronically track measurements (available on iOS and Android)

Technical Specification

Measurement Method	Oscillometric
Memory	60 measurements
Cuff Size	8.6-16.5" (22-42 cm)
Measurement Range	Pressure: 0-299 mmHg Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg Pulse: 40-180 beats/min
Measurement Accuracy	Pressure: ±3 mmHg Pulse: ±5%
Power Supply	4 x 1.5V batteries (R6P, LR6 or AA)
Battery Life	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries)
Monitor Dimensions	3.8 x 5.1 x 2.7" (9.6 x 13 x 6.8 cm)
Monitor Weight	8.5 oz (240 g) without batteries
Display Size	1.3 x 2.1" (3.4 x 5.3 cm)
Clinical Test	ANSI / AAMI SP-10 1992
Operating Conditions	+50°F to +104°F (10°C to 40°C) / 15% to 85% RH / 800 to 1060 hPa
Storage / Transport Conditions	-4°F to +140°F (-20°C to 60°C) / 10% to 95% RH / 700 to 1060 hPa

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Case Cube	UPC / GTIN
Individual Unit Box	4.9 x 3.5 x 7.1" (12.5 x 9 x 18 cm)	1.14 lb (0.52 kg)	0.07 ft ³ (0.002 m ³)	0-93764-60317-4
Case (10 units)	18.7 x 10.5 x 7.6" (47.5 x 26.7 x 19.4 cm)	12.5 lb (5.7 kg)	0.87 ft ³ (0.03 m ³)	1-00-93764-60317-1

What is In The Box

- 1 Monitor
- 1 Blood pressure cuff (8.6-16.5")
- 1 Instruction Manual
- 1 Quick Start Card