

WIDE BASE WEIGHT SCALE

UC-356BLE



*Scale tray sold separately

KEY FEATURES:

-  Wireless
-  550 lb/ 249 kg Capacity
-  Extra Wide Base
-  Anti-Tip
-  Motion Tolerance
-  Smart Connect
-  2 year warranty

A wide base wireless weight scale with capacity up to 550 pounds - automatically records and tracks measurements via mobile app to enable trending and sharing.

Trusted Accuracy

- 4 precision sensors to ensure accurate readings
- Motion technology for accurate measurements even with slight movement
- Measurements are recorded in pounds or kilograms
- High capacity of 550 lbs (249 kg)
- Optional scale tray to ensure accuracy on soft surfaces such as carpet

User Friendly

- Step on start – no buttons to push
- LCD display for easy reading
- Stores 30 measurements
- Auto off to save battery life
- Audio and visual indicator for weight completion

The A&D Promise

- Professional grade performance you can trust
- 2 year warranty



Heart
Track

POWERED BY MMHG



Technical Specifications

Maximum capacity	550 lb (249 kg)
Internal Memory	63 measurements
Minimum display	33 to 220 lb with an increment of 0.2 lb (15 to 100 kg with an increment of 0.1 kg) 220 to 550 lb with an increment of 0.5 lb (100 to 250 kg with an increment of 0.2 kg)
Wireless communication	Bluetooth Ver.4.0LE
Setting Items	Units: lb or kg
Power source	Four AA size alkaline batteries
Battery life	Approx. 1 year (two measurements per day when alkaline batteries are used)
Dimensions (HxWxD)	1 3/5 x 16 1/3 x 14 1/5" (40.6 x 414.8 x 360.7 mm)
Weight	7.7 lb (3.5 kg)
Data Outputs via Integration	Time, date, weight value, battery status, set memory

Packaging and Shipping Specifications

	Dimensions (WxDxH)	Weight	Individual Item Cube	UPC / GTIN
UC-356BLE Individual Unit (Box)	16 x 3 x 18" (41 x 8 x 43 cm)	8.4 lb (3.8 kg)	0.47 ft ³ (0.0133 m ³)	0-93764-60414-0
UC-356BLE Case (3 units)	19 x 17 x 8.5" (48 x 43 x 23 cm)	25.5 lb (11.6 kg)	3 ft ³ (0.085m ³)	1-00-937646-0414-7

Optional Accessories

Name	Order Code
Scale Tray For UC-356BLE (15.7x19.5x1.8")	AX-356BLE-ST

Included in the Box

- 1 Weight Scale
- 1 Instruction manual
- 4 AA batteries