

ESSENTIAL Manual Inflate Medium or Large CuffUpper Arm Blood Pressure Monitor



KEY FEATURES:

((C)) Irregular Heartbeat Detection

M 30

30 Reading Memory

(3)

Blood Pressure Rating Indicator

AA Battery Included

Extended 5 Year Warranty

F

UA-705V: Medium Cuff 9.4-14.2" (24-36 cm) UA-705VL: Large Cuff

14.2-17.7" (36-45 cm)

A classic manual inflation blood pressure monitor with accuracy at a touch of a button.

Quick and Easy Measurements

- Easy one button operation for quick and accurate measurement
- · Stores 30 measurements in memory
- · Displays average of all measurements

Trusted Accuracy

- · Clinically validated for accuracy
- · Indicates an irregular heartbeat (IHB) during a measurement
- Pressure Rating Indicator classifies blood pressure in accordance with medically approved standards

User Comfort

- · Flexible cuff options
 - UA-705V: Medium cuff included (9.4-14.2" (24-36 cm))
 - UA-705VL: Large cuff included (14.2-17.7" (36-45 cm))
- · Manual inflation enables a comfortable measurement

The A&D Promise

- Extended 5 year monitor warranty and 2 year cuff warranty
- · Cuffs not made with natural rubber latex
- · Battery life indicator
- Download the A&D Heart Track App for Free to manually log and electronically track measurements (available on iOS and Android)

Technical Specifications

Measurement Method	Oscillometric			
Memory	30 Reading Memory			
Cuff Size	UA-705V: Medium cuff 9.4-14.2" (24-36 cm) UA-705VL: Large cuff 14.2-17.7" (36-45 cm)			
Measurement Range	Pressure: 0 - 299 mmHg Systolic pressure: 60 - 279 mmHg Diastolic pressure: 40 - 200 mmHg Pulse: 40 - 180 beats/min			
Measurement Accuracy	Pressure: ±3mmHg or 2%, whichever is greater Pulse: ±5%			
Power Supply	1 type "AA" (1.5 volt) battery (included)			
Battery Life	Approx. 2,000 measurements			
Monitor Dimensions	3.2 x 4.2 x 2.2 in (8.1 x 10.5 x 5.4 cm)			
Monitor Weight	4.3 oz. (120g) without batteries			
Display Size	Digital: Pressure - 15mm character height; Pulse - 10mm character height. Pressure and pulse displayed simultaneously.			
Clinical Test	BHS Revised Protocol: 1993			
Operating Conditions	50°F to 104°F (10°C to 40°C)			
Storage/Transport Conditions	14°F to 140°F (-10°C to 60°C)			

This product meets the Validated Device Listing ("VDL") criteria (as listed on validatebp.org) for validation of clinical accuracy, based on the independent review and acceptance of documentation submitted by the manufacturer.

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
UA-705V Individual Unit (Box)	5.06 x 3.5 x 7.69" (13 x 9 x 19 cm)	1.32 lb (0.6 kg)	0.08 ft³ (0.0022 m³)	0 93764 60157 6
UA-705V Case (10 units)	18.5 x 11 x 8.5" (47 x 28 x 22 cm)	14.34 lb (6.5 kg)	1.00 ft³ (0.028 m³)	1 00 93764 60157 3
UA-705VL Individual Unit (Box)	5.06 x 3.5 x 7.69" (13 x 9 x 19 cm)	1.32 lb (0.6 kg)	0.08 ft³ (0.0022 m³)	0 93764 60158 3
UA-705VL Case (10 units)	18.5 x 11 x 8.5" (47 x 28 x 22 cm)	14.34 lb (6.5 kg)	1.00 ft³ (0.028 m³)	1 00 93764 60158 0

What is In The Box

- 1 Upper arm blood pressure monitor
- 1 Blood pressure cuff
 - UA-705V: Medium cuff 9.4-14.2" (24-36 cm)
 - UA-705VL: Large cuff 14.2-17.7" (36-45 cm)
- 1 Latex free bulb
- 1 Trilingual instruction manual
- · 1 Quick start card
- 1 AA battery

