







PREMIUM UPPER ARM BLOOD PRESSURE MONITOR



CLINICALLY
VALIDATED

KEY FEATURES:

-  Pre-Formed Cuff
9.0-17" (22-43 cm)
-  Reading Reminder
-  Irregular Heartbeat
Detection
-  Blood Pressure Rating
Indicator
-  280 Reading Memory
-  5 Year Warranty

Quick response blood pressure monitor with advanced technology and an easy to use pre-formed cuff.

Trusted Accuracy

- Clinically validated for accuracy
- Indicates an irregular heartbeat (IHB) during a measurement
- Pressure Rating Indicator classifies blood pressure in accordance with medically approved standards

User Friendly

- Pre-formed cuff for easy placement on the upper arm
- Easy, one button operation for quick and accurate measurement
- 3 programmable reminders to ensure user takes measurements at the same time each day
- Stores 280 measurements in memory
- Displays average of all measurements
- Clock display with time and date stamp for logging and trending measurements
- Fits most arm sizes (9-17 cm)
- Gentle inflation enables a more comfortable measurement

- Large display is easy to read
- Includes an AC adapter
- Carrying case included

The A&D Promise

- Professional grade performance you can trust
- 5 year warranty on monitor, 2 year warranty on cuff
- Cuffs are not made with natural rubber latex
- Battery life indicator
- Download the A&D Heart Track App for Free to manually log and electronically track measurements (available on iOS and Android)

Technical Specification

Measurement Method	Oscillometric
Memory	280 measurements
Cuff Size	9.0-17" (24-43 cm)
Measurement Range	Pressure: 0-280 mmHg Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg Pulse: 40-200 beats/min
Measurement Accuracy	Pressure: ±3 mmHg Pulse: ±5 %
Power Supply	4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TB-233)
Battery Life	Approx. 6 months with a daily measurement
Monitor Dimensions	6.4 x 4.4 x 2.4" (16 x 11 x 6 cm)
Monitor Weight	10.6 oz (300 g) without batteries
Display Size	2.4 x 2.8" (6.1 x 7.1 cm)
Clinical Test	ANSI / AAMI SP-10 1992
Operating Conditions	50°F to 104°F (+10°C to +40°C) / 15% to 85% RH / 800 to 1060 hPa
Storage/Transport Conditions	-4°F to 140°F (-20°C to +60°C) / 10% to 95% RH / 700 to 1060 hPa

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
UA-787EJCN Individual Unit (Box)	5.75 x 4.5 x 6.75" (15 x 19 x 15 cm)	2.38 lbs (1.1 kg)	0.15 ft ³ (0.004 m ³)	0-93764-60174-3
UA-787EJCN Case (10 units)	32 x 12 x 11" (81 x 30 x 28 cm)	28.4 lbs (13 kg)	2.62 ft ³ (0.074 m ³)	1-00-93764-60174-0

What is In The Box

- 1 Blood Pressure Monitor
- 1 Blood Pressure Cuff (9.0-17")
- 1 Manual
- 1 Quick Start Guide
- 1 Carrying Case
- AC Adapter