

WIRELESS WEIGHT SCALE



KEY FEATURES:



Wireless



7 User Compatible + Guest



Motion Tolerance



Step on Start



400 lb Capacity



Up to 30 Reading Memory per User



2 year warranty

A wireless weight scale with capacity up to 400 pounds - automatically records and tracks measurements.

Trusted Accuracy

- · 4 precision sensors to ensure accurate readings
- · Measurements are recorded in pounds or kilograms
- · High capacity of 400 lbs (180 kg)

User Friendly

- · Step on start no buttons to push
- · LCD display for easy reading
- · Register 7 users plus a guest
- Stores 30 measurements per user, plus unlimited readings on the A&D Heart Track app
- · Auto off to save battery life
- Audio and visual indicator for weight completion

 Free companion A&D Heart Track app automatically stores readings, tracks progress with simple graphs and allows you to share results with others

Convenience

- · Date and time stamp sync with the app to assist in trending
- · Easy to share with family or healthcare providers
- · Auto off to save battery life





Technical Specifications

Maximum capacity	400 lb (180 kg)
Internal Memory	7 users + Guest, 30 memory each
Minimum display	0.2 lb (0.1 kg)
Setting Items	Units: lb or kg
Power source	Four AAA size alkaline batteries
Dimensions (WxLxH)	12.6 x 12.6 x 0.9" (32 x 32 x 2.4 cm)
Weight	4.4 lb (1.98 kg) including batteries

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
UC-350BLE Individual Unit (Box)	14.8 x 2 x 14.5" (38 x 5 x 37 cm)	5.6 lb (2.5 kg)	0.25 ft³ (0.0071 cm³)	0-93764-60366-2
UC-350BLE Case (8 units)	18.3 x 16.3 x 16.1" (46 x 41 x 41 cm)	52 lb (24 kg)	2.78 ft³ (0.079 cm³)	1-00-93764-60366-9

What is In The Box

- 1 Weight Scale
- 1 Instruction manual
- 1 Quick start guide
- 4 AAA batteries

