

A&D Data Manager

A PC-based software platform to track, trend, and report the data generated by Ambulatory Blood Pressure Monitors. The software includes the following features:



Graphical Analysis

Maximize the value of different views of data points or adjust for your viewing preferences with histograms, trend analysis, correlation plots, and pie graphs.



Oscillometric Wave

Verify the accuracy of a given measurement with second by second waveform.



Advanced Analytics (TM-2441 only)

Combine environmental data such as temperature, air pressure, activity, and acceleration with blood pressure data to get a more comprehensive analysis.

A&D Data Manager is a good fit for:

- PC users that can install software on an up to date PC
- Users that need access to waveform data or environmental data often used in research
- Users that need to adjust parameters or do custom reporting

SYSTEM REQUIREMENTS

ABPM Hardware: TM-2440 or TM-2441

Computer CPU: 1GHz or higher processor

Memory: 2GB or more

Operating System: Windows 10 (x64), or 11 (x64)

SVA (Recommended): 800 x 600 pixels or more

Hard Disk: 24 GB available space

A&D Heart Track ABPM

A cloud-based ABPM software solution designed by hypertension experts. It increases efficiency of the studies.



Easy Patient Set Up

Patient Education and activity logging in the Patient Diary app fosters patient compliance for a more complete and accurate study.



No Manual Calculations or Entering Data

Automatically syncs patient data and calculates technical adequacy, removing the onerous processes of traditional ABPM studies.



Fast Reporting & Filing

Interpretation IQ – Automatic interpretation in seconds. Simply review, make extra notes, and file. Data is generated by the A&D Medical ABPM (TM-2441)

It includes a set up app to program the study, a patient diary app to record their activity, and web-based analysis and reporting program.

A&D Heart Track is a good fit for:

- Users that want to work in a cloud-based environment
- Users that are interested in streamlined studies, set up and reporting
- Users that need to work at multiple sites/workstations

SYSTEM REQUIREMENTS

ABPM Hardware: TM-2441

For app: Bluetooth enabled tablet - iPad (preferred) or Samsung A7 or A8 (Android 12 and below)

For portal: Internet browser, such as Chrome or Edge