

WIRELESS UPPER ARM BLOOD PRESSURE MONITOR UA-1200BLE

KEY FEATURES:

-  Wireless
-  Smart Connect
-  Unlimited memory in app + 500 device reading memory (100/user)
-  Multi-User Capability
-  Portable All-In-One Design
-  Fit & Movement Sensors



CLINICALLY
VALIDATED

The next generation of blood pressure monitoring, providing the ultimate user experience and accuracy.

Trusted Accuracy

- Clinically validated for accuracy*
- Indicates an irregular heartbeat (IHB) during a measurement
- Body movement sensor guides user to remain still during measurement
- Cuff fit error prevents inaccurate measurements by alerting user when the cuff is too loose

User Friendly

- Easy to use for up to 5 users and a guest
- Unlimited Reading Memory In-App (500, 100/user, On-Device)
- Free highly rated A&D Heart Track app automatically stores readings, tracks progress with simple graphs and allows you to share results with others

- Send your measurements wirelessly via Bluetooth® to your supported phone/tablet
- Smart Connect automatically detects and syncs measurements between the monitor and the Heart Track app
- Blood pressure measurements can be taken without your phone or tablet – allowing increased flexibility
- Compact, all-in-one device is small, light, and perfect for travel in both upper arm and wrist configurations

The A&D Promise

- Professional grade performance you can trust
- 5 year warranty
- Cuffs are not made with natural rubber latex



*This product meets the Validated Device Listing ("VDL" criteria (as listed on validatebp.org) for validation of clinical accuracy, based on the independent review and acceptance of documentation submitted by the manufacturer

Technical specifications

Measurement Method	Oscillometric
Memory	Unlimited on App
Cuff Size	8.6-16.5" (22-42 cm) (UA-1200BLE)
Measurement Range	Pressure: 0-299 mmHg Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg Pulse: 40-180 beats / minute
Measurement Accuracy	Pressure: ± 3 mmHg Pulse: $\pm 5\%$
Power Supply	3.7V Li-ion 325mAh
Battery Life	100 measurements
Monitor Dimensions	5 x 1 x 1.7" (13 x 3 x 4.3 cm) (UA-1200BLE)
Monitor Weight	0.45 lb (204 g) (including cuff (UA-1200BLE)
Clinical Test	According to ISO81060-2
Operating Conditions	+50°F to +104°F (10°C to 40°C) / 15% to 85% RH / 800 to 1060 hPa
Storage/Transport Conditions	-4°F to +140°F (-20°C to 60°C) / 10% to 95% RH / 700 to 1060 hPa
Ingress Protection	Device: IP20

Packaging and shipping specifications

	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
UA-1200BLE Individual Unit Box	5.3 x 2.8 x 8.1 in (13.4 x 7.1 x 21 cm)	1.52 lb (0.69 kg)	0.09 ft ³ (0.0025 m ³)	0-93764-60369-3 / 1-00-93764-60369-0
UA-1200BLE Case	10.7 x 8.2 x 13.9 in (27 x 21 x 35 cm)	11.69 lb (5.3 kg)	0.70 ft ³ (0.02 m ³)	

Included in the box

- 1 Blood pressure monitor
- 1 Quick start guide
- 1 Instruction manual
- 1 USB charging cable