

WIRELESS FINGERTIP PULSE OXIMETER UP-200BLE



KEY FEATURES:



Wireless



Auto Rotate



Pulse Rate Value Display



Auto off to conserve battery life



2 Year Warranty

Wireless pulse oximeter designed to measure both pulse rate and arterial hemoglobin oxygen saturation (SpO₂). An efficient, long-lasting battery makes it the ideal choice at home or on the go. Track and share measurements wirelessly in the free, highly rated A&D Heart Track App.

Trusted Accuracy

- Fingertip photoelectric sensor measures arterial hemoglobin saturation (SpO₂), pulse rate, and pulse waveform
- Identify SpO₂ measurement trends with graphs and averages

User Friendly

- Connect wirelessly with the A&D Heart Track app to record and save SpO₂ readings
- Record measurement with just one button
- View color display from any angle with auto screen rotation
- Record measurements anywhere with long lasting power
- Convenient carrying case for storage and travel

The A&D Promise

- Professional grade performance since 1977
- Established product quality supported with a two-year warranty



Technical Specifications

Storage Environment	Temperature: -40°C - +60°C (-40°F - +140°F) Relative humidity: <=95% Atmospheric pressure: 500hPa -1060hPa
Operating Environment	Temperature: 10°C - 40°C (50°F - +104°F) Relative humidity: <=75% Atmospheric pressure: 700hPa -1060hPa
Measurement of SpO2	Measurement Range: 0 -100% Accuracy: 70-100%, ±2%; Below 70% unspecified
Measurement of Pulse Rate	Measurement Range: 30 bpm -250 bpm Accuracy: ±2 bpm or +/- 2% (whichever is greater)
Resolution	SpO2: 1%, Pulse rate: 1 bpm
Optical Sensor	Red light (wavelength is 660nm, 6.65mW) Infrared (wavelength is 905nm, 6.75mW)
Power Requirements	2 x 1.5V AAA alkaline battery (3.0V)
Power Consumption	Smaller than 100 mA
Dimensions (WxLxH)	1.3 x 2.3 x 1.3" (3.2 x 5.8 x 3.4 cm)
Weight	2.0 oz (52 g) with the batteries

Packaging and Shipping Specifications

	Quantity	Dimensions (WxLxH)	Weight	Cube	UPC / GTIN
Individual Box	1	4.3 x 1.6 x 6.3" (11 x 4 x 16 cm)	5.8 oz (0.17 kg)	0.02 ft ³ (0.0007 m ³)	0-93764-60428-7
Inner Case	10	13 x 5 x 9" (34 x 13 x 23 cm)	4.0 lbs (1.8 kg)	0.36 ft ³ (0.01 m ³)	1-00-93764-60428-4
Outer Case	100	27 x 15 x 19" (68 x 37 x 49 cm)	42 lbs (19 kg)	4.35 ft ³ (0.12 m ³)	2-00-93764-60428-1

What is In The Box

- 1 Pulse Oximeter
- 1 Instruction Manual
- 1 Carrying Case
- 1 Lanyard
- 2 AAA Batteries