

TALKING + UPPER ARM BLOOD PRESSURE MONITOR



CLINICALLY
VALIDATED

KEY FEATURES:



TriCheck



Coaching



Irregular Heartbeat
Detection



Blood Pressure Rating
Indicator



90 Reading Memory



5 Year Warranty

A premium talking blood pressure monitor with advanced inflation technology that provides audible readings and instructions in multiple languages.

Trusted Accuracy

- Clinically validated for accuracy
- TriCheck™ mode averages 3 consecutive readings per medically recommended guidelines
- Irregular Heartbeat (IHB) alerts the user to irregular heartbeat patterns
- Pressure Rating Indicator classifies blood pressure in accordance with the medically approved standards
- Cuff fit error prevents inaccurate measurements by alerting user when the cuff is too loose
- Body movement sensor guides user to remain still during measurement

User Friendly

- Talking multi-lingual device gives verbal instructions and results
- Easy, one button operation for quick and accurate measurement
- Stores 90 measurements in memory

- Easy to use medium cuff (23-37 cm) that conforms to the user's arm to provide a comfortable experience
- Adjustable pressure settings
- Advanced inflation and deflation technology providing faster and more comfortable readings
- Clock display with date and time stamp for logging and trending measurements
- Large display is easy to read
- Includes an AC adapter
- Carrying case included

The A&D Promise

- Professional grade performance you can trust
- 5 year warranty on monitor, 2 year warranty on cuff
- Cuffs are not made with natural rubber latex
- Battery life indicator
- Download the A&D Heart Track App for Free to manually log and electronically track measurements (available on iOS and Android)

Technical Specification

Measurement Method	Oscillometric
Memory	90 measurements
Cuff Size	23-37 cm (9.0-14.6")
Measurement Range	Pressure: 0-280 mmHg Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg Pulse: 40-200 beats/min
Measurement Accuracy	Pressure: ±3 mmHg Pulse: ±5 %
Power Supply	4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TB-233)
Battery Life	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries)
Monitor Dimensions	14 x 10 x 6 cm (5.5 x 4.1 x 2.3")
Monitor Weight	300 g (10.6 oz) without batteries
Display Size	6.1 x 7.1 cm (2.4 x 2.8")
Clinical Test	ANSI / AAMI SP-10 1992
Operating Conditions	+10°C to +40°C (50°F to 104°F) / 15% to 85% RH / 800 to 1060 hPa
Storage/Transport Conditions	-20°C to +60°C (-4°F to 140°F) / 10% to 95% RH / 700 to 1060 hPa

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
UA-1030TCN Individual Unit (Box)	16.5 x 9.5 x 21 cm (6.5 x 3.8 x 8.2")	953 g (2.1 lbs)	0.004 m ³ (0.14 ft ³)	0-93764-60294-8
UA-1030TCN Case (10 units)	53 x 38 x 24 cm (21 x 15 x 9.5")	10.2 kg (22.6 lbs)	0.05 m ³ (1.73 ft ³)	1-00-93764-60294-5

What is In The Box

- 1 Blood Pressure Monitor
- 1 Blood Pressure Cuff (23-37 cm)
- 1 Manual
- 1 Quick Start Guide
- 1 Carrying Case
- 4 AA Batteries
- AC Adapter

Wholesaler

	McKesson	Kohl & Frisch	Jean Coutu	Familiprix	McMahon
Codes	17569	121175	560590	412023	27020901

Health Canada License - 87256