



WIRELESS UPPER ARM BLOOD PRESSURE MONITOR UA-1200CNBLE









KEY FEATURES:



Wireless



Smart Connect



Unlimited memory in app + 500 device reading memory (100/user)



Multi-User Capability



Portable All-In-One Design



Fit & Movement Sensors

The next generation of blood pressure monitoring, providing the ultimate user experience and accuracy.

Trusted Accuracy

- · Clinically validated for accuracy
- · Indicates an irregular heartbeat (IHB) during a measurement
- Body movement sensor guides user to remain still during measurement
- Cuff fit error prevents inaccurate measurements by alerting user when the cuff is too loose

User Friendly

- · Easy to use for up to 5 users and a guest
- Unlimited reading memory in-app (500, 100/user, on-device)
- Free highly rated A&D Heart Track app automatically stores readings, tracks progress with simple graphs and allows you to share results with others

- Send your measurements wirelessly via Bluetooth® to your supported phone/tablet
- Smart Connect automatically detects and syncs measurements between the monitor and the Heart Track app
- Blood pressure measurements can be taken without your phone or tablet – allowing increased flexibility
- Compact, all-in-one device is small, light, and perfect for travel

The A&D Promise

- · Professional grade performance you can trust
- 5 year warranty
- · Cuffs are not made with natural rubber latex







Technical specifications

Measurement Method	Oscillometric			
Memory	Unlimited on App			
Cuff Size	8.6-16.5" (22-42 cm)			
Measurement Range	Pressure: 0-299 mmHg Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg Pulse: 40-180 beats / minute			
Measurement Accuracy	Pressure: ±3 mmHg Pulse: ±5%			
Power Supply	3.7V Li-ion 325mAh			
Battery Life	100 measurements			
Monitor Dimensions	5 x 1 x 1.7" (13 x 3 x 4.3 cm)			
Monitor Weight	0.45 lb (204 g)			
Clinical Test	According to ISO81060-2			
Operating Conditions	+50°F to +104°F (10°C to 40°C) / 15% to 85% RH / 800 to 1060 hPa			
Storage/Transport Conditions	-4°F to +140°F (-20°C to 60°C) / 10% to 95% RH / 700 to 1060 hPa			
Ingress Protection	Device: IP20			

Packaging and shipping specifications

	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
UA-1200CNBLE Individual Unit Box	5.2 x 2.8 x 8.1 in (13.2 x 20.6 x 7.1 cm)	.99 lb (0.45 kg)	0.07 ft³ (0.002 m³)	0-93764-60369-3 / 1-00-93764-60369-0
UA-1200CNBLE Case	15.2 x 11.2 x 9.1 in (38.6 x 28.4 x 23.1 cm)	10.7 lb (4.85 kg)	0.90 ft³ (0.025 m³)	

Included in the box

- 1 Blood pressure monitor
- 1 Quick start guide
- 1 Instruction manual
- 1 USB charging cable

