

WIRELESS UPPER ARM BLOOD PRESSURE MONITOR UA-651BLE



KEY FEATURES:



Wireless



One button operation



Irregular Heartbeat
Detection



Unlimited memory on app



5 year warranty



CLINICALLY
VALIDATED

Simple and accurate blood pressure measurement at the touch of a button including a seamless Bluetooth® connection to the highly rated Heart Track app.

Trusted Accuracy

- Clinically validated* for accuracy
- Irregular Heartbeat (IHB) alerts the user to irregular heartbeat patterns

User Friendly

- Easy, one button operation for quick and accurate measurement
- Send your measurements wirelessly via Bluetooth to your supported phone/tablet
- Free companion A&D Heart Track app automatically stores readings, tracks progress with simple graphs and allows you to share results with others
- Data integration with Google Fit and Apple Health Kit
- Unlimited reading memory in the Heart Track app
- Easy to apply wide range cuff that comfortably fits most arms (8.6-16.5" / 22-42 cm)
- Compatible with small to large cuffs (not included)

The A&D Promise

- Professional grade performance you can trust
- 5 year warranty on monitor, 2 year warranty on cuff
- Cuffs are not made with natural rubber latex



Heart
Track

POWERED BY MMHG



Compatible with:



*This product meets the Validated Device Listing ("VDL") criteria (as listed on validatebp.org) for validation of clinical accuracy, based on the independent review and acceptance of documentation submitted by the manufacturer.

Technical Specifications

Measurement Method	Oscillometric
Memory	Unlimited in the Heart Track App
Cuff Size	8.6-16.5" / 22-42 cm
Measurement Range	Pressure: 0-299 mmHg, Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg, Pulse: 40-180 beats / minute
Measurement Accuracy	Pressure: ±3 mmHg, Pulse: ±5%
Power Supply	4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TB-233)
Number of Measurements	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries)
Monitor Dimensions	3.8 x 2.7 x 5.1" / 9.6 x 6.8 x 13 cm
Monitor Weight	8.8 oz (249 g) without batteries
Display Size	1.4 x 2.1" (3.6 x 5.3 cm)
Clinical Test	ISO81060-2: 2013
Operating Conditions	+50°F to +104°F (10°C to 40°C) / 15% to 85% RH / 800 to 1060 hPa
Storage/Transport Conditions	-4°F to +140°F (-20°C to 60°C) / 10% to 95% RH / 700 to 1060 hPa
Included In the Box	<ul style="list-style-type: none"> • 1 Upper arm blood pressure monitor • 1 Blood pressure cuff (8.6-16.5") • 1 Instruction manual • 1 Quick start guide • 1 Carrying case • 4 AA batteries

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
UA-651BLE Individual Unit (Box)	5.1 x 4 x 7.7" (13 x 10 x 19 cm)	1.52 lb (0.7 kg)	0.09 ft³ (0.0024 m³)	0-93764-60332-7
UA-651BLE Case (10 units)	21.7 x 11 x 8.5" (55 x 28 x 22 cm)	15.87 lb (7.2 kg)	1.18 ft³ (0.033 m³)	1-00-93764-60332-4

Optional Accessories

Name	Order Code	Weight	Dimensions	UPC / GTIN
AC Power Adapter	TB:233	2.4 oz	4.5 x 1.5 x 3.5"	0-93764-60235-1
Compatible Replacement Cuffs- Packaged in individual box				
Small Cuff 6.3-9.4" (16-24cm)	UA-289A	4.4 oz	4.7 x 1.5 x 7.5"	0-93764-60416-4
Medium Cuff 9.0-14.6" (23-37 cm)	UA-290A	6.3 oz	4.7 x 1.5 x 7.5"	0-93764-60417-1
Large Cuff 12.2-17.7" (31-45 cm)	UA-291A	8.7 oz	4.7 x 1.5 x 7.5"	0-93764-60418-8
Wide Range Cuff 8.6-16.5" (22-42 cm)	UA-420A	7.3 oz	4.7 x 1.5 x 7.5"	0-93764-60419-5