



ESSENTIAL UPPER ARM BLOOD PRESSURE MONITOR







KEY FEATURES:



One Button Operation



AFIB Indicator



Blood Pressure Rating Indicator



60 Reading Memory



5 Year Warranty



Cuff

22-42 cm (8.6-16.5")

A classic blood pressure monitor with the benefits of speed and accuracy at a touch of a button.

Trusted Accuracy

- · Clinically validated for accuracy
- Alerts the user to the presence of an irregular heartbeat (IHB) / Atrial Fibrillation (AFib)
- Pressure Rating Indicator classifies blood pressure in accordance with medically approved standards

User Friendly

- Easy one button operation for quick and accurate measurement
- · Stores 60 measurements in memory with recall
- Displays average of all measurements
- · Gentle Inflation enables a more comfortable measurement

The A&D Promise

- 5 year warranty on the monitor / 2 years on the cuff
- · Cuffs are not made with natural rubber latex
- · Battery life indicator
- Download the A&D Heart Track App for Free to manually log and electronically track measurements (available on iOS and Android)

Technical Specifications

Measurement Method	urement Method Oscillometric				
Memory	60 measurements				
Cuff Size	22-42 cm (8.6-16.5")				
Measurement Range	Pressure: 0-299 mmHg Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg Pulse: 40-180 beats / minute				
Measurement Accuracy	Pressure: ±3 mmHg Pulse: ±5%				
Power Supply	4 x 1.5V batteries (R6P, LR6 or AA) or AC Adapter sold seperately (TB:233)				
Battery Life	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries)				
Monitor Dimensions	9.6 x 13 x 6.8 cm (3.8 x 5.1 x 2.7")				
Monitor Weight	240 g (8.5 oz) without batteries				
Display Size	ze 3.4 x 5.3 cm (1.3 x 2.1")				
Clinical Test	Test ANSI / AAMI SP-10 1992				
Operating Conditions	+10°C to +40°C (50°F to 104°F) / 15% to 85% RH / 800 to 1060 hPa				
torage/Transport Conditions -20°C to +60°C (-4°F to 140°F) / 10% to 95% RH / 700 to 1060 hPa					

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
UA-651CN Individual Unit (Box)	12.5 x 9 x 18 cm (4.9 x 3.5 x 7.1")	0.6 kg (1.31 lb)	0.002 m³ (0.07 ft³)	0-93764-60313-6
UA-651CN Case (10 units)	47.5 x 26.7 x 19.4 cm (18.7 x 10.5 x 7.6")	6.47 kg (14.3 lb)	0.03 m³ (0.87 ft³)	1-00-93764-60313-3

What is In The Box

- 1 Blood pressure monitor
- 1 Blood pressure cuff (22-42 cm)
- 1 Instruction manual
- 1 Quick start guide
- 4 AA Batteries

Wholesaler

	McKesson	Kohl & Frisch	Jean Coutu	Familiprix	McMahon	LPG
Codes	058574	147149	560516	009672	29167801	035249

Health Canada License - 92815

