



WIRELESS UPPER ARM BLOOD PRESSURE MONITOR





KEY FEATURES:



Wireless



One button operation



Irregular Heartbeat
Detection



Unlimited memory on app



5 year warranty



A wireless blood pressure monitor that automatically records and tracks measurements via mobile app to enable trending and sharing.

Trusted Accuracy

- · Clinically validated for accuracy
- Irregular Heartbeat (IHB) alerts the user to irregular heartbeat patterns

User Friendly

- Easy, one button operation for quick and accurate measurement
- Send your measurements wirelessly via Bluetooth to your supported phone/tablet
- Free companion A&D Heart Track app automatically stores readings, tracks progress with simple graphs and allows you to share results with others
- · 30 reading memory in off-line mode, unlimited on app
- Easy to apply medium cuff that comfortably fits most arms (22-42 cm)
- · Gentle inflation enables a more comfortable measurement

The A&D Promise

- Professional grade performance you can trust
- Simple integration to telehealth systems with full service technical support
- 5 year warranty on monitor, 2 year warranty on cuff
- · Cuffs are not made with natural rubber latex
- · Battery life indicator







Technical Specifications

Measurement Method	Oscillometric				
Memory	Unlimited on App				
Cuff Size	22-42 cm (8.6-16.5")				
Measurement Range	Pressure: 0-299 mmHg Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg Pulse: 40-180 beats / minute				
Measurement Accuracy	Pressure: ±3 mmHg Pulse: ±5%				
Power Supply	4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TB-233)				
Battery Life	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries)				
Monitor Dimensions	9.6 x 13 x 6.8 cm (3.8 x 5.1 x 2.7")				
Monitor Weight	249 g (8.8 oz) without batteries				
Display Size	3.6 x 5.3 cm (1.4 x 2.1")				
Clinical Test	ISO81060-2: 2013				
Operating Conditions	+10°C to +40°C (50°F to 104°F) / 15% to 85% RH / 800 to 1060 hPa				
Storage/Transport Conditions	-20°C to +60°C (-4°F to 140°F) / 10% to 95% RH / 700 to 1060 hPa				

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
UA-651CNBLE Individual Unit (Box)	13 x 10 x 19 cm (5.1 x 4 x 7.7")	0.7 kg (1.52 lb)	0.0024 m³ (0.09 ft³)	0-93761-60348-8
UA-651CNBLE Case (10 units)	55 x 28 x 22 cm (21.7 x 11 x 8.5")	7.2 kg (15.87 lb)	0.033 m³ (1.18 ft³)	1-00-93764-60348-5

What is In The Box

- 1 Blood pressure monitor
- 1 Blood pressure cuff (22-42 cm)
- 1 Instruction manual
- 1 Quick start guide
- 1 Carrying case
- · 4 AA batteries

Wholesaler

	McKesson	Kohl & Frisch	Jean Coutu	Familiprix	McMahon	LPG
Codes	82959	152972	709096	037797	30525701	035250

Health Canada License - 94858

