

# PREMIUM UPPER ARM BLOOD PRESSURE MONITOR



CLINICALLY  
VALIDATED

## KEY FEATURES:



4 User Compatible



Snapshot Averaging



Irregular Heartbeat Detection



Pressure Rating Indication



60 Memory Recall per user



Date and Time Stamps



Body Movement Sensor



5 year warranty

Top rated, multi-user upper arm blood pressure monitor that is clinically validated for accuracy.

### Trusted Accuracy

- Clinically validated for accuracy
- Top-rated by several independent publications
- Irregular Heartbeat (IHB) alerts the user to irregular heartbeat patterns
- Pressure Rating Indicator classifies blood pressure in accordance with the medically approved standards

### User Friendly

- Simple operation for up to 4 different users with a unique 60 reading memory for each user
- Snapshot Averaging automatically calculates the average for all readings to determine the average baseline blood pressure

- AC adapter for operation without batteries
- Easy to apply wide-range cuff that comfortably fits most arms (22-42 cm)
- Gentle inflation enables a more comfortable measurement

### The A&D Promise

- Professional grade performance you can trust
- 5 year warranty on monitor, 2 year warranty on cuff
- Cuffs are not made with natural rubber latex
- Battery life indicator
- Download the A&D Heart Track App for Free to manually log and electronically track measurements (available on iOS and Android)

## Technical Specifications

<b>Measurement Method</b>	Oscillometric
<b>Memory</b>	60 measurements for four users
<b>Cuff Size</b>	22-42 cm (8.6-16.5")
<b>Measurement Range</b>	Pressure: 0-299 mmHg Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg Pulse: 40-180 beats / minute
<b>Measurement Accuracy</b>	Pressure: ±3 mmHg Pulse: ±5%
<b>Power Supply</b>	4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TB-233)
<b>Battery Life</b>	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries)
<b>Monitor Dimensions</b>	14 x 10 x 6 cm (5.5 x 4.1 x 2.3")
<b>Monitor Weight</b>	255 g (9.0 oz) without batteries
<b>Display Size</b>	4.3 x 4.8 cm (1.7 x 1.9")
<b>Clinical Test</b>	ISO81060-2: 2013
<b>Operating Conditions</b>	+10°C to +40°C (50°F to 104°F) / 15% to 85% RH / 800 to 1060 hPa
<b>Storage/Transport Conditions</b>	-20°C to +60°C (-4°F to 140°F) / 10% to 95% RH / 700 to 1060 hPa

## Packaging and Shipping Specifications

	<b>Dimensions (WxLxH)</b>	<b>Weight</b>	<b>Individual Item Cube</b>	<b>UPC / GTIN</b>
<b>UA-767FAM Individual Unit (Box)</b>	13 x 9.9 x 20 cm (5.1 x 3.9 x 7.7")	0.8 kg (1.7 lb)	0.0014 cm <sup>3</sup> (0.05 ft <sup>3</sup> )	0-93764-60347-1
<b>UA-767FAM Case (10 units)</b>	51 x 28 x 23 cm (20 x 11 x 9")	8.0 kg (17.6 lb)	0.031 cm <sup>3</sup> (1.08 ft <sup>3</sup> )	1-00-93764-60347-8

### What is In The Box

- 1 Blood pressure monitor
- 1 Blood pressure cuff (22-42 cm)
- 1 Instruction manual
- 1 Quick start guide
- 1 Carrying case
- 1 AC adapter

### Wholesaler

	<b>McKesson</b>	<b>Kohl &amp; Frisch</b>	<b>Jean Coutu</b>	<b>Familiprix</b>	<b>McMahon</b>
<b>Codes</b>	82964	152971	560084	038666	30504301

### Health Canada license - 94753