

# UA-789EJCN

# **PREMIUM UPPER ARM** BLOOD PRESSURE MONITOR



Quick response blood pressure monitor with advanced technology and an easy to use pre-formed cuff.

## **Trusted Accuracy**

- · Clinically validated for accuracy
- · Indicates an irregular heartbeat (IHB) during a measurement
- Pressure Rating Indicator classifies blood pressure in accordance with medically approved standards

#### **User Friendly**

- Pre-formed cuff for easy placement on the upper arm
- Easy, one button operation for quick and accurate
- measurement
  3 programmable reminders to ensure user takes
- measurements at the same time each dayStores 280 measurements in memory
- Displays average of all measurements
- Clock display with time and date stamp for logging and trending measurements

- Fits most arm sizes (22-43 cm)
- · Gentle inflation enables a more comfortable measurement
- · Large display is easy to read
- Includes an AC adapter
- · Carrying case included

#### **The A&D Promise**

- · Professional grade performance you can trust
- · 5 year warranty on monitor, 2 year warranty on cuff
- · Cuffs are not made with natural rubber latex
- · Battery life indicator
- Download the A&D Heart Track App for Free to manually log and electronically track measurements (available on iOS and Android)

# **UA-787EJCN**

## **Technical Specification**

Measurement Method	Oscillometric			
Memory	280 measurements			
Cuff Size	24-43 cm (9.0-17")			
Measurement Range	Pressure: 0-280 mmHg Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg Pulse: 40-200 beats/min			
Measurement Accuracy	Pressure: ±3 mmHg Pulse: ±5 %			
Power Supply	4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TB-233)			
Battery Life	Approx. 6 months with a daily measurement			
Monitor Dimensions	16 x 11 x 6 cm (6.4 x 4.4 x 2.4")			
Monitor Weight	300 g (10.6 oz) without batteries			
Display Size	6.1 x 7.1 cm (2.4 x 2.8")			
Clinical Test	ANSI / AAMI SP-10 1992			
Operating Conditions	+10°C to +40°C (50°F to 104°F) / 15% to 85% RH / 800 to 1060 hPa			
Storage/Transport Conditions	-20°C to +60°C (-4°F to 140°F) / 10% to 95% RH / 700 to 1060 hPa			

# **Packaging and Shipping Specifications**

	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
UA-787EJCN Individual Unit (Box)	15 x 19 x 15 cm (5.75 x 4.5 x 6.75")	1.1 kg (2.38 lbs)	0.004 cm³ (0.15 ft³)	0-93764-60174-3
UA-787EJCN Case (10 units)	81 x 30 x 28 cm (32 x 12 x 11")	13 kg (28.4 lbs)	0.074 cm³ (2.62 ft³)	1-00-93764-60174-0

#### What is In The Box • 1 Blood Pressure Monitor

• 1 Blood Pressure Cuff (22-43 cm)

### Wholesaler

	McKesson	Kohl & Frisch	Jean Coutu	Familiprix	McMahon
Codes	940767	87290	N/A	280081	22864901

- 1 Quick Start Guide
- 1 Carrying Case
- AC Adapter

• 1 Manual

#### Health Canada License - 30619



6695 Millcreek Drive, Unit 6 Missauga, ON L5N 5R8, Canada 1-800-461-0991 LifeSourceCanada.com