







BLOOD PRESSURE MONITOR WITH EXTRA LARGE CUFF



KEY FEATURES:

-  Extra Large Cuff
42-60 cm (16.5-23.6")
-  3 Memory Average
-  Reading Reminder
-  Irregular Heartbeat
Detection
-  60 Reading Memory
-  5 Year Warranty

The only home blood pressure monitor designed to provide accurate readings for people with larger arms.

Trusted Accuracy

- Clinically validated for accuracy
- Indicates an Irregular Heartbeat during a measurement
- Pressure Rating Indicator classifies blood pressure in accordance with medically approved standards

User Friendly

- Stores 60 measurements in memory
- Averages the last 3 readings
- Clock display with time and date stamp for logging and trending measurements
- 3 programmable reminders to ensure user takes measurements at the same time each day

- Extra large cuff to accommodate people with larger arms (42-60 cm)
- Fast and accurate readings with the high speed motor

The A&D Promise

- 5 year warranty on the monitor/ 2 years on the cuff
- Cuffs are not made with natural rubber latex
- Battery life indicator
- Download the A&D Heart Track App for Free to manually log and electronically track measurements (available on iOS and Android)

Technical Specification

Measurement Method	Oscillometric
Memory	60 measurements
Cuff Size	42-60 cm (16.5-23.6")
Measurement Range	Pressure: 0-280 mmHg Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg Pulse: 40-200 beats /min
Measurement Accuracy	Pressure: ±3 mmHg Pulse: ±5 %
Power Supply	AC adapter (included) or 4 x 1.5V batteries (R6P, LR6 or AA)
Battery Life	Approximately 120 measurements
Monitor Dimensions	16.3 x 11.2 x 6.2 cm (6.4 x 4.4 x 2.4")
Monitor Weight	350 g (12.3 oz)
Display Size	5.5 x 7 cm (2.2 x 2.8")
Clinical Test	According to ANSI/AAMI SP-10 2002
Operating Conditions	+10°C to +40°C (50°F to 104°F) / 15% to 85% RH / 800 to 1060 hPa
Storage/Transport Conditions	-20°C to +60°C (-4°F to 140°F) / 10% to 95% RH / 700 to 1060 hPa

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
UA-789AC Individual Unit (Box)	17 x 11 x 21 cm (6.5 x 4.0 x 8.3")	1.0 kg (2.2 lb)	0.0036 cm ³ (0.13 ft ³)	0-093764-60062-3
UA-789AC Case (10 units)	55 x 35 x 23 cm (22 x 14 x 9")	16.7 kg (36.8 lb)	0.045 cm ³ (1.59 ft ³)	1-00-93764-60062-0

What is In The Box

- 1 Blood pressure monitor
- 1 Blood pressure cuff (42-60 cm)
- 1 Instruction manual
- 2 Quick start guides
- 1 Carrying case
- 1 AC adapter

Wholesaler

	McKesson	Kohl & Frisch	Jean Coutu	Familiprix	McMahon	LPG
Codes	444489	95425	N/A	N/A	26008201	003631

Health Canada License - 73654