

ESSENTIALWRIST BLOOD PRESSURE MONITOR







KEY FEATURES:



One Button Operation



) AFIB Indicator



Blood Pressure Rating Indicator



60 Reading Memory



5 Year Warranty



Wrist Cuff 13.5-21.5 cm

A classic wrist blood pressure monitor with the benefits of speed and accuracy at a touch of a button.

Trusted Accuracy

- · Clinically validated for accuracy
- Indicates an irregular heartbeat (IHB) / Atrial fibrillation (AFib) during measurement
- Pressure Rating Indicator classifies blood pressure in accordance with World Health Organization (WHO) medically approved standards

User Friendly

- Easy one button operation for quick and accurate measurement
- · Stores 60 measurements in memory with recall

- · Displays average of all measurements
- Easy to use wrist cuff (13.5-21.5 cm)
- · Gentle inflation enables a more comfortable measurement
- · Compact size for storage and travel

The A&D Promise

- 5 year warranty
- · Cuffs are not made with natural rubber latex
- Battery life indicator
- Download the A&D Heart Track App for Free to manually log and electronically track measurements (available on iOS and Android)

Technical Specification

Measurement Method	Oscillometric			
Memory	60 measurements			
Wrist Circumference	13.5 -21.5 cm (5.3-8.5")			
Measurement Range	Pressure: 0-299 mmHg Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg Pulse: 40-180 beats /min			
Measurement Accuracy	Pressure: ±3 mmHg Pulse: ±5%			
Power Supply	2 x 1.5V alkaline batteries (LR03 or AAA)			
Battery Life	Approximately 250 measurements			
Monitor Dimensions	5.6 x 2.2 x 8.8 cm (2.2 x 0.9 x 3.5")			
Unit Weight	91 g (3.2 oz)			
Display Size	2.3 x 4 cm (0.9 x 1.6")			
Clinical Test	According to ISO81060-2:2013			
Operating Conditions	+10°C to +40°C (50°F to 104°F) / 15% to 85% RH / 800 to 1060 hPa			
-20°C to +60°C (-4°F to 140°F) / 10% to 95% RH / 700 to 1060 hPa				

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
UB-525CN Individual Unit (Box)	9.2 × 8.6 x 11.5 cm (3.6 x 3.4 x 4.5")	0.09 kg (0.19 lb)	0.00003 cm³ (0.001 ft³)	0-93764-60401-0
UB-525CN Case (20 units)	20 x 45 × 25 cm (7.9 x 17.9 x 9.9")	1.8 kg (4.05 lb)	0.00065 cm³ (0.023 ft³)	1-00-93764-60401-7

What is In The Box

- 1 Wrist blood pressure monitor
- 1 Instruction manual
- 1 Quick start guide
- · 2 AAA batteries

Wholesaler

	McKesson	Kohl & Frisch	Jean Coutu	Familiprix	McMahon
Codes	148748	163216	N/A	N/A	N/A

Health Canada License - 102407

