








WIRELESS WEIGHT SCALE



KEY FEATURES:

-  Wireless
-  7 User Compatible + Guest
-  Motion Tolerance
-  Step on Start
-  180 kg Capacity
-  Up to 30 Reading Memory per User
-  2 year warranty

A wireless weight scale with capacity up to 400 pounds - automatically records and tracks measurements.

Trusted Accuracy

- 4 precision sensors to ensure accurate readings
- Measurements are recorded in pounds or kilograms
- High capacity of 180 kg (400 lbs)

User Friendly

- Step on start – no buttons to push
- LCD display for easy reading
- Register 7 users plus a guest
- Stores 30 measurements per user, plus unlimited readings on the A&D Heart Track app
- Auto off to save battery life
- Audio and visual indicator for weight completion

- Free companion A&D Heart Track app automatically stores readings, tracks progress with simple graphs and allows you to share results with others

Convenience

- Date and time stamp sync with the app to assist in trending
- Easy to share with family or healthcare providers
- Auto off to save battery life



Technical Specifications

Maximum capacity	180 kg (400 lb)
Internal Memory	7 users + Guest, 30 memory each
Minimum display	0.1 kg (0.2 lb)
Setting Items	Units: kg or lb
Power source	Four AAA size alkaline batteries
Dimensions (WxLxH)	32 x 32 x 2.4 cm (12.6 x 12.6 x 0.9")
Weight	1.98 kg (4.4 lb) including batteries

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
UC-350BLE Individual Unit (Box)	38 x 5 x 37 cm 14.8 x 2 x 14.5"	2.5 kg (5.6 lb)	0.0071 cm ³ (0.25 ft ³)	0-93764-60377-8
UC-350BLE Case (8 units)	46 x 41 x 41 cm 18.3 x 16.3 x 16.1"	24 kg (52 lb)	0.079 cm ³ (2.78 ft ³)	

What is In The Box

- 1 Weight Scale
- 1 Instruction manual
- 1 Quick start guide
- 4 AAA batteries

Wholesaler

	McKesson	Kohl & Frisch	Jean Coutu	Familiprix	McMahon
Codes	106643	N/A	560518	072621	N/A

Health Canada License - XXXXX