


# UC-352CNBLE

# **PREMIUM WIRELESS** WEIGHT SCALE



## **KEY FEATURES:**



A premium wireless weight scale with capacity up to 450 pounds - automatically records and tracks measurements via mobile app to enable trending and sharing.

# **Trusted Accuracy**

- · 4 precision sensors to ensure accurate readings
- Motion technology for accurate measurements even with slight
  movement
- · Measurements are recorded in pounds or kilograms
- High capacity of 200 kg (450 lbs)

## **User Friendly**

- · Step on start no buttons to push
- · LCD display for easy reading
- Stores 90 measurements, plus unlimited readings on the A&D Heart Track app
- Send your measurements wirelessly via Bluetooth to your supported phone/tablet

 Free companion A&D Heart Track app automatically stores readings, tracks progress with simple graphs and allows you to share results with others

## Convenience

- · Date and time stamp sync with app to assist in trending
- · Easy to share with family or healthcare providers
- · Auto off to save battery life





# UC-352CNBLE

# **Technical Specifications**

Maximum capacity	200 kg (450 lbs)					
Internal Memory	90 measurements					
Minimum display	0.2 lb units from 0 to 220.0 lb 0.5 lb units from 220.0 to 450 lb 0.1 kg units from 0 to 100.0 kg 0.2 kg units from 100.0 to 200 kg					
Wireless communication	VZ (MURATA Manufacturing Co.Ltd.) Bluetooth Ver.4.0LE					
Setting Items	Units: kg or lb					
Power source	Four AAA size alkaline batteries					
Battery life	Approx. 1 year (two measurements per day when alkaline batteries are used)					
Dimensions	300 x 300 x 31 mm (11.8 x 11.8 x 1.2")					
Weight	2.1 kg (4.7 lb) including batteries					

# **Packaging and Shipping Specifications**

	Dimensions (WxDxH)	Weight	Individual Item Cube	UPC / GTIN
UC-352BLE Individual Unit (Box)	34.5 x 4.8 x 34.5 cm (13.6 x 1.9 x 13.6")	2.4 kg (5.3 lb)	0.0059 cm <sup>3</sup> (0.21 ft <sup>3</sup> )	
UC-352BLE Case (5 units)	36 x 26 x 37 cm (14 x 10.2 x 14.6")	12.5 kg (27.56 lb)	0.034 cm³ (1.2 ft³)	

## What is In The Box

- 1 Weight Scale
- 1 Instruction manual
- 1 Quick start guide
- 4 AAA batteries

#### Wholesaler

	McKesson	Kohl & Frisch	Jean Coutu	Familiprix	McMahon
Codes					

## Health Canada License - XXXXX



6695 Millcreek Drive, Unit 6 Missauga, ON L5N 5R8, Canada 1-800-461-0991 LifeSourceCanada.com