









PULSE OXIMETER FOR EASY AT-HOME MONITORING



KEY FEATURES:

-  Small Size and Lightweight
-  Auto Rotate
-  Simple Operation
-  Pulse Rate Value Display
-  Low Power Consumption
-  Low Battery Indicator
-  Auto off to conserve battery life
-  2 Year Warranty

Easy to use, reliable Pulse Oximeter

The A&D Medical Pulse Oximeter is easy to use with a one button operation. Place the Pulse Oximeter on your finger, press the button, and the results are clearly shown in any direction on the LCD. The device is portable making it easy to take measurements on the go and has power saving features to prolong battery life.

The Pulse Oximeter is a non-invasive device intended for checking oxygen saturation of arterial hemoglobin (SpO₂) and the pulse rate of adults at home.

Technical Specifications

Storage Environment	Temperature: -40°C - +60°C Relative humidity: <=95% Atmospheric pressure: 500-1060 hPa
Operating Environment	Temperature: 10°C - 40°C Relative humidity: <=75% Atmospheric pressure: 700-106 0hPa
Measurement of SpO2	Measurement Range: 0-100% Accuracy: 70-100%, ±2%; Below 70% unspecified
Measurement of Pulse Rate	Measurement Range: 30-250 bpm Accuracy: ±2 bpm
Power Requirements	2 x 1.5 V AAA alkaline battery (2.6-3.6 V)
Power Consumption	Smaller than 25 mA
Resolution	SpO2: 1%, Pulse rate: 1 bpm
Optical Sensor	Red light (wavelength is 660 nm, 6.65 mW) Infrared (wavelength is 905 nm, 6.75 mW)
Dimensions (LxWxH)	5.9 x 3.3 x 3.2 cm (2.3 x 1.3 x 1.3")
Weight	48 g (1.7 oz) with the batteries

Packaging and Shipping Specifications

	Quantity	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
Individual Box	1	11 x 4 x 16 cm (4.3 x 1.6 x 6.3")	0.16 kg (5.7 oz)	0.0007 m ³ (0.02 ft ³)	0-93764-60415-7
Inner Case	10	34 x 13 x 23 cm (13 x 5 x 9")	1.8 kg (4.0 lbs)	0.01 m ³ (0.36 ft ³)	1-00-93764-60415-4
Outer Case	100	68 x 37 x 49 cm (27 x 15 x 19")	19 kg (41 lbs)	0.12 m ³ (4.35 ft ³)	2-00-93764-60415-1

Included in the Box

- 1 Pulse Oximeter
- 1 Instruction Manual
- 1 Carrying Case
- 1 Lanyard
- 2 AAA Batteries

Wholesaler

	McKesson	LPG	Kohl & Frisch	Jean Coutu	Familiprix	McMahon
Codes	170585	035561	N/A	N/A	176959	N/A

Health Canada License - 102228