







# TALKING UPPER ARM BLOOD PRESSURE MONITOR UA-1030T



## KEY FEATURES:

-  Audible Coaching
-  TriCheck Advanced Averaging
-  Irregular Heartbeat Detection
-  Blood Pressure Rating Indicator
-  90 Reading Memory
-  5 Year Warranty



CLINICALLY  
VALIDATED

Trusted accuracy and reliable blood pressure measurements every time that provides audible readings and instructions in English, Spanish and French.

### Trusted Accuracy

- Clinically validated<sup>1</sup> for accuracy
- Advanced averaging that Averages the last three readings per clinical guidelines
- Irregular Heartbeat (IHB) alerts the user to irregular heartbeat patterns
- Pressure Rating Indicator classifies blood pressure in accordance with the American Heart Association's Blood Pressure Guidelines
- Cuff fit error prevents inaccurate measurements by alerting user when the cuff is too loose
- Body movement sensor guides user to remain still during measurement

### User Friendly

- Talking multi-lingual device gives verbal instructions and results in English, French, and Spanish
- Easy, one button operation for quick and accurate measurement

- Stores 90 measurements in memory
- Easy to use medium cuff (9.0-14.6") that conforms to the user's arm to provide a comfortable experience
- Adjustable pressure settings
- Advanced inflation and deflation technology providing faster and more comfortable readings
- Clock display with date and time stamp for logging and trending measurements
- Large display is easy to read

### The A&D Promise

- Professional grade performance you can trust
- 5 year warranty on monitor, 2 year warranty on cuff
- Cuffs are not made with natural rubber latex
- Download the A&D Heart Track App for Free to manually log and electronically track measurements (available on iOS and Android)

\*This product meets the Validated Device Listing ("VDL") criteria (as listed on [validatebp.org](http://validatebp.org)) for validation of clinical accuracy, based on the independent review and acceptance of documentation submitted by the manufacturer

## Technical Specification

<b>Measurement Method</b>	Oscillometric
<b>Memory</b>	90 measurements
<b>Cuff Size</b>	9.0-14.6" (23-37 cm)
<b>Measurement Range</b>	Pressure: 0-299 mmHg Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg Pulse: 40-180 beats/min
<b>Measurement Accuracy</b>	Pressure: $\pm 3$ mmHg Pulse: $\pm 5$ %
<b>Power Supply</b>	4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TB-233)
<b>Battery Life</b>	Approx. 1000 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries)
<b>Monitor Dimensions</b>	5.5 x 4.1 x 2.3" (14 x 10 x 6 cm)
<b>Monitor Weight</b>	10.6 oz (300g) without batteries
<b>Display Size</b>	2.4 x 2.8" (6.1 x 7.1 cm)
<b>Clinical Test</b>	According to ISO81060-2: 2013
<b>Operating Conditions</b>	+10°C to +40°C (50°F to 104°F) / 15% to 85% RH / 800 to 1060 hPa
<b>Storage/Transport Conditions</b>	-20°C to +60°C (-4°F to 140°F) / 10% to 95% RH / 700 to 1060 hPa

## Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
<b>UA-1030T Individual Unit (Box)</b>	6.5 x 3.8 x 8.3" (16.6 x 9.6 x 21.1 cm)	1.8 lbs (822 g)	0.12 ft <sup>3</sup> (0.003 m <sup>3</sup> )	0-93764-60293-1
<b>UA-1030T Case (10 units)</b>	14 x 20 x 9" (35 x 51 x 23 cm)	20.3 lbs (9.2 kg)	1.46 ft <sup>3</sup> (0.04 m <sup>3</sup> )	1-00-93764-60293-8

## What is In The Box

- 1 Blood Pressure Monitor
- 1 Blood Pressure Cuff (9.0-14.6")
- 1 Manual
- 1 Quick Start Guide
- 1 Carrying Case
- 4 AA Batteries
- AC Adapter

## Optional Accessories

Name	Order Code	Weight	Dimensions	UPC/ GTIN
AC Power Adapter	TB:233	2.4 oz	4.5 x 1.5 x 3.1"	0-93764-60235-1
<b>Compatible Replacement Cuffs- Packaged in individual box</b>				
Small Cuff 6.3-9.4" (16-24cm)	UA-289A	4.4 oz	4.7 x 1.5 x 7.5"	0-93764-60416-4
Medium Cuff 9.0-14.6" (23-37 cm)	UA-290A	6.3 oz	4.7 x 1.5 x 7.5"	0-93764-60417-1
Large Cuff 12.2-17.7" (31-45 cm)	UA-291A	8.7 oz	4.7 x 1.5 x 7.5"	0-93764-60418-8
Wide Range Cuff 8.6-16.5" (22-42 cm)	UA-420A	7.3 oz	4.7 x 1.5 x 7.5"	0-93764-60419-5