

A premium talking blood pressure monitor with advanced features and verbal guidance making it the ideal unit for the visually impaired.

Product Benefits

- TriCheck™ automated protocol takes 3 sequential measurements and averages in accordance with clinical guidelines
- Cuff fit error ensures proper cuff fit and body movement detection for accurate readings
- Irregular Heartbeat (IHB) indicates the presence of an irregular heartbeat

Comfort and Convenience

- Adjustable pressure settings to ensure comfort, customize the compression to users preference
- Advanced inflation and deflation technology provides faster and more comfortable readings
- Includes AC adapter and carrying case

The A&D Promise

- Clinical Grade Accuracy you can trust
- 5 year warranty on monitor, 2 year warranty on cuff
- Cuffs are not made with natural rubber latex



90 Measurement Memory

Irregular heartbeat symbol (IHB)


Cuff Fit Error

Movement Error

Pressure Rating Indicator

Battery Life Indicator

TriCheck™



Clock Display



Hypertension CANADA

Verbal instructions (EN, SP and FR)

One Button Operation



Unlock the Potential
Download the FREE
A&D Heart Track
Companion App



Technical Specification

Measurement Method	Oscillometric
Memory	90 measurements
Cuff Size	23-37 cm (9.0-14.6")
Measurement Range	Pressure: 0-280 mmHg, Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg, Pulse: 40-200 beats/min
Measurement Accuracy	Pressure: ± 3 mmHg, Pulse: ± 5 %
Power Supply	4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TB-233)
Battery Life	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries)
Monitor Dimensions	14 x 10 x 6 cm (5.5 x 4.1 x 2.3")
Monitor Weight	300 g (10.6 oz) without batteries
Display Size	6.1 x 7.1 cm (2.4 x 2.8")
Clinical Test	ANSI / AAMI SP-10 1992
Operating Conditions	+10°C to +40°C (50°F to 104°F) / 15% to 85% RH / 800 to 1060 hPa
Storage/Transport Conditions	-20°C to +60°C (-4°F to 140°F) / 10% to 95% RH / 700 to 1060 hPa

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
UA-1030TCN Individual Unit (Box)	16.5 x 9.5 x 21 cm (6.5 x 3.8 x 8.3")	.81 kg (1.8 lbs)	0.004 m ³ (0.14 ft ³)	0-93764-60294-8
UA-1030TCN Case (10 units)	50.5 x 35.3 x 23.1 cm (19.9 x 13.9 x 9.1")	9.2 kg (20.3 lbs)	0.05 m ³ (1.73 ft ³)	1-00-93764-60294-5

Optional Accessories

Name	Order Code
AC Power Adapter	TB:233

Compatible Replacement Cuffs- Packaged in individual box

Small Cuff 16-24 cm (6.3-9.4")	UA-289A
Medium Cuff 23-37 cm (9.0-14.6")	UA-290A
Large Cuff 31-45 cm (12.2-17.7")	UA-291A
Wide Range Cuff 21.8-41.9 cm (8.6-16.5")	UA-420A
Smooth Fit Cuff 23-37 cm (9.0-14.6")	UA-392A

Included in the Box

- 1 Blood Pressure Monitor
- 1 Blood Pressure Cuff 23-37 cm (9.0-14.6")
- 1 Manual
- 1 Quick Start Guide
- 1 Carrying Case
- 4 AA Batteries
- AC Adapter

Wholesaler

	Codes
McKesson	017569
Kohl & Frisch	121175
Jean Coutu	560590
Familiprix	412023
McMahon	27020901

Health Canada License - 87256