

A simple and accurate blood pressure monitor providing measurements at the touch of a button and includes a small cuff for petite arms.

**Product Benefits**

- Irregular Heartbeat (IHB) indicates the presence of an irregular heartbeat
- Displays averages of all measurements
- Small cuff fits petite arms 6.3-9.4" (16-24 cm)
- Compatible with Small to Large cuffs (not included)
- Stores 60 measurements in memory with recall

**Comfort and Convenience**

- One button operation for quick and easy functionality
- Advanced inflation and deflation technology provides faster and more comfortable readings
- Adjustable pressure settings to ensure comfort
- Includes AC adapter

**The A&D Promise**


- Clinical Grade Accuracy you can trust
- 5 year warranty on monitor
- Cuffs are not made with natural rubber latex



Measurement Memory

Pressure Rating Indicator

Battery Life Indicator



**Unlock the Potential**  
Download the FREE  
**A&D Heart Track**  
Companion App  
for manual entry



Irregular Heartbeat  
Symbol (IHB)

One Button  
Operation

## Technical Specifications

<b>Measurement Method</b>	Oscillometric
<b>Memory</b>	60 measurements
<b>Cuff Size</b>	6.3-9.4" (16-24 cm)
<b>Measurement Range</b>	<b>Pressure:</b> 0-299 mmHg, <b>Systolic pressure:</b> 60-279 mmHg <b>Diastolic pressure:</b> 40-200 mmHg, <b>Pulse:</b> 40-180 beats/min
<b>Measurement Accuracy</b>	<b>Pressure:</b> ±3 mmHg, <b>Pulse:</b> ±5%
<b>Power Supply</b>	AC Adapter (TB:233) included or 4 x 1.5V batteries (R6P, LR6 or AA) sold separately
<b>Battery Life</b>	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries)
<b>Monitor Dimensions</b>	3.8 x 5.1 x 2.7" (9.6 x 13 x 6.8 cm)
<b>Monitor Weight</b>	8.5 oz (240 g) without batteries
<b>Display Size</b>	1.3 x 2.1" (3.4 x 5.3 cm)
<b>Clinical Test</b>	ISO81060-2:2013
<b>Operating Conditions</b>	+50°F to +104°F (10°C to 40°C) / 15% to 85% RH / 800 to 1060 hPa
<b>Storage/Transport Conditions</b>	-4°F to +140°F (-20°C to 60°C) / 10% to 95% RH / 700 to 1060 hPa

## Packaging and Shipping Specifications

	<b>Dimensions (WxLxH)</b>	<b>Weight</b>	<b>Case Cube</b>	<b>UPC / GTIN</b>
<b>Individual Unit Box</b>	4.9 x 3.5 x 7.1" (12.5 x 9 x 18 cm)	1.6 lb (0.7 kg)	0.07 ft <sup>3</sup> (0.002 m <sup>3</sup> )	0-93764-60439-3
<b>Case (10 units)</b>	18.7 x 10.5 x 7.6" (47.5 x 26.7 x 19.4 cm)	15.1 lb (6.8 kg)	0.87 ft <sup>3</sup> (0.03 m <sup>3</sup> )	1-00-93764-60439-0

## Optional Accessories

Name	Order Code
AC Power Adapter	TB:233

### Compatible Replacement Cuffs- Packaged in individual box

Small Cuff 6.3-9.4" (16-24cm)	UA-289A
Medium Cuff 9.0-14.6" (23-37 cm)	UA-290A
Large Cuff 12.2-17.7" (31-45 cm)	UA-291A
Wide Range Cuff 8.6-16.5" (22-42 cm)	UA-420A

## What is In The Box

- 1 Monitor
- 1 Blood pressure cuff (6.3-9.4")
- 1 Instruction Manual
- 1 Quick Start Guide
- 1 AC Adapter