

The A&D Medical Professional Blood Pressure Kiosk is a compact powerhouse ideal for small pharmacies, mobile clinics, and workplace screenings. Designed for durability and frequent use, it allows patients to easily self-measure their blood pressure and obtain accurate results.

**Product Benefits**

- Experience fast, accurate and reliable results with the unique cuff inflation.
- Indicates the potential presence of an irregular heartbeat with an icon on the print out providing an early warning of possible cardiovascular issues.
- One button operation for quick and simple measurements.
- Anti-bacterial cuff cloth provides a sanitary solution. The cuff cloth is easily replaced and disposable covers are available.

**Convenience**

- Integrated printer that can be customized allows the patient to keep results as a record or provide to a healthcare professional. Measurements can be uploaded to the Heart Track app via a QR code on the result printout.
- The small footprint is ideal for integration in a health screening station or tight spaces.
- Seamlessly collect data via the optional Bluetooth® or RS232C communication modules.
- No installation required, just remove from box and plug in and turn on power.

**The A&D Promise**

- Clinical Grade Accuracy delivers precise measurements that have been validated in a clinical setting.
- Extended 2-year warranty on monitor provides long-term peace of mind and greater value.



## Technical Specifications

<b>Measurement Method</b>	Oscillometric
<b>Cuff Size</b>	Fits Arms 7.1-13.8" (18-35 cm)
<b>Measurement Range</b>	<b>Pressure:</b> 0-299 mmHg, <b>Systolic pressure:</b> 40-270 mmHg <b>Diastolic pressure:</b> 20-200 mmHg, <b>Pulse:</b> 30-240 beats/min
<b>Measurement Accuracy</b>	<b>Pressure:</b> ±3 mmHg or 2%, <b>Pulse:</b> ±5%
<b>Display</b>	3 Digit Display LED & LED Lamp
<b>Time/Date</b>	12 Hour Format, Year, Month, and Day
<b>Power Supply</b>	100-240V, 50-60 HZ
<b>AC Power Consumption</b>	50-80 VA
<b>Clinical Test</b>	EN1060-4:2004
<b>Wired Communication</b>	Optional RS-232C (TM-2657P-RS)
<b>Wireless Connection</b>	Optional Bluetooth® (TM-2657P-BT)
<b>Safety Standard</b>	IEC60601-1:2012
<b>EMD Compliance</b>	Complies with Electromagnetic Disturbances Standard IEC60601-1-2:2014
<b>Operating Conditions</b>	+50°F to +104°F (10°C to 40°C) / 15% to 85% RH (no condensation)
<b>Storage/Transport Conditions</b>	-4°F to +140°F (-20°C to 60°C) / 10% to 95% RH (no condensation)

## Packaging and Shipping Specifications

	<b>Dimensions (WxLxH)</b>	<b>Weight</b>	<b>UPC</b>
<b>Individual Unit Box</b>	13.3 x 22.8 x 18.1" (34 x 60 x 46 cm)	19.55 lb (8.9 kg)	0-93764-60350-1
<b>Kiosk Dimensions</b>	9.5 x 18.8 x 12.8" (24 x 48 x 33 cm)	12.4 lb (5.6 kg)	N/A
<b>Stand Dimensions</b>	4.63 x 25.38 x 22.38" (12 x 65 x 57 cm)	18.9 lb (8.6 kg)	N/A

## Optional Accessories

<b>Name</b>	<b>Order Code</b>
Arm Cuff Cover (Pack of 5)	AX-134005759-S
Dedicated Stand	TM-ST520
Adjustable Height Stool	TM-STA001
Printer Paper (Pack of 5)	AX:PP147-S
Disposable Kiosk Sleeves (Pack of 400)	AS-134010367

## What is In The Box

- 1 Kiosk monitor
- 1 Power cable
- 1 Printer paper roll
- 1 Instruction manual
- 1 Patient instruction panel