

# **ONE TOUCH**

Blood Pressure Monitor Upper Arm UA-660

# A traditional blood pressure monitor that delivers quick and accurate readings at the touch of a button.

#### **Product Benefits**

- Easily track blood pressure over time with 60 measurement memory eliminating the need for written log
- Displays the average of the last three measurements providing more consistent blood pressure data

#### **Comfort and Convenience**

- · One button operation for quick and easy functionality
- Advanced inflation and deflation technology provides faster and more comfortable readings
- Includes a carrying case and batteries for immediate use and portability

#### The A&D Promise

- Clinical Grade Accuracy delivers measurements that are as precise as those taken in a medical setting
- Extended 5-year warranty on monitor provides long-term peace of mind and greater value





## **Technical Specifications**

Measurement Method	Oscillometric	
Memory	60 measurements	
Cuff Size	8.6-16.5" (22-42 cm)	
Measurement Range	Pressure: 0-299 mmHg, Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg, Pulse: 40-180 beats/min	
Measurement Accuracy	Pressure: ±3 mmHg, Pulse: ±5%	
Power Supply	4 x 1.5V batteries (R6P, LR6 or AA)	
Battery Life	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries)	
Monitor Dimensions	3.8 x 5.4 x 2.6 in (9.6 x 13.8 x 6.7 cm)	
Monitor Weight	9.2 oz (260g) without batteries	
Display Size (WxH)	1.4 x 2.1" (3.5 x 5.3 cm)	
Clinical Test	ISO81060-2: 2020	
Operating Conditions	+50°F to +104°F (10°C to 40°C) / 15% to 85% RH / 800 to 1060 hPa	
Storage/Transport Conditions	-4°F to +140°F (-20°C to 60°C) / 10% to 95% RH / 700 to 1060 hPa	

## **Packaging and Shipping Specifications**

	Dimensions (WxLxH)	Weight	Case Cube	UPC / GTIN
Individual Unit Box	4.9 x 3.5 x 7.1 in (12.5 x 9 x 18 cm)	1.4 lb (620 g)	0.07 ft <sup>3</sup> (0.002 m <sup>3</sup> )	0-93764-60449-2
Case (10 units)	10.6 x 18.7 x 7.9" (27 x 47.5 x 20 cm)	14.9 lb (6.8 kg)	0.91 ft <sup>3</sup> (0.03 m <sup>3</sup> )	1-00-93764-60449-9

## **Optional Accessories**

#### Compatible Replacement Cuffs- Packaged in individual box

Small Cuff 6.3-9.4" (16-24 cm)	UA-289A	
Medium Cuff 9.0-14.6" (23-37 cm)	UA-290A	
Large Cuff 12.2-17.7" (31-45 cm)	UA-291A	
Wide Range Cuff 8.6-16.5" (22-42 cm)	UA-420A	
Cuffs are not made with natural rubber latex		

### What is In The Box

- 1 Blood pressure monitor
- 1 Blood pressure cuff 8.6-16.5" (22-42 cm)
- 1 Instruction manual
- · 1 Quick start guide
- 1 Carrying case
- 4 AA batteries

