

Experience trusted reliability and accuracy with our advanced cellular enabled weight scale. Featuring real-time data transmission, it enhances compliance in remote patient monitoring.

**Product Benefits**

- Enhanced data set provides access to complete, reliable and precise information maximizing compliance
- Motion tolerance to ensure patient is stable before measurement enhancing accuracy of the measurement
- Zero Weight operation ensures accurate measurements despite changes in location or ambient temperature
- Convenient step on start making it simple for the patient to use
- Low battery indicator alerts you when to replace the battery maintaining functionality
- Batteries included for immediate use out of the box
- Can be used as a standalone scale with 90-reading memory for tracking weight over time, expanding device usability
- Audible and visual notifications to indicate to the when a measurement starts and completes.


**Support and Integration**

- MVNO service leverages multiple cellular networks to find the best signal strength offering reliable connectivity to ensure seamless transmission of vital health data
- API securely delivers full unique data set to cloud
- Integration support with complete API and guide
- Expert team provides in-country technical support, customer service, and account management, enhancing your operational efficiency and ensuring ongoing support



**The A&D Promise**

- A&D designs and manufactures products ensuring unparalleled quality
- Devices offer clinical grade accuracy, providing precise and dependable measurements essential for remote patient monitoring



## Technical Specifications

<b>Capacity</b>	450 lb / 200 kg
<b>Memory</b>	90 measurements
<b>Minimum Resolution</b>	0 to 220 lb with an increment of 0.2 lb 220 to 450lb with an increment of 0.5 lb 0 to 100 kg with an increment of 0.1 kg 100 to 200 kg with an increment of 0.2 kg
<b>Wireless Communication</b>	LTE™-M (LTE™-Cat.M1) NB-IoT (LTE™Cat.NB1)
<b>Setting Units</b>	Units: lb or kg
<b>Power Supply</b>	4 AAA size alkaline batteries (included)
<b>Scale Dimensions (WxDxH)</b>	11.8 x 11.8 x 1.2" (300 x 300 x 31 mm)
<b>Scale Weight</b>	4.7 lb (2.1kg) including batteries
<b>Display Size (WxH)</b>	3.72 x 1.16" (9.44 x 2.94 cm)
<b>Operating Conditions</b>	+50°F to +104°F (+10°C to + 40°C) / 30% to 85% RH
<b>Storage/Transport Conditions</b>	14°F to 140°F (-10°C to +60°C / 30% to 85% RH

## Packaging and Shipping Specifications

	<b>Dimensions (WxLxH)</b>	<b>Weight</b>	<b>Case Cube</b>	<b>UPC / GTIN</b>
<b>Individual Unit Box</b>	13.6 x 1.9 x 13.6 in (34.6 x 4.8 x 34.6 cm)	5.3 lb (2.4 kg)	0.20 ft <sup>2</sup> (0.01 m <sup>2</sup> )	0-93764-60445-4
<b>Case (5 units)</b>	14 x 10.2 x 14.6 in (33.6 x 24.6 x 35 cm)	27.6 lb (12.5 kg)	1.21 ft <sup>2</sup> (0.03 m <sup>2</sup> )	1-00-93764-60445-1

## What is In The Box

- 1 Weight scale
- 1 Instruction manual
- 1 Quick start guide
- 4 AAA batteries